
ottobock.








10V39

DE	Gebrauchsanweisung	6
EN	Instructions for use	13
FR	Instructions d'utilisation	19
IT	Istruzioni per l'uso	26
ES	Instrucciones de uso	32
PT	Manual de utilização.....	39
NL	Gebruiksaanwijzing.....	46
SV	Bruksanvisning.....	53
DA	Brugsanvisning.....	59
NO	Bruksanvisning.....	65
PL	Instrukcja użytkowania.....	71
CS	Návod k použití.....	78
TR	Kullanma talimatı.....	84
EL	Οδηγίες χρήσης	90
RU	Руководство по применению	98

①

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English

Date of last update: 2020-04-27

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

1 Scope of delivery

10V39=* MovoWrist Flex (2)

647G375 Instructions for Use
11D105 Protective Cover (1)

1.1 Accessories (not included)

10A30=1/2"-20 (3)

10A30=M12×1.5 (3)

10A30=1

2 Description

2.1 Indications for use

The 10V39=* MovoWrist Flex is to be used solely for the fitting of upper limb prostheses in conjunction with the Ottobock system hands or Ottobock cable-activated hooks.

The wrist unit allows individual flexion and extension, locking in various positions in increments of 15° as well as rotation of the terminal device.

The wrist unit has been provided with a quick-release lock for an easy exchange of the terminal devices.

The 10A30=* Adapter necessary for the wrist unit depends on the type of connection of the terminal device and is therefore not included in the delivery of the wrist unit. The various movements of the wrist unit allow for a physiologically correct posture when using the prosthesis and reduce the necessity of unnatural compensating movements.

3 Technical data

The flexion wrist unit has 5 setting positions with the following angles:

Neutral or zero setting	0°	
Extension setting	-15°	
Flexion settings	15°/30°/45°	
Rotation	20 locked positions every 18°	
Weight	10V39=50	about 130 g
	10V39=45	about 110 g
	10A30=1/2"-20	about 20 g
	10A30=M12×1,5	about 20 g
	10A30=1	about 58 g
Overall length	33 mm, of which 12 mm are visible on the distal end of the forearm socket.	

With the 10A30= M12×1,5 und 10A30=1/2"-20 Adapter mounted, the visible length of the distal part increases to 17 mm.

4 Safety instructions

Attention!

- To guarantee a solid bond, clean the respective areas with acetone before applying the adhesive. Do not touch the areas after you have degreased them and when the acetone has dried up, apply the adhesive immediately.
- Follow all fabrication and safety considerations of the adhesive manufacturer.
- Please follow all recommended torque values.
- Use the 'Patient Information' section of this pamphlet (section 5) to inform your patient about the function and use of the Mo-voWrist Flex.

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water, and acids. Using this medical product under such conditions will render all claims against Otto Bock HealthCare null and void.

Please also inform your patient.

5 Patient information

The 10V39=* MovoWrist Flex has been developed for everyday life activities and must not be used for unusual activities such as extreme sports (free climbing, parachute jumping, paragliding, etc.). Careful handling of the prosthesis and its components not only increases the service life of the product but, above all, ensures the patient's safety. Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the prosthesis inspected for any damage. If necessary, the prosthetist will pass the prosthesis on to the Ottobock Service. Any damaged components may only be repaired by the Ottobock Service.

Attention!

- Please note that the position of terminal device may change when you operate the flexion or rotation push buttons unintentionally.

- Please note that the terminal device may open when you operate the unlock push button unintentionally.
- Before starting work and during work, check the wrist lock for proper function.

5.1 Operating a motor vehicle

An upper limb amputee's ability to drive a motor vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities. All persons are required to observe their country's national and state driving laws when operating vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test center. For maximum safety and convenience, Ottobock recommends that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork).

Caution!

Before starting the car, check the 10V39=* MovoWrist Flex for proper lock and secure seat.

6 Function

The 10V39=* MovoWrist Flex (fig. 1, item 2) allows for a fitting in conjunction with an Ottobock system hand or Ottobock cable-activated hook or another terminal device.

It features flexion and extension settings in different positions. This joint is also rotatable on the forearm socket in steps of 18°. To allow for an easy change of the terminal devices in every position, the 10V39=* MovoWrist Flex is equipped with a quick-release lock.

Three push-buttons are assigned to different functions (fig. 3). These are used for setting the rotation (1), exchanging system hands or a cable-activated hook (2) as well as for flexion and extension (3).

The wrist unit is available in two sizes and can be used on both sides:

- 10V39=45 for forearm socket with inner diameter for wrist unit connection of 45 mm
- 10V39=50 for forearm socket with inner diameter for wrist unit connection of 50mm

6.1 Rotation

The 10V39=* MovoWrist Flex unlocks for rotation by pressing the push button 1 (fig. 3). It can be freely positioned in 20 different engaging points. When the push button 1 is released in an engaging position, the wrist unit will lock.

6.2 Exchange

The push button 2 (fig. 3) unlocks the wrist unit for exchange of the system hands, cable-activated hook or terminal device. The system hand, hook or terminal device ejection is supported by a spring integrated in the 10A30=* Adapter.

Should a new terminal device be adapted, it must be slid into the wrist unit while locking is performed automatically.

6.3 Flexion

The push button 3 (fig. 3) releases flexion and extension of the 10V39=* MovoWrist Flex for positioning in 5 lock positions. After releasing the push button 3, the terminal device will lock in one of the flexion lock positions.

7 Mounting instructions MovoWrist Flex 10V39

7.1 Mounting MovoWrist Flex 10V39

The mechanics of the wrist unit must be separated from the lamination ring with Allen wrench size 3 (fig. 4). In order to avoid soiling the lamination ring in the subsequent steps, replace the included protective cover with slight pressure (fig. 5). The flexion level of the mechanics is indicated on the lamination ring by a notch (fig. 2). Bring the lamination ring in the forearm socket into the correct position according to the course of flexion. Mark the position of the notch on the forearm socket (fig. 6). To create a firm connection with the lamination ring, roughen the bonding surface on the inside of the forearm socket. Mark the 20-mm height of the lamination ring on the lower arm (fig. 7). Clean the bonding surfaces on the lamination ring and lower arm with acetone. After cleaning, let the acetone air dry and do not touch the bonding surfaces. Mask the outside of the lower arm off using Polyethylene Adhesive Tape 627B4 to avoid soiling it while gluing. Mix Orthocryl® Sealing Resin Compact Glue 636K18=1 and Hardener Paste 617H14 in accordance with the manufacturer's specifications. Apply sealing resin on the lamination ring and lower arm (fig. 8) and flush-align the lamination ring with the lower arm. Make sure that the notches in the lamination ring are flush with the marking on the lower arm made

previously. After hardening the sealing resin (for hardening time see the manufacturer's specifications), push out the protective cover (fig. 9). Make sure that the sealing resin is spread over the entire surface between the lamination ring and the forearm socket and that no sealing resin penetrates into the bore holes (fig. 10). After the removal of the protective cover, lightly apply the 636K14 Loctite® 601 on the lower third of the locking screw of the mechanics and screw into the lamination ring with a torque of 5 Nm (fig. 11). Remove excessive Loctite®.

7.2 Checking and adjusting the rotation function (fig. 14)

① Checking the rotation

The product is factory-set so that the lock pin ④ is drawn out of the lamination ring ② when the push button is pressed ③, thus allowing unimpeded rotation.

Since the adjustment of the locking pin may lead to loss of functionality, an inspection must be carried out following installation in the prosthesis.

Check the proper function of the rotation as shown in the illustration (① Pressing the push button / ② Checking the rotation).

⑤ Adjusting the lock pin

Adjust the lock pin as shown in the illustration ③. In doing so, make only a small change (1 rotation) before checking whether the rotation is functioning properly ④. The lock pin must be adjusted so that the lamination ring engages and disengages securely.

7.3 Mounting the adapter 10A30

Mounting the adapter 10A30=M12×1.5/10A30=1/2"-20

Clean the thread surfaces of the terminal device and 10A30=* Adapter with acetone. Let the acetone air dry and do not touch the thread surfaces. Apply the 636K14 Loctite® 601 on the lower third of the thread surface of the terminal device (fig. 12). Screw the terminal device completely into the 10A30=* Adapter. Remove excessive Loctite®. Mount the terminal device with adapter into the wrist unit and align the terminal device to the lower arm, if necessary, by unscrewing from the adapter (by a maximum of ¼ turn) (fig. 13). Remove the terminal device with the 10A30=* Adapter from the wrist unit again and remove excessive Loctite® from the terminal device or wrist unit. Let the adhesive air dry according to the manufacturer's specifications.

Mounting the adapter 10A30=1

Adapter for using the 10V39 in conjunction with terminal devices with standard pins.

Insert the standard pin of the terminal device into the standard connection of the 10A30=1 and secure the closure by turning the clamp sleeve to the right. Mount the terminal device with adapter into the wrist unit.

7.4 Fabricating a laminated socket

Use the 743A18 Alignment Tool with wrist unit diameter 45 mm or 50 mm during the fabrication of the laminated socket. For assembly of the 10V39 MovoWrist Flex, please follow the instructions under section 7.1.

8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.3 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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