

## 9X18

<b>DE</b> Gebrauchsanweisung (Fachpersonal) .....	3
<b>EN</b> Instructions for Use (Qualified Personnel) .....	7
<b>FR</b> Instructions d'utilisation (Personnel spécialisé) .....	11
<b>IT</b> Istruzioni per l'uso (Personale tecnico specializzato) .....	16
<b>ES</b> Instrucciones de uso (Personal técnico especializado) .....	21
<b>PT</b> Manual de utilização (Pessoal técnico) .....	25
<b>NL</b> Gebruiksaanwijzing (Vakmensen) .....	30
<b>SV</b> Bruksanvisning (Fackpersonal) .....	34
<b>DA</b> Brugsanvisning (Faguddannet personale) .....	38
<b>NO</b> Bruksanvisning (Fagpersonell) .....	43
<b>EL</b> Οδηγίες χρήσης (Τεχνικό προσωπικό) .....	47
<b>JA</b> 取扱説明書 (有資格担当者) .....	52

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---

English

## 1 Timeliness of the Instructions for Use

### INFORMATION

Date of the last update: 2013-01-22

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

## 2 Product Description

### 2.1 Description

The 9X18 Cable Pull Switch is suitable for incorporation in the definitive socket. The 9X18 Cable Pull Switch is suitable for patients who cannot generate enough muscle signal to control a myoelectric prosthesis but who want to take advantage of the high rehabilitation value of myoelectric prostheses.

If the cable pull switch is mounted in the provisional socket, the patient can practise the control or switching as early as the first stage of the rehabilitation. This will ensure a smooth transition to the use of the definitive prosthesis.

The cable pull switch is used to control the prosthetic components with digital control and to switch between the prosthetic components. It can be

used in combination with MYOBOCK electrodes, linear transducer and 4-step transducer.

## 2.2 Function

The cable pull switch features 2 control points that are controlled independently and one after the other. A function for control of the prosthetic components is assigned to each control point. The control point 1 should always be assigned to the 'close hand' function.

Pulling the cable of the cable pull switch wire will reach the first control point (S1) and activate control signal A. When the cable is pulled again, the second control point (S2) is reached, control signal A drops off and control signal B is activated. Now the control switch B is solely active.

## 3 Application

### 3.1 Indication

The 9X18 Cable Pull Switch is to be used **solely** in exoprosthetic fittings of the upper extremities and for myoelectric control of arm prostheses.

### 3.2 Application

The 9X18 Cable Pull Switch is suitable for both right or left arm amputees as well as bilateral amputees.

### 3.3 Conditions of Use

The 9X18 Cable Pull Switch is suitable for indoor and outdoor use.

### 3.4 Qualification

The fitting of patients with the 9X18 Cable Pull Switch may only be carried out by a prosthetist.

## 4 Safety Instructions

Failure to follow the safety instructions given below can lead to faulty control or malfunction of the 9X18 Cable Pull Switch and result in a risk of injury for the patient.

## 5 General/Patient-Specific Safety Instructions

### 5.1 General Safety Instructions

- Check the pull-cable regularly for wear and tear, and the cable pull switch housing for damage.
- When mounted into a provisional socket, the cable pull switch may only be worn by an amputee within an orthopaedic workshop or as an inpatient in a rehabilitation facility.
- Please note that the cable travel of the cable pull switch is limited. To prevent excessive strain on the cable pull switch, make sure the harness is correctly fitted and accommodates the patient's movements.
- The cable pull switch may only be loaded in the pulling direction.
- When mounting a cable with cable housing, make sure there are no tight angles. Avoid kinks. Make sure the cable runs smoothly.
- The cable pull switch may not be used for bathing prostheses.
- Instruct the patient in the proper handling of the cable pull switch according to the section "Patient-Specific Safety Instructions".

### 5.2 Patient-Specific Safety Instructions

- Do not increase your tension on the cable once you have felt the cable come to a stop within the cable pull switch.
- Handle the cable pull switch with care and regularly check the product for damage (e.g. cracks or breakage) at the pull-cable or housing. Should the cable pull switch be damaged or not function properly, please send it immediately to your Ottobock service facility.
- Driving vehicles of any kind or operating machinery is strictly prohibited while wearing a provisional socket fitted with a cable pull switch.
- Do not let foreign particles or liquids get into the cable pull switch. The cable pull switch should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- The cable pull switch and damaged components may only be opened or repaired by certified Ottobock Myo-Service technicians.
- Regularly clean the cable pull switch with a cloth **slightly** dampened with a soap solution. Do not allow moisture to enter into the cable pull switch or into the contacts.
- Since all moving parts are subject to wear and tear, annual service is required.

## 6 Delivery

### 6.1 Scope of Delivery

- 1 pc. 9X18 Cable Pull Switch
- 2 pcs. 501S46=M3x5 Oval head screw
- 2 pcs. 501S46=M3x8 Oval head screw
- 1 pc. 647G401 Instructions for Use

## 7 Preparation for Use

### 7.1 Connecting the Electrode Cable

The cable pull switch is delivered without connecting cable. It must be ordered separately. Appropriate cables are 13E50=\*, 13E97=\*, 13E98=1200 and 13E99=1200 connecting cables. Abundantly apply 633F11 Silicone Grease before connecting the contacts of the connecting cable. The grease will protect the connectors from sweat and humidity.

### 7.2 Mounting the Pull-Cable to the Stirrup

The cable pull switch is delivered with pre-assembled pull-cable. When mounting the cable pull switch to a provisional or definitive socket you must connect the free end of the pull-cable to the stirrup.

The pull-cable of the cable pull switch can also be fed through a cable housing to the shoulder harness. When positioning the cable housing, make sure there are no tight angles, edges or kinks and that the cable runs smoothly.

- 1) Determine the necessary cable length, but do not yet shorten the cable.
- 2) Form the pull-cable into a loop and feed it into the stirrup (Fig. 1).
- 3) Place the wedge into the stirrup and tighten the pull-cable (Fig. 2).

Then check the length of the cable and adjust if necessary. You must first remove the wedge before doing this. After the correct cable length is set, cut the pull-cable to length.

## 8 Disposal

This product must not be disposed of with regular domestic waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

## 9 Legal Information

### 9.1 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this device may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise. The device may only be opened and repaired by authorised Ottobock technicians.

### 9.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

### 9.3 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

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