








## 3R80, 3R80=ST



<b>DE</b> Gebrauchsanweisung .....	7
<b>EN</b> Instructions for use .....	19
<b>FR</b> Instructions d'utilisation .....	32
<b>IT</b> Istruzioni per l'uso .....	45
<b>ES</b> Instrucciones de uso .....	58
<b>PT</b> Manual de utilização .....	71
<b>NL</b> Gebruiksaanwijzing .....	84
<b>SV</b> Bruksanvisning .....	97
<b>PL</b> Instrukcja użytkowania .....	109
<b>HU</b> Használati utasítás .....	122
<b>CS</b> Návod k použití .....	134
<b>RO</b> Instrucțiuni de utilizare .....	147
<b>TR</b> Kullanma talimatı .....	159
<b>RU</b> Руководство по применению .....	171
<b>JA</b> 取扱説明書 .....	185
<b>ZH</b> 使用说明书 .....	197

## INFORMATION

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### INFORMATION

Date of last update: 2021-06-28

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.

- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

## 1.1 Construction and Function

The 3R80\* Prosthetic Knee Joint is a monocentric prosthetic knee joint with rotation hydraulics and a mechanically controlled stance phase stabilisation system. It is waterproof and equipped with a manually operated lock. The prosthetic knee joint has a pyramid adapter (3R80) or threaded connector (3R80=ST) as the proximal connection.

The integrated rotation hydraulics control the swing phase and secure the stance phase. Damping for flexion, extension and the stance phase are individually adjustable. This allows comfortable walking over a broad range of speeds.

Stance phase damping is activated when bearing weight on the prosthesis. The activation threshold is adjustable. It provides for stance phase flexion up to 4° while bouncing, and more than 4° while yielding. Walking on slopes and down stairs step-over-step is therefore possible.

An extension assist spring supports the swing phase by generating an extension moment. The additional spring included in the scope of delivery can be installed to increase the extension moment.

The lock keeps the prosthetic knee joint at full extension. This function is suitable for situations where flexion is to be prevented (e.g. when standing for extended periods of time or wading through water).

## 1.2 Combination possibilities

### INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

Prosthetic components that are especially well suited for combination with the product are listed below.

Designation	Reference Number
Prosthetic Feet	1A30, 1C30, 1C40, 1D35, 1E56, 1C60, 1C61, 1C62 <b>Waterproof:</b> 1C63, 1C64
Foam Cover	3S107

In a waterproof walking aid, the product can be combined with prosthetic components from the Ottobock Aqualine product line.

### Unallowable combinations

Prosthetic hip joints	7E10*
-----------------------	-------

## 2 Intended use

### 2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

## 2.2 Area of application

3R80\*

Recommended for mobility grades **3 and 4** (unrestricted outdoor walker and unrestricted outdoor walker with particularly high demands). Approved for a body weight of up to **150 kg**.

## 2R57=16-285-WF

Approved for a body weight of up to **150 kg**.

## 2.3 Environmental conditions

### Allowable environmental conditions

Temperature range for use: -10 °C to +60 °C

Storage and transport conditions: -10 °C to +60 °C, 250 mbar to 1013 mbar, condensing

Acid resistance: salt water, chlorine water, soapsuds, perspiration, urine

Moisture: fresh water, salt water, chlorine water, soapsuds, steam

Allowable immersion depth: 3 m

### Unallowable environmental conditions

Mechanical vibrations or impacts

Dust, sand, highly hygroscopic particles (e.g. talcum)

Sauna

## 2.4 Lifetime

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

## 3 Safety

### 3.1 Explanation of warning symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

#### **CAUTION**

#### **Excessive strain on the product**

Risk of injury due to breakage of load-bearing components

- ▶ Use the prosthetic components according to the MOBIS classification.

#### **CAUTION**

#### **Improper combination of prosthetic components**

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

**⚠ CAUTION**

**Use under restricted environmental conditions**

Risk of injury due to damage to the product

- ▶ Do not expose the product to restricted environmental conditions.
- ▶ If the product has been exposed to restricted environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

**⚠ CAUTION**

**Exceeding the service life**

Risk of injury due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded.

**⚠ CAUTION**

**Reuse on other patients and improper maintenance**

Fall due to loss of functionality or damage to the product

- ▶ Only use the product for a single patient.
- ▶ Maintain the product regularly to extend the service life of the product (see the section "Maintenance").

**⚠ CAUTION**

**Reaching into the area of the joint mechanism**

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

**⚠ CAUTION**

**Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

**Signs of changes in or loss of functionality during use**

Noticeable changes in functionality may include poor response, incomplete extension, diminished swing phase control or stance phase stability, abnormal noises, etc.

**4 Scope of delivery**

The scope of delivery for the product is shown on page 2 (see figure 1).

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

Fig.		Item no.	Quantity	Designation	Reference Number
-	■	-	1	Instructions for Use	647G403=1

Fig.		Item no.	Quantity	Designation	Reference Number
–	■	–	1	Quick Start Guide	647H515=1
1	–	1	1	Prosthetic Knee Joint	–
1	■	2	1	Tube Adapter, titanium, waterproof	2R57=WF
1	■	3	1	Adjustment Wrench	710H10=2x3
1	■	4	1	Compression Spring	513D59=1.3X8.8X86
1	■	5	1	Pyramid Dome Protektor Kit	2Z11=KIT

## 5 Preparing the product for use

### 5.1 Information on fabrication of a prosthesis

#### CAUTION

##### **Incorrect alignment, assembly or adjustment**

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

#### CAUTION

##### **Initial use of the prosthesis by the patient**

Fall due to lack of patient experience, incorrect alignment or incorrect adjustment of the prosthesis

- ▶ For the safety of the patient, use a suitable device (e.g. parallel bars, handrail) during initial standing and walking.
- ▶ Do not use devices with insufficient stability (e.g. crutches or walking sticks).

Fabricating a prosthesis that is optimally tailored to the patient demands technical knowledge of orthopaedic technology. The alignment and adjustments may only be carried out by a prosthetist.

The alignment and adjustments are mainly completed in the following steps:

1. Bench alignment
2. Static alignment
3. Dynamic trial fitting

These steps are first carried out with a test prosthesis in order to determine the best possible combination and positioning of the prosthetic components relative to each other. Once all settings have been tailored to the patient, the definitive prosthesis is produced. The same steps are repeated for this purpose.

In dynamic trial fitting, the patient must be given sufficient time to become familiar with the functions of the prosthesis, thereby learning safe use for everyday life.

### 5.2 Information on Using a Cosmetic Cover

#### CAUTION

##### **Use of talcum**

Risk of injury, damage to the product due to lack of lubrication

- ▶ Do not use talcum on the product or other prosthesis components.

- ▶ Only use the foam cover approved by the manufacturer (see section "Combination Possibilities" – see page 20).
- ▶ To reduce friction and to eliminate noise, apply 519L5 Silicone Spray directly onto the contact surfaces of the cosmetic foam cover.
- ▶ After fabricating the cosmetic cover, repeat the dynamic trial fitting since the cover affects the adjustments of the prosthesis (e.g. damping).

### 5.3 Bench Alignment

#### Use 2Z11=KIT

**INFORMATION:** The connection area of the prosthetic knee joint can be protected against scratches with the protective film of the 2Z11=KIT during alignment in the workshop and testing in the trial fitting area.

- ▶ Use the protective film as illustrated in the accompanying document for the 2Z11=KIT.
- ▶ Remove the protective film before the patient leaves the fitting area.

<b>Bench Alignment Process</b>	
<b>+ = anterior placement / - = posterior placement (in relation to the alignment reference line)</b>	
<b>No.</b>	<b>see fig. 4</b>
	<p><b>Required equipment:</b> Alignment apparatus (e.g., 743L200 L.A.S.A.R. Assembly or 743A220 PROS.A. Assembly)</p> <p>Position the prosthetic foot in the alignment apparatus in accordance with the following values:</p>
<b>①</b>	a-p positioning, middle of the prosthetic foot to the alignment reference line ①: <b>+30 mm</b>
<b>②</b>	<p>Heel height (see section "Determining the Effective Heel Height" – see page 25): <b>effective heel height (x) + 5 mm</b></p> <p>External foot rotation: <b>approx. 5°</b></p>
	<p>Check the prosthetic knee joint according to the section "Checking the Factory Settings" (see page 27).</p> <p>Position the alignment reference point ② (rotation axis) of the prosthetic knee joint in the alignment apparatus according to the following values:</p>
<b>③</b>	<p>Height of the alignment reference point in the alignment apparatus: <b>MTP-floor measurement + 20 mm</b></p> <p>a-p positioning of the alignment reference point to the alignment reference line: <b>0 mm</b></p> <p>External rotation of the prosthetic knee joint: <b>approx. 5°</b></p>
<b>④</b>	<p><b>INFORMATION: If the length-adjustable 2R45=34 Tube Adapter (max. body weight 125 kg) is used in the alignment of the test prosthesis, the required length for the definitive prosthesis is easier to determine.</b></p> <p>When using the included tube adapter, note the section "Adapting the Length of the Tube Adapter" (see page 25).</p> <p>Connect the prosthetic knee joint to the prosthetic foot and chosen adapters according to their instructions for use and the section "Installing the Tube Adapter" (see page 25).</p>
<b>⑤</b>	<p>Mark the socket reference point and tuberosity reference point on the prosthetic socket according to the section "Determining the Socket Reference Point" (see page 26).</p> <p>Position the tuberosity reference point ④ of the prosthetic socket at the height of the tuberosity-floor measurement.</p>
<b>⑥</b>	<p><b>CAUTION! The correct position and flexion of the prosthetic socket is crucial for the safety of the patient and functionality of the prosthetic knee joint, and avoids premature wear.</b></p> <p><b>INFORMATION: The 4R118 Adapter Plate (max. body weight 125 kg) can be used for further posterior placement.</b></p> <p>Position the prosthetic socket so the socket reference point ③ lies on the alignment reference line ①.</p> <p>Set the determined socket flexion angle <math>\beta</math> (see section "Determining Socket Flexion" – see page 26).</p>
<b>⑦</b>	<p>Position the lamination anchor or socket adapter so that all prosthetic components between the socket and prosthetic knee joint can be connected to each other correctly.</p> <p>Join the prosthetic socket and the prosthetic knee joint using the adapters you selected.</p> <p>Be sure to follow the adapter instructions for use during adjustment and installation.</p>

### 5.3.1 Determining the Effective Heel Height

> **Required measuring equipment:**

743S12 Heel Height Measuring Device (see fig. 2)

- 1) Measure the posterior heel height ( $x_1$ ) of the shoe.
- 2) Measure the anterior heel height ( $x_2$ ) of the shoe.
- 3) Determine the effective heel height ( $x$ ) using the formula ( $x = x_1 - x_2$ ).

### 5.3.2 Shortening the Tube Adapter

**⚠ CAUTION**

**Incorrect processing of tube**

Fall due to damage to the tube

- ▶ Do not clamp the tube into a vice.
- ▶ To shorten the tube, use only a tube cutter or a cutting device.

> **Recommended tools and materials:**

719R5 Tube Cutter, 718R1 Tube Deburrer

- 1) **CAUTION! The length of the tube must not be less than the required length for the patient, since length compensation in the insertion area of the prosthetic knee joint is prohibited.**

Use the tube cutter to shorten the tube to the required length (see fig. 7).

- 2) Carefully deburr the inside and outside of the tube with the tube deburrer (see fig. 8).
- 3) Dispose of the remainder.

### 5.3.3 Installing the Tube Adapter

**⚠ CAUTION**

**Incorrect mounting of the tube**

Risk of injury due to breakage of load-bearing components

- ▶ Slide the tube all the way to the stop in the intended prosthetic component when mounting.

**⚠ CAUTION**

**Improper assembly of the screw connections**

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

> **Recommended tools and materials:**

710D4 Torque Wrench, degreasing cleaner (e.g. 634A3 Acetone)

- 1) Loosen the cap screw by **approx. 2 turns** (see fig. 10).
- 2) Clean the contact surfaces of the prosthetic knee joint and tube adapter using a degreasing cleaning agent.
- 3) Slightly rotate the tube as you insert it into the prosthetic knee joint to the stop (see fig. 9).
- 4) Using the torque wrench, tighten the cap screw (**10 Nm**) (see fig. 10).

### Installing the Prosthetic Foot

> **Recommended tools and materials:**

710D4 Torque Wrench, 636K13 Loctite®

- 1) Position the pyramid receiver of the tube adapter on the pyramid adapter of the prosthetic foot.
- 2) **For final assembly:**  
Use Loctite® to secure the set screws of the pyramid receiver.

- 3) Screw in the set screws.
- 4) Use the torque wrench to tighten the set screws (**10 Nm**).
- 5) Replace any set screws that are protruding or recessed too much with suitable ones (see selection table).
- 6) **For final assembly:**  
Use the torque wrench to retighten the set screws (**15 Nm**).

Selection table for set screws	
Reference number	Length (mm)
506G3=M8X12-V	12
506G3=M8X14	14
506G3=M8X16	16

### 5.3.4 Determining Socket Flexion

> **Required measuring equipment and materials:**

662M4 Goniometer, Pin

- 1) On the lateral side of the residual limb, mark the course of the thigh bone.
- 2) With the patient lying down, instruct the patient to fully bend the contralateral leg and hold it with the hands.
- 3) Slide the flat hand ① under the patient's back (see fig. 3).
- 4) Bring the residual limb to full flexion.
- 5) Slowly move the residual limb in the extension direction until pressure on the hand is reduced.
- 6) Slowly move the residual limb back in the flexion direction until the pressure reduction stops.
- 7) In this residual limb position, measure the angle  $\alpha$  (hip flexion contracture) with the goniometer.
- 8) **INFORMATION: Correct socket flexion ( $\beta$ ) during bench alignment is crucial for optimum walking. Take the individual situation of the patient into account for the tolerance.**  
If the angle  $\alpha = 0^\circ$ , the optimum value of the angle  $\beta = 3$  to  $5^\circ$  to the alignment reference line.  
If the angle  $\alpha > 0^\circ$ , the optimum value of the angle  $\beta = \alpha + 5$  to  $10^\circ$ .

### 5.3.5 Determining the Socket Reference Point

> **Required measuring equipment and materials:**

743A80 50:50 Gauge, Pin

- 1) Position the 50:50 gauge on the lateral side of the prosthetic socket in the distal area and mark a centre point ① (see fig. 5).
- 2) Position the 50:50 gauge in the proximal area and mark a centre point ②.
- 3) Using these points, mark the centre line of the prosthetic socket.
- 4) Clearly mark the tuberosity height as the tuberosity reference point ④ on the centre line.
- 5) Clearly mark the socket reference point ③ **30 mm** proximal to the tuberosity reference point.

## 5.4 Static Alignment

Static Alignment Process	
+ = anterior placement / - = posterior placement (in relation to the load line)	
No.	see fig. 6
	<b>Required equipment:</b> 743L100 L.A.S.A.R. Posture, 3D L.A.S.A.R. Posture 743L500
①	To determine the load line ①, position the patient on the L.A.S.A.R. Posture as follows: <ul style="list-style-type: none"> <li>• Prosthetic foot (with shoe) on force measurement plate (sufficient load: &gt; <b>35%</b> body weight)</li> <li>• Other foot (with shoe) on height compensation plate</li> <li>• The tips of the shoes are in one line</li> </ul>

Static Alignment Process	
+ = anterior placement / - = posterior placement (in relation to the load line)	
No.	see fig. 6
②	Optimise the static alignment solely by changing the plantar flexion. Adjustment is performed only with the anterior and posterior set screws of the tube adapter on the prosthetic foot. a-p positioning of the alignment reference point ② (rotation axis of the prosthetic knee joint) to the load line: <b>-35 mm</b>

## 5.5 Dynamic Trial Fitting

 <b>CAUTION</b>
<p><b>Adjusting the Settings</b></p> <p>Fall due to incorrect or unfamiliar settings</p> <ul style="list-style-type: none"> <li>▶ Only adapt the settings to the patient gradually.</li> <li>▶ Explain the effects of the adjustments on the use of the prosthesis to the patient.</li> </ul>

During dynamic trial fitting, the alignment and settings of the prosthesis are checked and adapted for optimum walking according to the needs and abilities of the patient.

The patient has to learn the safe use of the prosthesis through intensive training.

The following subsections describe the adjustment possibilities for adapting the product to the patient.

The list below provides an overview of the sequence for working through the subsections:

- Checking the Factory Settings
- Setting Flexion Damping
- Setting Extension Damping
- Setting Stance Phase Damping
- Setting the Stance Phase Damping Threshold

The adjustment is illustrated graphically in the supplied Quick Start Guide.

### 5.5.1 Opening the Covers

> **Tools required:**

Flat screwdriver

- 1) Open the extension damping and flexion damping covers (see fig. 11).
- 2) Open the stance phase damping cover (see fig. 12).

### 5.5.2 Checking the Factory Settings

Function	Means of adjustment	Factory setting	Meaning
Flexion damping	Adjustment valve "F" (see fig. 13)	Stop "-"	Minimum damping
Extension damping	Adjustment valve "E" (see fig. 14)	Stop "-"	Minimum damping
Stance phase damping	Adjustment ring "S" (see fig. 15)	Stop "+"	Maximum damping
Stance phase damping threshold	Adjustment ring "!" (see fig. 16)	Bumper without play, pre-stressed	Low threshold (stance phase damping is activated at low load)

> **Tools required:**

710H10=2x3 Adjustment Wrench

- 1) Check whether the settings match the factory settings (see table).
- 2) In case of deviations, restore the factory settings using the adjustment wrench.

### 5.5.3 Setting Flexion Damping

Function	Means of adjustment	Adjustment	Meaning
Flexion damping	Adjustment valve "F" (see fig. 13)	Turn to the right "+"	Increase damping
		Turn to the left "-"	Reduce damping

> **Tools required:**

710H10=2x3 Adjustment Wrench

▶ **CAUTION! Even at slow walking speeds, the prosthetic foot has to swing through adequately.**

Adjust flexion damping using the adjustment wrench so that the prosthetic foot does not swing through too far, according to the patient's walking speed.

### 5.5.4 Setting Extension Damping

Function	Means of adjustment	Adjustment	Meaning
Extension damping	Adjustment valve "E" (see fig. 14)	Turn to the right "+"	Increase damping
		Turn to the left "-"	Reduce damping

> **Tools required:**

710H10=2x3 Adjustment Wrench

▶ **CAUTION! The prosthetic knee joint has to reach full extension, even at slow walking speeds. The extension assist force can be increased if required (see section "Installing the Compression Spring" – see page 29).**

Using the adjustment wrench, adjust extension damping so that the prosthetic knee joint does not swing too hard against the extension stop.

### 5.5.5 Setting Stance Phase Damping

Function	Means of adjustment	Adjustment	Meaning
Stance phase damping	Adjustment ring "S" (see fig. 15)	Turn to the right "+"	Increase damping
		Turn to the left "-"	Reduce damping

> **Required training location:**

Stairs and ramp with handrail (e.g. 758G8=1 Test Platform)

**Tools required:**

710H10=2x3 Adjustment Wrench

- 1) Place the prosthetic foot on the bottom step so that the middle of the foot goes slightly beyond the step's edge.
- 2) Briefly lock the prosthetic knee joint by extending the residual limb.
- 3) Initiate flexion while bearing weight.
- 4) **INFORMATION: In case of problems switching between the swing phase and stance phase, readjust the stance phase damping threshold (see section "Setting the Threshold" – see page 28).**

Adjust stance phase damping so that the prosthetic knee joint does not flex too quickly or slowly while walking down stairs.

- 5) Once stance phase damping and the threshold are correctly adjusted for walking down stairs, check the settings while walking down a ramp and adjust if necessary.
- 6) Harmonise the settings for walking down stairs and walking down a ramp according to the needs of the patient.

### 5.5.6 Optional: Setting the Threshold

Function	Means of adjustment	Adjustment	Meaning
Stance phase damping threshold	Adjustment ring "I" (see fig. 16)	Turn to the right "+" (tension is applied to bumper)	The threshold is increased – a higher load is required to activate stance phase damping

Function	Means of adjustment	Adjustment	Meaning
Stance phase damping threshold	Adjustment ring "!" (see fig. 16)	Turn to the left "-" (tension on the bumper is released)	The threshold is decreased – a smaller load is required to activate stance phase damping

> **Tools required:**

710H10=2x3 Adjustment Wrench

- ▶ Use the adjustment wrench to adapt the stance phase damping threshold to the needs of the patient.

**5.5.7 Optional: Installing the Compression Spring**

The additional compression spring is installed in the lower joint section if the patient needs more support in order to bring the prosthetic knee joint to full extension.

> **Required tools and materials:**

Torque wrench (2 Nm), 513D59=1.3X8.8X86 Compression Spring

- 1) Unscrew the spring retainer from the lower joint section (see fig. 17).
- 2) Insert the compression spring into the larger compression spring in the lower joint section (see fig. 18).
- 3) Screw the spring retainer into the lower joint section (installation torque: **2 Nm**).
- 4) Check the settings of the prosthetic knee joint and reaching extension with a walking test.

**5.5.8 Closing the Covers**

- 1) Close the extension damping and flexion damping covers (see fig. 19).
- 2) Close the stance phase damping cover (see fig. 20).

**5.6 Finishing the Prosthesis**

> **Required tools:**

Torque wrench

- 1) Tighten all screws of the prosthetic components to the specified installation torque values.
- 2) Verify proper functioning of the prosthesis.

**6 Use**

**6.1 Information for use**

**⚠ CAUTION**

**Overheating of the hydraulics due to excessive strain**

Burns, injuries caused by falling due to changes in functionality and damage to prosthetic components

- ▶ Do not touch overheated prosthetic components.
- ▶ Reduce all activities if there are changes in functionality, in order to allow the overheated prosthetic components to cool down.
- ▶ If overheating or changes in functionality occur, have the prosthesis checked by authorised qualified personnel.

**⚠ CAUTION**

**Use in wet areas**

Fall due to slipping or altered functioning of the prosthetic knee joint

- ▶ Before using the prosthesis in wet areas, lock the prosthetic knee joint (see section "Locking the Prosthetic Knee Joint" – see page 30).
- ▶ Adjust the walking speed to the local conditions.
- ▶ On wet surfaces, only use prosthetic feet with a non-slip sole.

► **Inform the patient.**

**NOTICE**

**Mechanical overload**

Impaired functionality due to mechanical damage

- Check the product for damage prior to each use.
- Do not use the product if its functionality has been impaired.
- Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

**INFORMATION**

**Leak in the product's hydraulic system**

Environmental damage or skin irritation due to leaking hydraulic oil

- Before each use, check the product for leaking hydraulic oil.
- In case of leaks, do not continue using the product and have it repaired promptly.

**6.1.1 Locking the Prosthetic Knee Joint**

**Locking the Prosthetic Knee Joint**

- 1) Bring the prosthetic knee joint to full extension under load.
- 2) Push down the lock lever (see fig. 21).
- 3) Check to verify that the prosthetic knee joint is locked.

**Unlocking the Prosthetic Knee Joint**

- 1) Bring the prosthetic knee joint to full extension under load.
- 2) Pull up the lock lever (see fig. 22).
- 3) Check to verify that the prosthetic knee joint is unlocked.

**6.2 Cleaning**

**⚠ CAUTION**

**Use of unsuitable cleaning agents or disinfectants**

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- Only clean the product with the approved cleaning agents.
- Only disinfect the product with the approved disinfectants.
- Observe the instructions for cleaning and care.

**INFORMATION**

- Clean the product after contact with liquids or if it gets dirty.
- For disinfection, only use disinfectants that do not affect the materials of the product. For further information, please consult the manufacturer.
- Observe the cleaning instructions for all prosthetic components.

- 1) Clean the product with a damp cloth and soap.
- 2) Rinse the product with clean, fresh water.
- 3) Dry the product with a soft cloth.
- 4) **NOTICE! Do not store the product while moist.**  
Allow to air dry in order to remove residual moisture.

## 7 Maintenance

### CAUTION

#### Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

► Observe the following maintenance instructions.

- Arrange regular maintenance intervals with the patient depending on the level of use.
- Following an individual period for the patient to get accustomed to the prosthesis, check the settings of the prosthetic joint and adapt them to the patient's requirements again as needed.
- The prosthetic components should be inspected after the first 30 days of use.
- Inspect the entire prosthesis for wear during normal consultations.
- Conduct annual safety inspections.
- As part of the safety inspections, inspect the prosthetic joint for wear and proper functionality. Special attention should be paid to movement resistance, bearings and abnormal noises. Full flexion and extension must be ensured at all times. Readjust as required.

## 8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

## 10 Technical data

Reference Number	3R80	3R80=ST
Weight [g]	1240	1255
System height [mm]	163	179
Proximal system height up to the alignment reference point [mm]	28	44
Distal system height up to the alignment reference point [mm]	135	
Knee flexion angle	150°	
Proximal connection	Pyramid adapter	Threaded connector
Distal connection	Tube clamp (Ø 34 mm)	
Max. body weight [kg]	150	
Mobility grade	3 + 4	

<b>Reference Number</b>	<b>2R57=16-285-WF</b>
<b>Weight [g]</b>	240
<b>Diameter [mm]</b>	34
<b>Min. system height [mm]</b>	77
<b>Max. system height [mm]</b>	317
<b>Distal connection</b>	Pyramid receiver
<b>Max. body weight [kg]</b>	150

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