ottobock.

8K22, 8K23

DE	Gebrauchsanweisung	3
EN	Instructions for use	8
FR	Instructions d'utilisation	13
IT	Istruzioni per l'uso	18
ES	Instrucciones de uso	23
PT	Manual de utilização.....	29
NL	Gebruiksaanwijzing.....	34
SV	Bruksanvisning	40
DA	Bruksanvisning.....	45
NO	Bruksanvisning	50
PL	Instrukcja użytkowania	55
HU	Használati utasítás	60
CS	Návod k použití	65
RO	Instrucțiuni de utilizare.....	70
TR	Kullanma talimatı.....	75
EL	Οδηγίες χρήσης	80
RU	Руководство по применению	85
JA	取扱説明書.....	91
ZH	使用说明书.....	96

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Explanation of symbols

CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting/use.

INFORMATION

Last update: 2020-10-01

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

The scope of delivery is shown on the front page.

1 Single components

1.1 Single components

For available single components refer catalogue.

1.2 Accessories

Determine glove size according to the size of the hand.

8S4=*	Prosthetic Glove
8S5=*	Prosthetic Glove
8S6=*	Prosthetic Glove
8S4N=*	MovoSkin Natural
8S5N=*	MovoSkin Natural
8S6N=*	MovoSkin Natural
640F12	Special Cleaner
640F13	Pump Sprayer for 640F12

Accessories for long below elbow and through-wrist residual limbs

9S96=40	Chassis	(for hand size L/R6 ¾)
9S184=*	Chassis	(for hand sizes L/R7 ¼, L/R7 ¾ and L/R8)
506G4=M3x5	Set Screw	(for hand size L/R6 ¾)
506G4=M4x5	Set Screw	(for hand sizes L/R7 ¼, L/R7 ¾ and L/R8)
9E94=*	Lamination Ring	(for hand size L/R6 ¾, L/R7 ¼, L/R7 ¾ and L/R8)

The accessories are not included. They must be ordered separately.

2 Description

2.1 Indications for use

The Ottobock System Hands – voluntary opening – are suited solely for the exoprosthetic fitting of the upper extremity in combination with Ottobock arm components.

2.2 Field of application

The System Hands can be used for any residual limb length. The System Hands – voluntary opening – are opened by a harness control system.

2.3 Safety instructions

Please forward the following safety instructions to your patients:

CAUTION

The Ottobock prosthetic components described in this manual have been developed for everyday life activities and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety!

Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the prosthesis inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Service.

CAUTION

An upper extremity amputee's ability to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities.

All persons are required to observe their country's national and state driving laws when operating motor vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test center.

For maximum safety and convenience, Ottobock recommends that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork, automatic shift). Risk-free driving should be ensured even when the prosthesis is not functioning.

NOTICE

Please avoid exposing prosthetic components to surroundings that corrode metal parts, for example, freshwater, saltwater, acids and other liquids.

Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

NOTICE

The prosthesis and its components should not be subjected to intense smoke, dust, vibrations, shocks or high temperatures. Do not allow debris or liquids to get into the prosthesis and its components.

Failure to follow these instructions can lead to malfunctioning and damage of the prosthesis.

NOTICE

Do not use any aggressive cleaning agents since they could cause damage to bearings, seals and plastic parts.

2.4 Function

The System Hands have an outer pull. The connection to the body harness is dorsal.

The Ottobock System Hands are available in the following designs:

8K22=L/R* with metric thread M12x1.5 with outer pull

8K23=L/R* with thread measured in inches ½"-20 with outer pull

For various sizes see table in section 3.

The System Hands – voluntary opening – belong to the cable-activated prostheses. They are opened by pulling the cable (active) and close under spring tension. In closed position the hands lock automatically.

The System Hands carry out the three-finger grip with index and middle fingers as well as thumb.

To ensure a good pattern of opening or closing movements of the whole hand, the ring and little fingers are connected with the middle finger through the inner hand.

The grip force can be changed by adjusting the coil spring. In this case the System Hand should be sent to Ottobock Service.

3 Technical data

Size	6 ¾	7 ¼	7 ¾	8
Weight [g/oz]	215/7.6	300/10.6	330/11.6	340/12
Total height [mm/in]	125/4.92	152/5,98	155/6.1	157/6.18
Thum length [mm/in]	105/4.13	132/5.19	135/5.31	137/5.39

4 Handling

4.1 Assembly

The System Hands are to be mounted to an Ottobock wrist unit.

For long below elbow and through-wrist residual limbs, the chassis with threaded stud should be replaced by the chassis with 9E94=* Lamination Ring as follows:

9S96=40 for hand size L/R 6 ¾

9S184=* for hand sizes L/R 7 ¼, L/R 7 ¾ and L/R 8

The components are attached with the 506G4=* Set Screw (Fig. 1).

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To avoid damaging, use the 711M1 Hand Holding Tool for holding the System Hand while pulling on the glove. The 711M2 Mounting Plate is used for attaching the Hand Holding Tool to a workbench (Fig. 2).

4.2 Combination possibilities

The Ottobock System Hands – voluntary opening – can be combined with the adapters suitable for modular arm prostheses and with the Ottobock wrist units.

If necessary, you must select a suitable adapter from the wrist unit accessories.

5 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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