








ottobock.







12K19, 12K5, 12K20

| | | |
|----|----------------------------------|-----|
| DE | Gebrauchsanweisung | 6 |
| EN | Instructions for use | 13 |
| FR | Instructions d'utilisation | 21 |
| IT | Istruzioni per l'uso | 29 |
| ES | Instrucciones de uso | 36 |
| PT | Manual de utilização | 43 |
| NL | Gebruiksaanwijzing | 52 |
| SV | Bruksanvisning | 59 |
| DA | Bruksanvisning | 66 |
| NO | Bruksanvisning | 73 |
| PL | Instrukcja użytkowania | 80 |
| HU | Használati utasítás | 87 |
| CS | Návod k použití | 94 |
| RO | Instrucțiuni de utilizare | 101 |
| TR | Kullanma talimatı | 109 |
| EL | Οδηγίες χρήσης | 116 |
| RU | Руководство по применению | 124 |
| JA | 取扱説明書 | 132 |
| ZH | 使用说明书 | 138 |

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Explanation of symbols

English

CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION

Last update: 2022-09-06

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

The scope of delivery is shown on the front page.

1 Single components

1.1 Single components

For available single components refer to the catalogue.

2 Description

2.1 Indications for use

The 12K19, 12K5 and 12K20 Ottobock Elbow Set-Ups are to be used **solely** for the exoprosthetic fitting of the upper extremities in combination with Ottobock arm components.

2.2 Field of application

The 12K19, 12K5 and 12K20 Ottobock Elbow Set-Ups can be used for all cosmetic, body-powered and myoelectric arm prostheses.

2.3 Lifetime

Expected lifetime: 5 years

The design, manufacturing and requirements for the intended use of the product are based on the expected lifetime.

2.4 Safety instructions

Please make your patients aware of the following safety instructions:

CAUTION

Mechanical overload

Injuries due to breakage of load-bearing components, impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

The Ottobock prosthetic components described in these instructions for use have been developed for activities of daily living and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.). Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety!

Should prosthetic components be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the components inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Service.

2.5 Environmental conditions

CAUTION

Contact with cold, heat, embers or fire

Injuries (e.g. burns) and damage to the product due to unallowable environmental conditions

- ▶ Keep the product away from open flames, embers and other sources of heat.
- ▶ Do not expose the product to unallowable environmental conditions (see the section „Environmental conditions“).

NOTICE

Contact with salt water or water containing chlorine as well as abrasive substances (e. g. sand)

Damage and premature product wear

- ▶ Following any contact with the substances identified above, promptly clean the product.

Prosthetic components must not be exposed to corrosive elements such as fresh water, salt water and acids.

Using this medical device in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

Prosthetic components must not be subjected to intense smoke, dust, vibrations, shocks or high temperatures. Do not allow debris or liquids to get into the prosthesis and its components. Failure to follow these instructions can lead to malfunctioning and damage to the prosthesis.

2.6 Function

The Ottobock Elbow Set-Ups are available in the following versions:

- 12K19=40 Elbow Set-Up for Children
Colour: caucasian skin colour
- 12K5=* Elbow Set-Up with lamination ring
Colour: caucasian skin colour

- **12K20=* Elbow Set-Up with lamination ring**
Colour: dark brown
Similar to 12K5=*

For available sizes, please refer to section 3.

The elbow set-ups have a passive elbow lock. The different locking positions can be set with a lock button (Fig. 1, item A). The 12K19=40 set-up has 10 possible locking positions in 8° increments and 12K5=* and 12K20=* set-ups have 13 possible locking positions in 8° increments.

The lock button (Fig. 1, item A) is used to manually release and lock the elbow lock.

12K19=40

The friction (restraint) setting of the internal/external rotational movement (sickle movement) in the elbow set-up can be changed by tightening the adjustment screws (Fig. 2, item A).

12K5=* / 12K20=*

The friction (restraint) setting of the rotational movement (sickle movement) in the elbow set-up can be changed by tightening the 13G11/ 13G39 Set Screw (Fig. 3).

For myoelectric arm prostheses, the optional 9E185=50-1 extension cable can be used and routed through the elbow joint.

3 Technical data

| Article number | 12K19=40 | 12K5=45 12K20=45 | 12K5=50 12K20=50 |
|--------------------------------|----------|---------------------|---------------------|
| Weight (g) | 260 | 370 | 330 |
| Upper arm connection Ø (mm) | 54 | 67 | |
| Wrist connection Ø (mm) | 40 | 45 | 50 |

4 Handling

NOTICE

Incorrect alignment, assembly or adjustment

Product damage due to incorrectly installed or adjusted prosthetic components

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

NOTICE

Mechanical damage to the product

Change in or loss of functionality due to damage

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

4.1 Preparing for assembly

For assembling the elbow set-up you need a prefabricated inner socket and a shaped foam part for the outer socket. Inner socket and shaped foam part need to be connected to each other.

4.2 Assembly

- In order to determine the position and length of the elbow set-up, at first temporarily attach all components (including the wrist unit and terminal device) to the shaped foam part.
- In case of unilateral fittings, the distal end of the thumb (Fig. 4) and the centre of rotation of the elbow joint (Fig. 5) are the determining factors for the length of the components. In case of bilateral fittings, priority is given to the prostheses' functionality. When adapting the length, make sure that the prosthesis will offer the patient a high degree of freedom of movement.

- Minor length adaptations are possible by shortening the shaped foam part and the forearm component. It is important, however, that the diameter for the wrist unit remains unchanged in the forearm socket.
- Mark the position of the lamination ring on the shaped foam part (Fig. 6). After determination of the prosthesis length remove the elbow set-up from the shaped foam part.
- Only use Orthocryl sealing resin for attaching the lamination ring, hand adapter etc.; pour from inside into the grooves avoiding the formation of bubbles.

4.2.1 12K5=* / 12K20=*

Disassembling the lamination ring

- To remove the 13G8=67 Lamination Ring (Fig. 1, item B) from the elbow set-up, you first have to separate the 13Z2=4 Lock Ring (Fig. 7, item A) from the 13G11 / 13G39 Set Screw. After loosening the set screw (Fig. 3) you can remove the lock ring. Then remove the lamination ring with attached parts from the elbow ball (Fig. 8).

Assembly

- Isolate the inner contours of the lamination ring with silicone grease and transfer it to the shaped foam part with light putty. This will ensure exact positioning of the lamination ring (Fig. 9). Pull a PVA bag over the shaped foam part and inner contours of the lamination ring and tie off distally. Isolate the interior of the lamination ring with silicone grease and degrease the exterior. Attach the lamination ring over the PVA bag to the shaped foam part (Fig. 10). Laminate the outer socket over the shaped foam part and lamination ring.
- Once the laminate has cured, lay bare the lamination ring by grinding (Fig. 11 – without damaging the anodised layer on the lamination ring) and remove the shaped foam part.
- Mount the elbow ball on the lamination ring again.
- Once again, temporarily affix the wrist unit and terminal device to the elbow set-up to verify the length of the prosthesis. Shorten the forearm socket a bit more, if necessary. It is important, however, that the diameter for the wrist unit remains unchanged in the forearm socket.

Optional (12K5=*/12K20=*): installing the 9E185=50-1 extension cable

For better representation of the installation, the elbow component is shown without the prosthetic socket in the corresponding illustrations.

- Disassemble the lamination ring with upper arm rotator as illustrated and remove the blind plugs (Fig. 12).
- After the prosthetic socket has been laminated, pull the extension cable through the elbow joint and connect it to the 13E190=150 distributor (Fig. 13).
- Once the extension cable is correctly positioned in the upper arm rotator, seal the proximal opening with silicone bonding agent (e.g. 617H46) (Fig. 14, item A). Mount the outer socket on the elbow component and position the compensating plug in the distal opening of the elbow joint (Fig. 14).

To avoid damaging the cable, the elbow joint must not be rotated more than 360° after installing the extension cable.

4.2.2 12K19=40

Disassembling the lamination ring (12K19=40)

- After loosening the 506G1=M8x15 Adjustment Screws (Fig. 2, item A), the lamination ring with 13G9=47 Upper Arm Rotation Piece (Fig. 2, item B) can be separated from the elbow ball.

Assembly

Refer to Assembly of 12K5=* / 12K20=* (Section 4.2.1).

4.2.3 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Do not use any aggressive cleaning agents since they could cause damage to bearings, seals and plastic parts. The product is to be cleaned with a damp cloth and mild soap only.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

9E185=50-1 extension cable

The product meets the requirements under the RoHS Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

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