


ottobock.

16X12, 16X13, 16X14, 16Y27, 16Y31=N

DE	Gebrauchsanweisung	4
EN	Instructions for use	12
FR	Instructions d'utilisation	20
IT	Istruzioni per l'uso	31
ES	Instrucciones de uso	42
PT	Manual de utilização	53
NL	Gebruiksaanwijzing	63
SV	Bruksanvisning	74
DA	Brugsanvisning	84
NO	Bruksanvisning	94
PL	Instrukcja użytkowania	104
HU	Használati utasítás	114
CS	Návod k použití	124
RO	Instrucțiuni de utilizare	134
TR	Kullanma talimatı	145
EL	Οδηγίες χρήσης	155
RU	Руководство по применению	165
JA	取扱説明書	176
ZH	使用说明书	186

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Explanation of symbols

⚠ WARNING Warnings regarding possible risks of severe accident or injury.

⚠ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting/use.

INFORMATION

Last update: 2015-04-10

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

The scope of delivery is shown on the front page.

1 Single components

1.1 Single components

For available single components refer to the catalogue.

1.2 Accessories

16Y26	Pull Cable	for 16Y27
21Y79	Webbing Plate	for 16Y27
16Y30	Modular Connection Piece	for 16Y27
16Y32	Plastic Cover	for 16Y27

The accessories are not included. They must be ordered separately.

2 Description

2.1 Indications for use

The 16X12 elbow joint bars, 16X13 elbow joint bar with cable lock, 16X14 elbow joint bar without lock as well as the 16Y27 joint piece with cable

lock and 16Y31=N joint piece without lock are to be used **solely** for the exoprosthetic fitting of the upper limbs.

2.2 Field of application

Not all patients are able to fully control the cable controlled arm prosthesis with these components. This should be examined in each individual case. However, each patient can use a passive prosthesis.

The 16X12 Elbow Joint Bars are used for very long above-elbow residual limbs or an elbow disarticulation when no fitting with an above-elbow component cannot be performed.

In case of a unilateral cosmetic fitting, the 16X14 Elbow Joint Bar without Lock can be used, too.

2.3 Safety instructions

Please forward the following safety instructions to your patients:

⚠ WARNING

Risk of accident when driving a motor vehicle. An upper extremity amputee's ability to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities.

All persons are required to observe their country's national and state driving laws when operating motor vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test center.

For maximum safety and convenience, Ottobock recommends that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork, automatic shift). Risk-free driving should be ensured even when the prosthesis is not functioning.

⚠ CAUTION

Risk of injury as a result of excessive strain. The Ottobock prosthetic components described here have been developed for activities of daily living and must not be used for unusual activities such as extreme sports (free climbing, paragliding etc.).

Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety!

Should prosthetic components be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the components inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Service.

NOTICE

Risk of corrosion. Prosthetic components must not be exposed to corrosive elements such as fresh water, salt water and acids.

Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

NOTICE

Damage due to incorrect environmental conditions. Prosthetic components must not be subjected to intense smoke, dust, vibrations, shocks or high temperatures. Do not allow debris or liquids to get into the prosthesis and its components. Failure to follow these instructions can lead to malfunctioning and damage of the prosthesis.

NOTICE

Damage caused by improper cleaning agents. Do not use any aggressive cleaning agents since they could cause damage to bearings, seals and plastic parts. The product may only be cleaned with a damp piece of cloth and neutral soap.

2.4 Function

16X12

The 16X12 Elbow Joint Bars consist of the 16X13 Elbow Joint Bar with Cable Lock and 16X14 Elbow Joint Bar without Lock.

The joint bars are made of stainless steel of a flat profile and can be used on both sides.

16X13

The 16X13 Elbow Joint Bar is lockable. It is equipped with the 16Y27 Joint Piece with Cable Lock. The lock is operated via cable.

16X14

The 16X14 Elbow Joint Bar is freely movable. It is equipped with the 16Y31=N Joint Piece without Lock.

3 Technical Data

Article number	16X12	16X13	16X14
Weight [g]	350	185	165

The 16Y27 Joint Piece with Cable Lock has 18 locking positions in 7° increments.

4 Handling

4.1 Installation preparation

In order to determine the position and length of the 16X12 elbow joint bars, temporarily attach all components (including the adapter and hand) to the shaped foam part. To ensure easy access to the cable of the 16Y27 joint piece with cable lock, the joint piece should be adjusted on the lateral side (Fig. 1).

For a unilateral fitting, the distal thumb end and the pivot point of the elbow joint are crucial for determining the length of the components.

For a bilateral fitting, the focus is on the functionality of the prosthesis. When adapting the length, make sure the prosthesis allows the patient a high degree of freedom of movement.

16X12

A pre-fabricated inner socket and above-elbow harness system must be prepared as a prerequisite for the assembly with the 16X12 Elbow Joint Bars. Assemble both components (fig. 3). You will also need a foam part for fabrication of the outside socket.

Before assembling the 16X12 Elbow Joint Bars, check whether the 16Y27 Joint Piece with Cable Lock must be installed on the medial or lateral side. This should be decided according to the movement possibilities of the pa-

tient. Lateral positioning of the joint piece with lock allows easier control with the body harness and is recommended.

It may be necessary to change the position of the joints through positioning the joint plate towards the forearm (Fig. 4, arrow):

- For this purpose, first loosen both freely accessible screws and remove the 16Y29 Bar from the guide.
- Remove the plastic cover (fig. 5/6).
- Turn the dismantled bar by 180° along its longitudinal center line and reassemble the bar (fig. 7).
- Tighten the screws alternately until the bar is guided in the joint piece without play. If the torque is too strong, friction will occur in the joint. This would cause higher wear and tear.
- Subsequently, loosen the second 16Y29 Bar (fig. 7), turn it by 180° – as described above – and reassemble.

The position of the other joint bar should be changed in the same manner.

INFORMATION

The correct torque of the screws is **2.5 Nm to 3 Nm** (fig. 7/8). Use 710D6 Torque Wrench to tighten the screws.

NOTICE

Malfunction of the joint. Do not use threadlockers for securing screw. Failure to follow this instruction can lead to malfunction of the joint.

4.2 Assembly

16X12

- Temporarily attach all components (incl. wrist unit and hand) to the inner socket for positioning and determining the length of the elbow joint bars. Make sure the axis of rotation of the 16Y27 Joint Piece with Cable Lock and 16Y31=N Joint Piece without Lock are in one line. This minimizes the wear of joints through transversal forces which would occur.
- Mount the elbow joint bars so that the joint plate is positioned towards the forearm. More space is available in the forearm for accommodation of the joint plate, which facilitates subsequent shaping of the bars.

- Mark the necessary length on the elbow joint bars both in the above-arm and forearm areas and shorten the bars. Adjust the length of the bars to the physiological conditions as necessary. The longer the bar the better the force transmission. The minimum length of the bars from the axis of rotation must not be less than 10 cm. Carefully debur the cutting edges.
- To avoid the excessive increase of the forearm, they are counter-sunk into the foam part. Transfer the contours for this purpose (fig. 14). Remove the material from the contour corresponding to the bar thickness.
- Adjust the bars in the above-elbow area to the physiological conditions by bending using the 711S4 Bending Iron. To prevent the damage of the mechanism of the 16Y27 Joint Piece with Cable Lock, do not start bending nearer than from the distal area of the first hole (fig. 15).
- Subsequently attach the bars temporarily once again together with the forearm foam part to the inner socket. Verify the length relationships and movement possibilities and carry out necessary adjustments.
- Mark the bar contours on the inner socket (fig. 16). Place a dummy (e.g. 616T3 Thermolyn Trolen), thickness 2 mm, into the marked area.
- Laminate the above-elbow and forearm sockets individually.
- Attach the elbow joint bars. For this purpose, drill at least two holes (\varnothing 4.5 mm) for each bar into the outside above-elbow socket. Transfer the hole position on the bars and tap the M4 threads in the bars. Attach the bars with 503F3 Truss Head Screw (fig. 17).
- Bond the bars for the forearm socket.

 **CAUTION**

Processing the bars. Clamp the bars for drilling holes. Use a sharp drill. The material must not become discolored. Carefully debur the cutting edges.

4.3 Adjustment and final assembly

NOTICE

Risk of damage to the joint piece with cable lock during unlocking/locking. The cable must be integrated into the body harness so it can be pulled upwards as much vertically as possible.

In case of improper cable installation requiring from the patient to pull the cable too much in direction to the center of the torso, there is the risk of damage to the joint piece with cable lock.

4.3.1 Adjustment of the control cable to the harness

16X12

- We recommend to use the 10Y19=1 Coupling Screw and 21A12 Threaded Fitting, short (for version with Perlon cable) or the 10Y25 Clamp Sleeve, short (for version with a steel cable) for attachment of the control cable to the body harness.
- Cut the steel cable of the joint piece with cable lock on a suitable place. We recommend to cut the cable approx. 12 cm behind the rotation axis.
- Remove the isolating tubing and steel cable from the joint (fig. 19). Shorten the isolating tubing according to the steel cable length and pull it on the steel cable again.
- Clamp the coupling screw of the 10Y19=1 Screw Coupling on the steel cable using the 736Y6 Clamping Tool.
- We recommend to re-mount the steel cable of the lock and to don the prosthesis once again for determination of the necessary length of the body harness pull-cable. This allows you to determine the length easily and shorten the body harness pull-cable accordingly.
- Pull the coupling sleeve of the screw coupling over the pull-cable and screw the 21A12 Threaded Fitting, short on the Perlon cable (fig. 20). Verify the function of the locking mechanism by pulling the body harness.
- When using a body harness with steel cable, use a 10Y25 Clamp Sleeve, short instead of the 21A12 Threaded Fitting, short. Clamp the clamp sleeve with the 736Y6 Clamping Tool.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH