




12S4, 12S7

DE	Gebrauchsanweisung	6
EN	Instructions for use	12
FR	Instructions d'utilisation	18
IT	Istruzioni per l'uso	26
ES	Instrucciones de uso	31
PT	Manual de utilização.....	38
NL	Gebruiksaanwijzing	44
SV	Bruksanvisning	51
DA	Brugsanvisning.....	57
NO	Bruksanvisning	63
PL	Instrukcja użytkowania	69
HU	Használati utasítás	75
CS	Návod k použití	81
RO	Instrucțiuni de utilizare.....	91
TR	Kullanma talimatı.....	98
EL	Οδηγίες χρήσης	100
RU	Руководство по применению	107
JA	取扱説明書.....	114
ZH	使用说明书.....	120

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Explanation of symbols

English

⚠ WARNING Warnings regarding possible risks of severe accident or injury.

⚠ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

INFORMATION

Last update: 2016-01-14

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

The scope of delivery is shown on the front page.

1 Single components

1.1 Single components

For available single components refer to the catalogue.

1.2 Accessories

10R2=M12x1.5 Adapter for 12S7 Ottobock Ball Shoulder Joint

The accessories are not included in the delivery. They must be ordered separately.

2 Description

2.1 Indications for use

The 12S4 Ottobock shoulder joint and the 12S7 Ottobock ball shoulder joint are to be used **exclusively** for the exoprosthetic fitting of the upper limbs in combination with Ottobock arm components.

2.2 Field of application

The 12S4 Ottobock Shoulder Joint and the 12S7 Ottobock Ball Shoulder Joint can be used on both sides for shoulder disarticulation and shoulder girdle amputation.

These Ottobock Shoulder Joints can be used as component in an Ottobock modular prosthesis.

2.3 Safety instructions

Please make your patients aware of the following safety instructions:

WARNING

Risk of accident when driving a motor vehicle. An upper extremity amputee's ability to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities.

All persons are required to observe their country's national and state driving laws when operating motor vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test centre.

For maximum safety and convenience, Ottobock recommend that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork, automatic shift). Risk-free driving should be ensured even when the prosthesis is not functioning.

CAUTION

Risk of injury as a result of excessive strain. The Ottobock prosthetic components described in these instructions for use have been developed for activities of daily living and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety.

Should prosthetic components be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the components inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Service.

NOTICE

Risk of corrosion. Prosthetic components must not be exposed to corrosive elements such as fresh water, salt water, and acids.

Using this medical device in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

NOTICE

Risk of damage due to incorrect environmental conditions. Prosthetic components must not be subjected to intense smoke, dust, vibrations, shocks or high temperatures. Do not allow debris or liquids to get into the prosthetic components. Failure to follow these instructions can lead to malfunctioning and damage to the prosthesis.

NOTICE

Damage caused by improper cleaning agents. Do not use any aggressive cleaning agents since they could cause damage to bearings, seals and plastic parts. The product is to be cleaned with a damp cloth and mild soap only.

2.4 Function

The Ottobock Shoulder Joints are available in the following versions:

12S4 Ottobock Shoulder Joint with lamination ring

12S7 Ottobock Ball Shoulder Joint with threaded stud

The 12S4 Ottobock shoulder joint as well as the 12S7 Ottobock ball shoulder joint have two degrees of freedom.

They allow for anteversion/retroversion and abduction/adduction movements. The Ottobock Shoulder Joints are connected to the elbow component with bars.

Connection of an Ottobock modular prosthesis to the 12S7 Ottobock Ball Shoulder Joint is possible with the 10R2 Adapter.

3 Technical data

Article number	12S4	12S7
Weight [g]	135	80

4 Handling

4.1 Preparing for assembly

For assembling the Ottobock Shoulder Joints you need a well fitting inner socket. Attach a shaped part which may be needed for shoulder compensation to the inner socket.

For position and length adjustment, temporarily attach all single components to the inner socket.

In case of unilateral fittings, the distal end of the thumb and the centre of rotation of the elbow joint are the determining factors for the length of the components (Fig. 1).

In case of bilateral fittings, priority is given to the prosthesis' functionality. When adapting the length, make sure that the prosthesis will offer the patient a high degree of freedom of movement.

With the 12S4 Ottobock shoulder joint, length adaptations are possible by shortening the upper arm bars. Make sure, however, that there will be a distance of at least 6 mm between the end of the bar and the fixing hole centre. This will ensure full surface contact of the fixing screw.

With the 12S7 Ottobock Shoulder Joint, only make length adaptations via the 12R4=* Elbow Joint.

Pay particular attention to the free swinging range of the shoulder joint. The posterior stop should be at approximately 10° of retroversion (Fig. 2).

4.2 Assembly

INFORMATION

For patients with subtotal shoulder disarticulation, we recommend to mount the shoulder joint distally at a little bit lower shoulder height for reasons of symmetry and optical reasons.

12S4

- 1) Transfer the position of the lamination ring to the shaped part (Fig. 2).
- 2) Apply 640Z5=* Plaster Insulating Cream to the inner side of the lamination ring and fill with light putty (Fig. 3/4). Using the light putty, transfer the inner contours of the lamination ring to the shaped foam part (Fig. 5)

12S7

- 1) Individually adapt the shoulder support by bending and shortening. To avoid damages to the shoulder support during bending, we recommend using the 711S4 and 711S5 Bending Irons (Fig. 6/7).
- 2) Transfer the position of the shoulder support to the shaped part (Fig. 8).
- 3) Apply 640Z5=* Plaster Insulating Cream to the shoulder support and fill with light putty (Fig. 9/10). Using the light putty, transfer the inner contours of the shoulder support to the shaped foam part (Fig. 11/12).

Fabrication of the outer socket

The fabrication of the outer socket is described in the following at the example of the 12S4 Ottobock Shoulder Joint:

- 1) Prior to lamination, the free swinging range and length adaptation must be checked.
- 2) Disassemble the lamination ring or shoulder support from the Ottobock Shoulder Joint.
- 3) Laminate the outer socket over the shaped part and lamination ring / shaped part and shoulder support. Use 616G12 Carbon Fibre Cloth in the transition area from outer socket to lamination ring / outer socket to

shoulder support in order to achieve secure force transmission between the outer socket and modular arm component (Fig. 13). In the area of the fixing holes, use the 507S15 Steel Washer with hole.

- 4) Once the laminate of the outer socket has cured, laminate the cap for later attachment of the cosmetic foam cover. This will allow easy disassembly of the cosmetic foam cover when changing the friction setting on the Ottobock Shoulder Joint. For this purpose, pull three layers of 623T3=* Perlon Stockinette over the outer socket, which serve as spacer between the outer socket and cap.
- 5) Lay a PVA bag over the Perlon Stockinette for laminating the cap onto. Again, use 616G12 Carbon Fibre Cloth and the 507S15 Steel Washer with hole for better force transmission (Fig. 14).
- 6) Once the laminate has cured, drill the required fixing holes between the outer socket and inner socket. Then remove the inner socket and the shaped foam part. Next, drill the bore holes through the outer socket and cap.
- 7) In the area of the joint connection, remove the lamination ring or shoulder support from the laminate (Fig. 15). Create the outer shape of the outer socket and cap by cutting and grinding (Fig. 16).

Subsequently, assemble the outer socket, cap and shoulder joint.

4.3 Adjustments and final assembly

12S4

The T-screw (Fig. 3; item 1) is used to adjust the friction of anteversion and retroversion movement. Adjust it on the inner side of the outer socket.

2S7

Change the movement friction setting by adjusting the clamp of the ball on the 501Z2=M4x18 Cap Screw (Fig. 12; item 1).

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard

of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH