
ottobock.

2R119

DE	Gebrauchsanweisung	5
EN	Instructions for use	10
FR	Instructions d'utilisation	15
IT	Istruzioni per l'uso	20
ES	Instrucciones de uso	25
PT	Manual de utilização.....	30
NL	Gebruiksaanwijzing.....	35
SV	Bruksanvisning.....	40
DA	Brugsanvisning.....	45
NO	Bruksanvisning.....	50
PL	Instrukcja użytkowania.....	55
HU	Használati utasítás	60
CS	Návod k použití.....	65
RO	Instrucțiuni de utilizare.....	70
TR	Kullanma talimatı.....	75
EL	Οδηγίες χρήσης	80
RU	Руководство по применению	85
JA	取扱説明書.....	90
ZH	使用说明书.....	95

INFORMATION

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INFORMATION

Last update: 2020-10-01

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

Figures 5, 6 and 8 apply to both, test sockets and definitive sockets.

Explanation of symbols

English

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting/use.

1 Description

1.1 Intended Use

The Ottobock Vacuum Connector is to be used **exclusively** for prosthetic fittings with Harmony. The vacuum connector is used for the connection between the socket and Harmony pump.

1.2 Safety instructions

Please forward the following safety instructions to your patients:

NOTICE

Risk of corrosion. Please avoid exposing prosthetic components to surroundings that corrode metal parts, for example, freshwater, saltwater, acids and other liquids.

Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

NOTICE

Risk of damage if using powder. Any kind of powder on the inner side of the socket or on the skin can restrict or disable the functionality of the vacuum connector during use. Avoid using powder in the prosthetic socket or on the liner.

1.3 Function

The vacuum connector has no thread (click connection). It is opened and closed by pressing together two levers. It serves to conduct air from the socket to the Harmony pump through a hose.

2 Components

2R119 Vacuum Connector (Fig. 1)
(1) Lower screw ring
(2) Upper screw ring
(3) Upper connector part
(4) Hose connector
(5) Vacuum forming dummy
(6) Lamination dummy
(7) 501S86=M5x12 Countersunk Head Screw
(8) 503S3 Screw
(9) Pin wrench
(10) Hose

3 Handling

The dummy provided is used during vacuum forming or lamination and allows for precise fitting into the socket made of different materials. Please always make sure that the full dummy surface rests on the plaster positive.

3.1 Setup

NOTICE

Risk of damage to the vacuum connector. The vacuum connector is not intended for the fixation of inner socket and frame socket. It has no holding function. Attach the vacuum connector only to the flexible inner socket.

3.1.1 Installation Preparation for Sockets Made of ThermoLyn

The enclosed vacuum forming dummy (5) provides for that the lower screw ring is flush on the inner side of the socket.

- Soak 99B25 Nylon Stockinette with 85F1 Plaster Parting Agent and pull it over the plaster positive.
- Secure the vacuum forming dummy to the plaster positive with the enclosed 503S3 Screw (Fig. 2).

3.1.2 Completion of a Test Socket

- After vacuum forming of the test socket of ThermoLyn, grind away the material covering the vacuum forming dummy until it is flush (Fig. 3).
- Remove the 503S3 Screw and vacuum forming dummy.
- Remove the vacuum connector and apply 617H46 Bonding Agent to the contact surface of the lower screw ring if need be (Fig. 4). Position the lower screw ring in the socket (Fig. 5). Tighten the upper screw ring. The pin wrench can be used for tightening the lower and upper screw rings (Fig. 6).

3.1.3 Completion of a Flexible Inner Socket with Carbon Frame

The enclosed lamination dummy (6) provides for the space required for the upper screw ring in the laminated socket.

- After vacuum forming of the definitive inner socket of ThermoLyn, grind away the material covering the vacuum forming dummy until it is flush (Fig. 3).
- Remove the 503S3 Screw.
- Put the lamination dummy onto the vacuum forming dummy and secure with the enclosed countersunk head screw (Fig. 7). The 5Z1 Lamination Kit with the aid of the Technical Information 2.3.5 or the 5Z4 Lamination Kit for CAT-CAM can be used for reinforcement of the frame.
- After laminating, grind away the material covering the lamination dummy until it is flush and remove the countersunk head screw.

- Remove the carbon frame with inner socket from the plaster positive.
- Remove the vacuum connector and apply 617H46 Bonding Agent to the contact surface of the lower screw ring if need be (Fig. 4). Position the lower screw ring in the socket (Fig. 5). Tighten the upper screw ring. The pin wrench can be used for tightening the lower and upper screw rings (Fig. 6).
- Position the flexible inner socket in the carbon frame. Make sure that it is easy to install and remove the upper connector part into/from the base. Remove more excessive material if necessary.

3.1.4 Installation Preparation for Laminated Sockets

The enclosed vacuum forming dummy provides for that the lower screw ring is flush with the inner side of the socket.

- Pull 99B25 Nylon Stockinette over the plaster positive.
- Secure the vacuum forming dummy to the plaster positive with the enclosed 503S3 Screw (Fig. 2).
- Pull 99B81 PVA Bag over the plaster positive and reinforce the laminated socket adequately.

3.1.5 Completion of a Laminated Socket without Flexible Inner Socket

- After laminating the lamination resin socket, grind away the material covering the vacuum forming dummy until it is flush.
- Remove the 503S3 Screw and vacuum forming dummy.
- Remove the vacuum connector and apply 617H46 Bonding Agent to the contact surface of the lower screw ring if need be (Fig. 4). Position the lower screw ring in the socket (Fig. 5). Tighten the upper screw ring. The pin wrench can be used for tightening the lower and upper screw rings (Fig. 6).
- The lamination dummy is not needed here.

3.1.6 Final assembly

Using the enclosed hose (10) establish the connection between the upper connector part and the Harmony pump.

To do so, attach one end of the hose to the hose connector (4) and attach the other end either to the one-way valve of the mechanical Harmony pump or to the filter of the electric pump (Fig. 8).

3.2 Instructions for use

Position the upper connector part in the base by simply pressing it in place. An audible feedback is emitted once the upper connector part has been pressed in properly. To remove the upper connector part, lift the two wings with your index finger and thumb and pull the upper connector part out of the base.

3.3 Maintenance instructions

The vacuum connector should be cleaned regularly with clear, warm water.

4 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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