









9E420=L, 9E420=R

DE Gebrauchsanweisung	5
EN Instructions for use	18
FR Instructions d'utilisation	30
IT Istruzioni per l'uso	43
ES Instrucciones de uso	56
NL Gebruiksaanwijzing	69
SV Bruksanvisning	82
JA 取扱説明書	95

INFORMATION

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INFORMATION

Date of the last update: 2016-09-06

- Please read this document carefully.
- Follow the safety instructions and the precautions specified in this accompanying document.
- If assistance is required during the start-up, use or maintenance of the product, as well as in the event of unexpected operating behaviour or circumstances, contact the manufacturer (see manufacturer's address on the back of this accompanying document).

Explanation of Symbols

English

 **WARNING** Warnings regarding possible risks of severe accident or injury.

 **CAUTION** Warnings regarding possible risks of accident or injury.

 **NOTICE** Warnings regarding possible technical damage.

 **INFORMATION** Additional information on the fitting / use.

1 Application

Only put the product into use in accordance with the information contained in the accompanying documents supplied. Instruct the patient in the proper use and care of the product. The product may not be transferred to the patient without prior instruction.

1.1 Indications of use

The 7in1 Controller 9E420 is intended **exclusively** for exoprosthetic fittings of the upper limbs and performs control and system-internal coordination of all active prosthetic components.

1.2 Use / Application

The 7in1 Controller 9E420 can be used for unilateral or bilateral amputees or in case of dysmelia. The 9E420 7in1 Controller may only be used for the Ottobock Electric Hand 2000 and the VASI Hand VV* in combination with the 757B25=1 and 757B35=* MyoEnergy Integral.

1.3 Conditions of Use

The 7in1 Controller 9E420 has been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

1.4 User / Qualifications

Fitting a user with the controller may only be performed by prosthetists who have been trained and certified by Ottobock in the use of the prosthetic components of the prosthesis system.

2 Safety Instructions

WARNING

Non-observance of safety instructions. Failure to follow the safety instructions given below can lead to damage to or malfunction of the product. Follow the safety instructions and the precautions specified in this accompanying document.

⚠ WARNING

Fitting of children. When fitting children, there is the risk that they might swallow small system components. This can result in a risk of injury.

Never leave children unattended during the fitting.

⚠ CAUTION

Adjusting the electrodes. If the electrode settings are excessively sensitive, powerful electromagnetic radiation (e.g. in the vicinity of high-voltage power lines, transmitters, transformers, department store security systems) may cause malfunctions of the 7in1 Controller 9E420 and the connected components. This can result in a risk of injury to the patient.

Ensure that the electrode settings are at the lowest possible sensitivity.

⚠ CAUTION

Penetration of liquids. The penetration of liquids may lead to malfunction or defects of the 7in1 Controller 9E420 and result in a risk of injury to the patient.

The 7in1 Controller 9E420 may not be used for bathing prostheses.

⚠ CAUTION

Installation of cables. Improper installation may cause damage to the cables. This can lead to faulty control or malfunction of the 7in1 Controller 9E420 and result in a risk of injury to the patient.

When installing the cables, make sure there are no tight angles or kinks.

⚠ CAUTION

Use of unsuitable prosthetic components. Using unsuitable prosthetic components may lead to faulty control or malfunction of the 7in1 Controller 9E420 and result in a risk of injury to the patient.

Only combine the 7in1 Controller 9E420 with prosthetic components approved by Otto Bock (see Instructions for Use and catalogues).

⚠ CAUTION

Incorrect assembly. Incorrect assembly of the prosthetic components may lead to faulty control or malfunction of the 7in1 Controller 9E420 and result in a risk of injury to the patient.

For proper assembly, pay special attention to Section 9 „Connecting the 7in1 Controller“ in these Instructions for Use.

⚠ CAUTION

Improper use. Any type of improper use can lead to faulty control or malfunction of the 7in1 Controller 9E420 and result in a risk of injury to the patient.

Use Section 3 „Patient Information“ to inform your patient about the proper use of the 7in1 Controller 9E420.

NOTICE

Improper care of the housing

- Damage to the housing through the use of acetone, benzene or similar solvents.
- Only use a damp cloth and mild soap to clean the housing.

INFORMATION

The 7in1 Controller 9E420 can be individually adapted to the needs of the patient. Please be aware that only a carefully adjusted prosthesis will improve the patient's quality of life.

3 Patient Information

WARNING

Operating the product near active implantable systems

- Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the product.
- When operating the product in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

CAUTION

Penetration of dirt and humidity. The penetration of dirt and humidity may lead to faulty control or malfunction of the 7in1 Controller 9E420 and result in a risk of injury to the patient.

Ensure that neither solid particles nor liquids can penetrate inside the 7in1 Controller 9E420.

CAUTION

Mechanical overloading. External mechanical influences or loads, such as impacts and vibration, can lead to faulty control or malfunction of the 7in1 Controller 9E420 and result in a risk of injury to the patient.

The 7in1 Controller 9E420 should not be subjected to mechanical vibrations or impacts.

CAUTION

Thermal overloading. Extreme temperature conditions can lead to faulty control or malfunction of the 7in1 Controller 9E420 and result in a risk of injury to the patient.

Avoid exposure to temperatures outside the allowable range (see Section 12 „Technical Data“).

CAUTION

Penetration of dirt and humidity. The penetration of dirt and humidity may lead to faulty control or malfunction of the 7in1 Controller 9E420 and result in a risk of injury to the patient.

Ensure that neither solid particles nor liquids can penetrate inside the 7in1 Controller 9E420.

⚠ CAUTION

Distance to communication devices. If the distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices) is too small, interference with internal data communications can result in unexpected

prosthesis behaviour. Therefore, keeping the following minimum distances to these HF communication devices is recommended:

- Mobile phone GSM 850/GSM 900: 0.99 m
- Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
- DECT cordless phones incl. base station: 0.35 m
- WiFi (routers, access points,...): 0.22 m
- Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

INFORMATION

This product may not be disposed of with regular domestic waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

4 Adjusting the MyoBock electrodes

INFORMATION

- Allow the patient to rest during the adjustment of the electrodes, otherwise muscle fatigue will lead to inconsistent results and the therapist will, as a result, tend to adjust the electrodes with too high a sensitivity.
- The electrodes should be placed on intact skin only and with as much electrode-skin contact as possible.
- In the case of strong interference from electric devices, the position of the electrodes should be verified and adjusted if necessary. If interference cannot be eliminated, please contact Otto Bock Myo-Service.
- Before connecting, apply sufficient Silicone Grease 633F11 to the electrode pin contacts and to the connecting cable contacts. The grease will protect the connection from humidity.

5 Scope of Delivery

1 pc. 7in1 Controller 9E420

1 pc. Contact Plate with cable 9E414

1 pc. Lamination Dummy 9E421

6 Accessories

1 pc. Lamination Ring 10S16

1 pc. Installation Tool 711M50-1

1 pc. Adapter Cable 757VAK2 for VASI Hand

1 pc. MyoEnergy Integral 757B25=1, 757B35=*

1 pc. MyolinoLink 60X6

1 pc. MyolinoSoft 560X3=* Software

7 Description

The 7in1 Controller can be used to control all sizes of the Ottobock Electric Hand 2000. When combined with the Adapter Cable 757VAK2, all hand sizes of the VASI Hand can be controlled as well.

7.1 Control Modes

7 control modes and the option to make patient-specific adjustments are available for optimal adjustment to the patient.

Settings can only be changed via Bluetooth data transfer and the 560X3=* MyolinoSoft software. In order to do so, a wireless connection must be established between the 7in1 Controller and the PC with the 60X6 MyolinoLink, which is connected to the charging receptacle of the 757B25=1, 757B35=* MyoEnergy Integral. Please consult the instructions for use enclosed with the respective products for detailed information.

Upon delivery, the DMC plus program for 2 electrodes is activated. See Section 10.

7.2 Control components

The following control components can be used and combined with one another.

- MyoBock Electrode 13E200=*
- Suction Socket Electrode 13E202=*
- Harness Pull Switch 9X14
- Cable Pull Switch 9X18
- Rocker Switch 9X25
- Pressure Switch 9X37
- Linear Transducer 9X50/9X52

7.3 Power Supply

Power to the 7in1 Controller and the Ottobock Electric Hand 2000 is supplied exclusively by the 757B25=1 or 757B35=* MyoEnergy Integral. Please consult the instructions for use enclosed with the product for detailed information.

7.4 Safety Function

In an emergency, it is possible to open the electric hand without Myo signals from the patient. Read Section 11 of these Instructions for Use.

7.5 Battery Manager

An intelligent „battery manager“ alerts you to low voltage conditions - the hand will become slower and the grip force will decrease. Finally, when there is very little battery capacity remaining, the hand will switch off automatically to protect the battery against harmful total discharge.

8 Recommended Installation Instructions for the 7in1 Controller

The 7in1 Controller is positioned in the forearm socket. The flat, slightly rounded design allows it to be secured in the outer socket or between the outer and inner socket depending on space availability. When little space is available, the 7in1 Controller is laminated into the outer socket. If there is sufficient space between the outer and inner socket, the 7in1 Controller can be secured using a hook-and-loop closure strap.

8.1 Laminating instructions (Fig. 1 through Fig. 15)

If there is insufficient space between the inner and outer socket, then the 7in1 Controller can be laminated into the outer socket.

INFORMATION

Before handling and processing the Pedilen Rigid Foam 200 617H12=* and the Pedilen Hardener 617P21, please read the corresponding material safety data sheet.

Before handling and processing the C-Orthocryl Lamination Resin 617H55=* or the Orthocryl Lamination Resin 617H19=* and the Hardener 617P37=*, please read the corresponding material safety data sheet.

Ensure sufficient ventilation of working rooms and spaces.

For your own protection wear appropriate respirator, protective gloves and goggles when sanding the socket.

INFORMATION

Please note that for the following lamination technique, the 7in1 Controller lamination dummy must not be fixed to the inner socket **with** double-sided tape or the like. It is only held in place by the Perlon Stockinette 99B25. Modelling with Plastilin Clay is therefore not required.

1. Fabricate the inner socket. Cut Perlon Stockinette 99B25 to twice the length of the socket positive. Pull over inner socket and tie off distally. Push the 7in1 Controller lamination dummy under the stockinette and position, but do not fix in place (Fig. 1).
2. Pull the Perlon Stockinette 99B25 down, removing all folds, and tie off proximally. Verify the position of the lamination dummy. Warm the PVC Bag 99B71 and pull it over the stockinette. Tie the bag off and evacuate (Fig. 2).
3. Pour foam positive using Pedilen Hard Foam 200 617H12=* and Pedilen Hardener 617P21 (Fig. 3).
4. Remove foam positive from the inner socket. The contour of the lamination dummy can be seen on the inside of the foam positive (Fig. 4). Grind the foam positive distally first.
5. Place lamination ring and lamination cover distally and mark contour (Fig. 5). Fasten the 7in1 Controller lamination dummy to the foam positive with double-sided tape in order to obtain the contour during subsequent grinding (Fig. 6). Finish grinding foam positive (Fig. 7).
6. Make a hole with a diameter of approx. 15 mm and a depth of 10 mm in the distal area of the foam core using countersink drill or knife (Fig. 8). Apply talcum powder to foam core.
7. Pull PVA Bag 99B81 over the foam core. Tie off the PVA bag proximally and distally. Cut off excess bag in the distal area.
8. Switch vacuum on and place knot into the hole. Place adhesive tape over the knot (Fig. 9).
9. Mount lamination ring and lamination cover distally. Cut Perlon Stockinette 623T3=* to twice the length of the socket positive, pull over the socket positive and tie off distally (Fig. 10).
10. In order to improve the socket's stability, it is reinforced with pieces of Carbon Fibre Cloth 616G12 placed around and along the length of the socket – these must be measured and cut to fit the shape of the socket. To do so, the contour of the carbon fibre cloth pieces is determined with a paper template and then transferred to the carbon fibre cloth.

INFORMATION

Note that the lamination ring and the 7in1 Controller must not be covered with the carbon fibre cloth for reasons of electrical safety (Fig. 11).

11. Tie the Perlon Stockinette off in the lowermost groove of the lamination ring with a string (Fig. 12). Pull Perlon Stockinette down without wrinkles and tie off proximally (Fig. 13).
12. Cut second layer of Perlon Stockinette 623T3=* to twice the length of the socket positive and pull it over the socket positive. Tie Perlon Stockinette off in the uppermost groove of the lamination ring using a string. Pull Perlon Stockinette down without wrinkles and tie off proximally.
13. Pull PVA Bag 99B81 down over the positive model without wrinkles and tie off proximally. Apply vacuum. Laminate the outer socket using known methods and allow it to harden.
14. Remove lamination cover (Fig. 14) and finish socket as usual (Fig. 15).

9 Connect 7in1 Controller (Fig. 16 to Fig. 20)

To end up with a functional 7in1 Controller, the following sequence must be observed:

1. Turn the power supply off by pressing the MyoEnergy Integral 757B25=1 or 757B35=* charging receptacle (Fig. 21, arrow).
2. Lubricate the contact bushings of the cables and the connector plugs of the 7in1 Controller with Silicone Grease 633F11.
3. Plug MyoEnergy Integral 757B25=1 or 757B35=* (2-pole and 3-pole connector), electrode (at least 1 on input 1) and the 3-pole connector of the contact plate into the 7in1 Controller all the way (Fig. 16).
4. Fix the contact plate to the lamination ring. To do so, lay the contact plate on the installation tool – the integrated magnet holds the contact plate in place (Fig. 17). Insert the installation tool in the socket from the proximal direction (Fig. 18) and allow the contact plate to engage in the lamination ring O-ring.

INFORMATION

For demonstration purposes, the socket has been cut open.

5. If the contact plate must be pressed out, then this is also done using the installation tool (Fig. 19).
6. When using the VASI Hand, plug the Adapter Cable 757VAK2 into the 7in1 Controller and connect to the VASI Hand.
7. Lay the 7in1 Controller into the laminated trough of the outer socket (Fig. 20) or fix it in the outer socket using a hook-and-loop closure strap.
8. Connect the terminal device and turn the system on again by pressing the charging receptacle (Fig. 21, arrow).

INFORMATION

If the components were connected while the MyoEnergy Integral was turned on, then a system reset (switching the MyoEnergy Integral off and on) must be performed.

10 Control Options

10.1 DMC plus control with 2 electrodes

The grip speed or grip force are determined by the strength of the electrode signal (resulting from muscle tension), making the gripping process significantly more physiological and the grasping of small, fragile objects possible without a problem. The grip force can be increased up to the maximum grip force at any time through a stronger electrode signal.

The gripping strength can be adjusted via the 560X3=* MyolinoSoft software.

Opening: Proportionally controlled by the OPENING electrode.

Closing: Proportionally controlled by the CLOSING electrode.

Program	OPEN	CLOSE	Indications
DMC plus	Sustained electrode signal Speed proportional	Sustained electrode signal Grip force: proportional After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the prosthetic component with unintended electrode signals is prevented. Speed proportional	For patients with 2 strong electrode signals.

10.2 DMC LowInput control with 2 electrodes

The grip speed and grip force are determined by the strength of the electrode signal (resulting from muscle tension). As a result, the gripping process is significantly more physiological and the grasping of small, fragile objects possible without a problem.

In contrast to the DMC plus control, the maximum grip speed and grip force are reached more quickly. This control should only be used for children with low muscle potential.

The gripping strength can be adjusted via the 560X3=* MyolinoSoft software.

Opening: Proportionally controlled by the OPENING electrode.

Closing: Proportionally controlled by the CLOSING electrode.

Program	OPEN	CLOSE	Indications
DMC LowInput	Sustained electrode signal Reduced proportionality range: maximum speed once LOW-threshold has been reached Speed proportional	Sustained electrode signal Grip force: proportional Reduced proportionality range: maximum speed once LOW threshold has been reached Speed proportional	For patients with 2 strong electrode signals.

10.3 AutoControl LowInput control with 2 electrodes

The prosthetic component closes at proportional speed, and the grip force increases proportionally with the duration of gripping.

The gripping strength can be adjusted via the 560X3=* MyolinoSoft software.

Opening: at constant speed by a muscle signal at any level exceeding the ON threshold of the OPENING electrode.

Closing: Proportionally controlled by the CLOSING electrode.

Program	OPEN	CLOSE	Indications
AutoControl LowInput	Sustained electrode signal Speed constant	Sustained electrode signal Grip force: Time proportional Reduced proportionality range: maximum speed once LOW threshold has been reached. Speed proportional	For patients with 2 weak electrode signals.

10.4 Digital control with 2 electrodes

The prosthetic component closes at constant speed. The grip force increases proportionally to the duration of gripping. The OPEN/CLOSE speed and the grip force build-up speed can be adjusted using the MyolinoSoft 560X3=* software.

The gripping strength can be adjusted via the 560X3=* MyolinoSoft software.

Opening: Digitally; signal on the OPENING electrode opens the prosthetic component at constant speed.

Closing: Digitally; signal on the CLOSING electrode closes the prosthetic component at constant speed.

Program	OPEN	CLOSE	Indications
Digital	Sustained electrode signal Speed constant	Sustained electrode signal Grip force: Duration of the signal Speed constant	For patients with 2 weak electrode signals.

10.5 VarioControl using 1 electrode

INFORMATION

For this control mode, an electrode must be connected to input1.

With this control mode, the OPENING speed is determined by the speed and strength of the muscle contraction. The CLOSING speed and grip force depend on the reduction of muscle tension. The gripping strength can be adjusted via the 560X3=* MyolinoSoft software.

- Opening:** proportional. The OPENING speed is determined by the speed and strength of the musclecontraction.
- Closing:** proportional. The CLOSING speed is determined by the speed and strength of the musclerelaxation.
- Stopping:** Through very slow muscle relaxation to the electrode, the electric hand stays in the opened position.

OPEN	CLOSE	Indications
Increasing electrode signal through muscle contraction. Speed and strength of muscle contraction to the electrode. Speed proportional	Declining electrode signal through muscle relaxation. Grip force: Proportional to the decline of the electrode signal. After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the prosthetic components with unintended electrode signals is prevented. Speed proportional	For patients with 1 strong electrode signal or tendency to co-contraction.

10.6 DoubleChannel control with 1 electrode

INFORMATION

For this control mode, an electrode must be connected to input1.

With this control mode, the prosthetic component is opened with a quick and strong signal and closed with a slow, gentle signal.

The gripping strength can be adjusted via the 560X3=* MyolinoSoft software.

Opening: Digitally; fast strong signal on the electrode opens the prosthetic component at constant speed.

Closing: Digitally; slow, gentle signal on the electrode closes the prosthetic component at constant speed.

Program	OPEN	CLOSE	Indications
DoubleChannel	Quick, strong electrode signal that reaches the upper threshold within 80 ms and is maintained above the upper threshold for at least 40 ms. Speed constant	Slow, gentle electrode signal that does not reach the upper threshold within 80 ms. Grip force: Duration of the signal Speed constant	For patients with 1 strong electrode signal.

10.7 EVO Digital control with 1 electrode

INFORMATION

For this control mode, an electrode must be connected to input1.

This control mode is used for easier control of the prosthetic component for toddlers or children which cannot differentiate the different muscle signals.

Opening and closing of the prosthetic component is controlled by only one electrode. When sufficient muscle signal is generated (ON threshold exceeded), the electric hand opens at a constant speed, and closes again at a constant speed when the muscle signal disappears.

The gripping strength can be adjusted via the 560X3=* MyolinoSoft software.

Opening: Digitally; signal on the electrode opens the prosthetic component at constant speed.

Closing: Digitally; the prosthetic component closes at constant speed without a signal on the electrode.

Program	OPEN	CLOSE	Indications
EVO digital	Sustained electrode signal Speed constant	No electrode signal Grip force: maximum Speed constant	For patients with 1 weak electrode signal.

11 Opening the Electric Hand in Emergencies (Fig. 21)

The 7in1 Controller is equipped with an „EMERGENCY OPEN Switch“. This safety function makes it possible to quickly open the electric hand without myo signals.

Pressing the MyoEnergy Integral 757B25=1 or 757B35=* charging receptacle button (Fig. 21 arrow) opens the electric hand and the system turns off. This requires pressing the button for approx. 3 seconds (Fig. 21, arrow).

12 Cleaning and Care

Clean the product with a damp, soft cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed. Make sure the product is not damaged and that no liquid can penetrate into system components.

Then, dry the prosthesis component with a soft cloth.

13 Technical Data

Voltage	7.4 V
Weight	5 grams
Power Supply	757B25=1, 757B35=*
Average current consumption	4 mA
Operating temperature	0 °C to +60 °C (- 4 °F to +140 °F)
Storage temperature	-20 °C to +60 °C (- 4 °F to +140 °F)
Relative humidity	to 80%, non-condensing
Dimensions (LxWxH)	28 x 26,5 x 5 mm
Storage and transport temperature in packaging	-20° – +60°C
Storage and transport temperature without packaging	-20° – +60°C
Service life	5 years

14 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise.

The product may only be opened and repaired by authorised Otto Bock technicians.

15 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective-applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

16 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices.

This product has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to Appendix VII of the guidelines.

Bluetooth is a registered trademark of the respective holder.

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