Otto Bock®

QUALITY FOR LIFE

7772 Malleo Arexa

Ⓓ	Sprunggelenkorthese	4
ⒼⒷ	Ankle Orthosis	9
Ⓔ	Órtesis para articulación astrágalo-calcánea	13

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⚠ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

INFORMATION

Please read these instructions for use carefully and pay special attention to the safety information!

1 Description

1.1 Intended Use

The 7772 Malleo Arexa Ankle Orthosis is to be used **exclusively** for orthotic fitting of the lower limb.

1.2 Indications and Effects

The Malleo Arexa ankle orthosis is used for functional treatment of ligament instabilities, post-operatively after ligament ruptures, after ankle distortions, as well as for functional-conservative therapy of lateral collateral ligament ruptures. Indication must be determined by the physician.

The orthosis stabilizes the ankle. Eversion and inversion of the foot joint are restricted without negatively affecting plantar and dorsal flexion of the foot.

1.3 Contraindications

1.3.1 Absolute Contraindications

None.

1.3.2 Relative Contraindications

Skin diseases /injuries, inflammatory symptoms, prominent scars on body parts to which the orthosis is to be applied that are swollen, reddened and feel excessively warm; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the orthosis is to be applied; sensory and circulatory disorders.

1.4 Safety Instructions

INFORMATION

The patient is to be instructed in the proper use /care of the product.

The initial fitting and application of the product must be carried out by trained, qualified personnel only. The daily duration of use and period of application are dependent on medical indication by the physician.

CAUTION

Risk of injury as a result of improper use. The product is designed for use on one patient only. Parts to be fitted and those parts that come directly into contact with the skin can cause functional and hygienic risks if the orthosis is used by another person.

An orthosis / support applied too tightly to the body can cause local pressure and, in some cases, even restrict adjacent vessels or nerves. Do not apply the product too tightly. Consult a physician immediately if you experience unusual changes (such as increase in pain).

Improper changes to the product are not permitted.

NOTICE

Risk of damage due to incorrect environmental conditions. This product is not flame-resistant. Keep the product away from flames or other heat sources.

The product should not come into contact with grease or acidic agents, unguents and lotions, as this can reduce the useful life – especially of pads.

2 Handling

2.1 *Adapting and Application of the product*

The orthosis is universally applicable for the left and right sides. Sizing is according to shoe size.

The initial fitting and application of the orthosis must be carried out by qualified personnel.

The orthosis is delivered ready for trial fitting.

1. Open all closures on the brace (Fig. 1). Place it under the insole of the shoe. (If the insoles can't be taken out, see instruction 5). The back of the brace shall be placed against the heel counter of the shoe.
2. Close the upper and lower closures (Fig. 2). Tighten as hard as possible to still maintain a comfortable feeling.
3. Pull the diagonal strap through the D-ring and fasten it with the velcro (Fig. 3). The strap can be pulled through the D-ring either from the inside or outside the brace.
 - a. From the outside of the foot to the inside of brace, when the injury is on the outside of the foot (most common injury).
 - b. From the foot inside to the outside of brace, when the injury is on the inside foot.
4. Finally close the wide strap by pulling it around the ankle and tighten it so the brace feels stable. (Be careful not to tighten too hard thus stopping the blood flow).
5. If you have a shoe where you can not remove the insole, so use the adhesive EVA pad, which is included in the package (Fig. 5). Remove the protective film and place the pad on the plastic part that is in contact with the foot (Fig. 6).

2.2 *Instruction for use and care*

NOTICE

Damage caused by wear and tear. Before each use, the orthosis must be checked for functional reliability and for possible wear or damage. If a closure or any other part of the orthosis shows any signs of wear, the orthosis should no longer be used.

Improper changes to the product are not permitted.

INFORMATION

Neoprene products, time of use should not exceed 3-4 hours.

Products in other materials, time of use should not exceed 8 hours, unless otherwise is recommended by prescriber.

INFORMATION

If you are allergic to heat, we advise you to use a product that not contains neoprene.

The orthosis must be used in strict conformity with the indication. The physician's instructions must be strictly observed.

Material: Neopren / SBR, plastic (EVA)

Cleaning textile components:

1. Pull the textile component off the plastic shell.
2. Close all hook-and-loop closures.
3. Hand wash in 40 °C warm water with a standard mild detergent. Rinse thoroughly to remove all soap residues since they can cause skin irritations.
4. Stretch the product while still wet.
5. Allow to air dry. Do not expose to direct heat (from the sun, an oven or radiator) while drying.
6. Pull the dry SBR / neoprene sock over the plastic shell.

Cleaning plastic components: Wipe the outside with a slightly damp cloth.

Disposal after usage: The material used in the product is recyclable and should be disposed of accordingly.

3 Further Usage Restrictions

The product is designed for use on one patient only. Re-use of the product is not permitted.

The daily duration of use and period of application are dependent on the medical indication.

4 Liability

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

5 CE Conformity

This product meets the requirements of the 93/42/EEC guidelines for medical devices. This product has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH