616S5=1, 616S6=1 BetaSil

DE Gebrauchsanweisung	4
EN Instructions for use	7
FR Instructions d'utilisation	11
IT Istruzioni per l'uso	14
ES Instrucciones de uso	18
PT Manual de utilização	22
NL Gebruiksaanwijzing	25
SV Bruksanvisning	29
DA Brugsanvisning	32
NO Bruksanvisning	36
FI Käyttöohje	39
PL Instrukcja użytkowania	43
HU Használati utasítás	46
CS Návod k použití	50
RO Instrucțiuni de utilizare	53
HR Upute za uporabu	57
SL Navodila za uporabo	60
SK Návod na používanie	63
BG Инструкция за употреба	67
TR Kullanma talimatı	70
EL Οδηγίες χρήσης	74
RU Руководство по применению	78
JA 取扱説明書	82
ZH 使用说明书	85
KO 사용 설명서	88

INFORMATION

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INFORMATION

Date of the last update: 2015-01-14

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

616S5=1, 616S6=1 BetaSil is a silicone material with the following properties:

- 2-component silicone material that cross-links at room temperature
- Permanent bond with new silicone and silicone gel liners
- **No permanent bond with polyurethane liners or used (contaminated) silicone liners**
- High viscosity that prevents flowing down, even on vertical surfaces with application thicknesses less than 5 mm
- No irritating or allergenic substances of content, so that application on unbroken skin and closed scar tissue is possible (no open wounds)
- Easy to remove (does not adhere to the skin)

2 Application

2.1 Application

BetaSil is designed for use in orthopaedics technology. It may only be used by qualified personnel.

2.2 Area of Application

616S5=1, 616S6=1 BetaSil is used for filling up scars and concavities and for adding cushions. This achieves a complete form fit with the residual limb for a standard silicone liner.

The product is suitable for direct application to the skin. It must not be applied to open wounds. The product should not come into contact with body orifices such as the eyes, mouth, nostrils, ears or genitals. For safety reasons, a test should be performed prior to actual application, for example on the inside of the wrist. Do not use if erythema or skin rash occur.

3 Safety

3.1 Explanation of warning symbols

Warnings regarding possible risks of accident or injury.

Warnings regarding possible technical damage.

3.2 General Safety Instructions

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

Use of the product on damaged skin

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product on healthy skin.
- ▶ **Inform the patient.**

Improper processing of the material

Injuries due to careless processing

- ▶ Follow the processing instructions.

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

4 Scope of Delivery

The scope of delivery for the product is shown on page 1 (cover picture).

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

616S5=1, 616S6=1 BetaSil					
Fig.	Item no.		Quantity	Designation	Reference number
-	-	■	1	Instructions for Use	647G637=1
-	-	-	1	Dispenser	-
-	-	●	1	Refill Kit	616S6=1
				Consisting of:	
	-	-	3	BetaSil Cartridge	-
	-	-	15	Static Mixer	-

5 Usage

5.1 Notes on processing and storage

Processing

At a room temperature of 23°C, the processing time of BetaSil is about 5 minutes. Lower temperatures increase the processing time while higher temperatures reduce it. The processing time begins as soon as the two components mix. BetaSil cannot be shaped after the end of the processing time. Vulcanisation is complete after approximately 15-20 minutes at a room temperature of 23°C.

After vulcanisation, BetaSil will crosslink with itself as long as the surface is free of grease. New BetaSil can simply be applied to the old, with a recommended layer overlap of 1 cm. The shorter the time between the first and second application, the better the bond will be.

Once fully vulcanised, BetaSil is permanently stable and not altered by outside influences (temperature, air pressure, relative humidity). Unprocessed BetaSil can be stored for later use, provided the cartridge is sealed with the cover cap.

Storage

BetaSil in sealed cartridges can be stored at room temperature until the specified expiry date. At a storage temperature in excess of 25°C, the material may thicken and become unusable. Storage in a refrigerator is optimal, separate from unprocessed lamination resin.

5.2 Adapting a Liner

> Required tools and materials:

634A58 isopropyl alcohol, 749Y16 sanding sleeve 40 grit, 749Z12=160X16 sanding brush, 749F7 foam polishing roller

- 1) Determine the patient's residual limb dimensions to establish the liner size.
- 2) Check the fit of the liner on the patient and trim the liner if necessary (Fig. 9).
- 3) **INFORMATION: Do not use soap for cleaning under any circumstances. Soap prevents good cross-linking since it acts as a parting agent. Acetone, thinner or similar solvents irritate the skin and destroy the silicone.**
Clean the entire area of the liner where BetaSil will be applied with a degreasing cleaner.
- 4) Allow the liner to air dry for at least 5 minutes.
- 5) Load the BetaSil cartridge with a static mixer and insert it in the dispenser (see Page 10).
- 6) Discharge the first drop of BetaSil onto a cloth. This ensures that the components are mixed. At the same time, the drop serves as a reference sample to assess the state of vulcanisation.
- 7) Apply the BetaSil to the desired area of skin, e.g. scar tissue (Fig. 10).
- 8) Remove excess material with your finger (Fig. 11).
- 9) Pull on the liner.
- 10) **If necessary:** apply additional BetaSil (Fig. 12)

- 11) After 15 min to 20 min, remove the liner and turn it inside-out.
- 12) **NOTICE! Work carefully when sanding in order to avoid damage to the silicone material of the liner.**
 Trim the BetaSil with scissors or sand it with a sanding sleeve (see fig. 13).
 Then sand with the sanding brush (see fig. 14).
 Finally smooth with the polishing roller (see fig. 15).

5.3 Preparing the Dispenser

For homogenous mixing of the two silicone components, a static mixer is threaded onto the cartridge. Then the cartridge is inserted into the dispenser that pushes out the silicone.

Vulcanisation of the silicone begins as soon as the two components mix in the static mixer. When work interruptions (e.g. to change the cartridge) are shorter than the processing time of the silicone, the static mixer can continue to be used. Partly emptied cartridges can be resealed with the cover cap. The silicone in the cartridge remains unmixed and can still be used. The static mixer must be disposed of after the work is complete.

> **Required tools and materials:**

Dispenser, Static Mixer, BetaSil Cartridge

- 1) Rotate the cover cap on the cartridge anticlockwise by 90° and pull it off (Fig. 5 and 6).
- 2) Apply the static mixer to the cartridge and lock it with a 90° clockwise rotation (Fig. 7).
- 3) Pull the dispenser plunger back to the stop.
- 4) Put the cartridge into the dispenser with the notch at the bottom (Fig. 8).
- 5) Lock the dispenser.
- 6) Push the plunger into the cartridge by operating the handle until the static mixer fills with silicone.
- 7) If the cartridge is not empty after completing the work, reseat it with the cover cap.

6 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return and collection processes.

7 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

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Otto Bock HealthCare GmbH