
AxonRotation 9S503

DE Gebrauchsanweisung (Fachpersonal)	3
EN Instructions for use (qualified personnel)	7
FR Instructions d'utilisation (Personnel spécialisé)	12
IT Istruzioni per l'uso (Personale tecnico specializzato)	16
ES Instrucciones de uso (Personal técnico especializado)	21
PT Manual de utilização (Pessoal técnico)	26
NL Gebruiksaanwijzing (Vakmensen)	30
SV Bruksanvisning (Fackpersonal)	35
DA Brugsanvisning (Faguddannet personale)	40
NO Bruksanvisning (Fagpersonell)	44
FI Käyttöohje (Ammattihenkilöstö)	49
CS Návod k použití (Odborný personál)	53
JA 取扱説明書 (有資格担当者)	58

1

INFORMATION

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INFORMATION

Last update: 2017-03-01

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

The "9S503 AxonRotation" Axon-Bus system component will be referred to simply as the product below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Function

In combination with an Axon-Bus gripping component, e.g. the 8E500 Michelangelo Hand, the product facilitates fitting with components in the modular Ottobock Axon-Bus prosthesis system.

Pronation/supination of the Axon-Bus gripping component are possible from the neutral position and are carried out actively.

These functions support the patient during his or her daily activities and provide the highest rehabilitation value.

The product promotes a physiologically correct body posture and reduces the need for unnatural, compensating body movements.

The product can only be used in combination with the Axon-Bus system component "13E500 AxonMaster" version 1.4 or higher.

Power is supplied to the product by the battery which comes integrated in the socket.

Control settings for the product can only be carried out when all Axon-Bus components are connected. For further information, see the instructions for use provided with the Axon-Bus components.

To prevent mechanical damage, the coupling (connection between the motor and gripping component) slips when rotary movement is obstructed.

3 Application

3.1 Indications for use

The product is to be used exclusively for exoprosthetic fittings of the upper limbs and enables active pronation and supination of the Axon-Bus gripping component.

3.2 Use/Field of Application

In combination with myoelectrically controlled Axon-Bus gripping components in the modular Ottobock Axon-Bus prosthesis system, the product can be used from the transradial or transhumeral amputation level or, in case of dysmelia, for forearm or upper arm fittings.

3.3 Conditions of use

The product was developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

Permissible ambient conditions are described in the technical data (see Page 11).

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

3.4 Qualification

The product may be fitted only by qualified personnel authorised by Ottobock after completing the corresponding training.

4 Safety

4.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

CAUTION

The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

CAUTION

Independent manipulation of the product

Injury due to malfunction and resulting unexpected prosthesis actions.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

CAUTION

Distance to HF communication devices is too short (e.g. mobile phones, Bluetooth devices, WiFi devices)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
 - Mobile phone GSM 850/GSM 900: 0.50 m
 - Mobile phone GSM 1800/GSM 1900/UMTS: 0.35 m
 - DECT cordless phones incl. base station: 0.18 m
 - WiFi (routers, access points,...): 0.11 m
 - Bluetooth devices (third-party products not approved by Ottobock): 0.11 m

CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

► Clean the product and the contacts with a damp cloth and mild soap only.

5 Scope of Delivery and Accessories

Scope of Delivery

- 1 pc. 9S503 AxonRotation (active rotation unit)
- 4 pcs. M3 x 8 Torx screws
- 1 pc. Torx offset screwdriver
- 1 pc. 647G648 Instructions for Use (qualified personnel)

6 Combination possibilities

This product is to be combined exclusively with components in the Axon-Bus prosthetic system. Components of the Ottobock MyoBock system or components from other manufacturers cannot be used with this product.

7 Preparation for use

- 1) Use silicone grease to seal the plug connection.
- 2) Connect the product with the connecting cable.
- 3) Position the product so that the arrow mark points in the direction of the ulna (outside of the arm). (see fig. 1)
- 4) Screw the product to the lamination ring. (see fig. 1)

INFORMATION: Tighten the screws hand-tight and do not use a thread lock such as Loctite.

8 Maintenance

Since all moving mechanical parts are subject to wear, regular service inspections are recommended. These service inspections include an inspection of the product and, when necessary, the replacement of wear parts.

9 Legal information

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

9.3 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

The product meets the requirements under the RoHS Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

10 Appendix

10.1 Symbols Used

Declaration of conformity according to the applicable European directives

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

SN YYYW NNN

Serial number

Legal manufacturer

IP33

Protection against penetration of solid foreign objects with a diameter > 2.5 mm, protection against water spray

10.2 Technical data

Environmental conditions	
Storage and transport in original packaging	-20 °C/-4 °F to +40 °C/+104 °F
Storage and transport without packaging	-20 °C/-4 °F to +40 °C/+104 °F max. 80 % relative humidity, non-condensing
Operation	-10 °C/+14 °F to +60 °C/+140 °F max. 80 % relative humidity, non-condensing

General	
Reference number	9S503
Product service life	5 years
Active pronation from the neutral position	160°
Active supination from the neutral position	160°
Maximum torque	1.5 Nm
Idle speed	25 RPM
Rotation angle	150°/sec.
Weight	approx. 140 g/4.94 oz

Load limits	
Maximum static vertical load	700 N

Load limits	
Maximum dynamic vertical load	200 N

10.3 Axon

The term "Axon" stands for **A**daptive **eX**change **O**f **N**europlacement data. The Axon-Bus is an Ottobock innovation for the field of exoprosthetics: a data transmission system, derived from safety-related bus systems in the aviation and automobile industries. For the user this means enhanced safety and reliability because of a considerably reduced sensitivity to electromagnetic interference in comparison with conventional systems.

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Otto Bock Healthcare Products GmbH