
## 50R30 Lumbo TriStep

<b>DE</b> Gebrauchsanweisung .....	5
<b>EN</b> Instructions for use .....	9
<b>FR</b> Instructions d'utilisation .....	13
<b>IT</b> Istruzioni per l'uso .....	18
<b>ES</b> Instrucciones de uso .....	22
<b>PT</b> Manual de utilização .....	26
<b>NL</b> Gebruiksaanwijzing .....	31
<b>SV</b> Bruksanvisning .....	35
<b>NO</b> Bruksanvisning .....	39
<b>FI</b> Käyttöohje .....	43
<b>PL</b> Instrukcja użytkowania .....	47
<b>SK</b> Návod na používanie .....	51
<b>TR</b> Kullanma talimatı .....	56
<b>EL</b> Οδηγίες χρήσης .....	60
<b>RU</b> Руководство по применению .....	65
<b>ZH</b> 使用说明书 .....	70

**INFORMATION**

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### INFORMATION

Last update: 2015-09-30

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

## 1 Indications for use

The 50R30 Lumbo TriStep lumbar spine bridging brace with mobilizing function is intended **solely** for the orthotic fitting of the lumbar spine and is **designed for contact with intact skin only**.

## 2 Indications and effects

Post-operative stabilization of the lumbar spine (e.g. spondylodesis), post-operative correction of lordosis of the lumbar spine (e.g. after intervertebral disk surgery), facet syndrome, lumbar spinal canal stenosis, symptomatic stenoses of the lumbar intervertebral foramina, degenerative instability, spondylolysis, spondylolisthesis, herniated intervertebral disk, follow-up treatment of stable vertebral body fractures of the lower lumbar spine without neurological damage.

**The indication must be diagnosed by a physician.**

The 50R30 Lumbo TriStep is a bridging brace with a mobilizing function. By allowing for the removal of individual elements, the Lumbo TriStep enables systematic rehabilitation of the lumbar spine in three steps.


In the first stage of therapy, stabilization and relief are achieved through correction of lordosis of the lumbar spine.

In the second stage of therapy, the bridging element (5) is removed to make the bridging brace become a lumbar spine support brace. At this point, physiological lordosis will be nearly restored.

In the third stage of therapy, the support element (4) is removed from the brace. The remaining lumbar spine brace stabilizes and relieves the lumbar region.

## 3 Size selection

The brace is available in three sizes. Size selection is based on the pelvis circumference as follows:

	<i>Size</i>	<i>Pelvis circumference</i>
	<b>50R30=S</b>	65 – 80 cm / 25 – 31 inches
	<b>50R30=M</b>	80 – 100 cm / 31 – 39 inches
	<b>50R30=L</b>	100 – 125 cm / 39 – 49 inches

## 4 Scope of delivery and design (Fig. 1)

### 4.1 Scope of delivery

- |  |                       |
|--|-----------------------|
| (1) Back part  | (8) Clip rivet        |
| (2) Side parts attachable by hook and loop closure (width: 5 cm / 2 in; 7.5 cm / 3 in; 10 cm / 4 in) | (9) Strap tensioner   |
| (3) Abdominal closure with grip pockets  | (10) Perlon cord      |
| (4) Support element  | (11) Strap            |
| (5) Four-link bridging element   | (12) Mounting adapter |
| (6) Strap guide loop   | (13) Lumbar pad       |
| (7) Strap holder   |                       |

### 4.2 Design

The Lumbo TriStep is delivered ready for trial fitting.

The Lumbo TriStep is a three-component lumbar spine bridging brace that can be reduced in function by removing components in three steps. A support element (4) and an adaptable bridging element (5) are attached to the outside of the lumbar spine brace.

### System:

#### **Therapy stage 1 – lumbar spine bridging brace:**

Lumbar spine brace + support element (4) + bridging element (5)

#### **Therapy stage 2 – lumbar spine support brace:**

Lumbar spine brace + support element (4)

#### **Therapy stage 3 – lumbar spine brace:**

Lumbar spine brace + optional lumbar pad (13)

The 50R30 Lumbo TriStep consists of a back part (1), side parts of different widths, attachable by hook and loop closure (2), and an abdominal closure (3). The support element (4) is made of plastic and consists of vertical splints which decrease in length from the dorsal to the ventral region. The individual splint ends are attached to the brace with hook and loop closures.

Strap guide loops (6) connect the Perlon cord (10), which runs through the frame, to the wide strap (11) which ends on the abdominal closure (3). Strap holders (7) and strap tensioners (9) secure the support element (4) and the bridging element (5) on the abdominal closure (3). The strap tensioners (9) are also used to apply tension on both sides on the support element (4) and the bridging element (5).

Two mounting adapters (12) attach the bridging element (5), which consists of four vertical splints, to the support element (4).

In the third stage of therapy, a lumbar pad (13) can be attached to the inside of the brace by hook and loop connection, if desired.

## 5 Adapting

The initial fitting and application must be carried out by qualified personnel.

### **Therapy stage 1:**

Open the belt tensioners (8) and the abdominal closure (3) and apply the brace. For optimum adjustment of brace circumference and compression, the individual side parts (2), which are attached to the brace by hook and loop closure, can be removed. To do so, remove the unneeded side parts (equally wide) (2) on both sides (Fig. 7) and reattach the abdominal closure parts (3) (Fig. 8) (the red line serves as a guide).

If necessary, the length of the strap (11) can be adjusted at the strap tensioner (9) (Fig. 9).

### **Therapy stage 2:**

Open the mounting adapters (12) (Fig. 10) and remove the bridging element (5) (Fig. 11).

### **Therapy stage 3:**

To remove the support element (4), release the strap holders (7), which are attached to the abdominal closure parts (3) with clip rivets (8) (Fig. 12). Then remove the support element (4) with strap (11) from the brace (Fig. 13).

Once the brace has been applied to the patient, check for proper placement, compression and fit.

## 6 Application

Apply the lumbar brace with the red line on the inner side down such that the support element (4) is positioned over the entire width of the back (Fig. 4).

To close the brace, put your fingers into the pockets of the abdominal closure (3), pull the brace with equal force on both sides to the front, and attach the right part of the abdominal closure (3) to the left part (Fig. 5). Next, apply equal tension to the support element (4) on both sides using the strap tensioners (9) and attach the strap tensioners (9) to the upper right part of the abdominal closure (3) with the hook and loop connection (Fig. 6). Make sure that the attachment surfaces of the strap tensioners are completely on the abdominal closure and not in contact with clothing or skin.

Explain to the patient how to apply the brace and ensure its correct fit. Practice the application and removal of the brace with the patient.

## 7 Side effects

### **Attention!**

All aids – supports and braces – applied too tightly to the surface of the body may result in local pressure spots or, in some cases, even restrict the adjacent vessels or nerves.

## 8 Contraindications

### **Attention!**

In case of the following clinical pictures, the application and use of this type of aid is indicated only after consultation with the physician:

- Skin diseases and injuries of body parts to be fitted, especially in case of inflammatory symptoms. This also applies to prominent scars that are swollen, reddened, and feel excessively warm when touched.
- Lymphatic flow disorders, including indistinct soft tissue swellings in the region of the lower extremities.

## 9 Instructions for use

The brace must be used in strict conformity with the indication. The patient's acceptance and cooperation play an important part. It is essential to instruct the patient in the proper use of the brace!

**The physician's instructions must be strictly observed.**

The product should not come into contact with fatty or acidic agents, ointments, and lotions.

Improper changes to the brace are not permitted.

It is recommended that you always wear a T-shirt or undershirt under the brace.

The daily duration of use and total period of application must be determined by the attending physician.

Patients should be instructed to consult a physician immediately if they experience unusual changes.

## 10 Care instructions

Upon delivery of the product, the patient must be informed that only the textile part of the brace can be washed (after removal of the support element (4) and the bridging element (5), if attached). The support element (4) and the bridging element (5) are not washable.

The textile part of the brace can be washed at 30 °C (86 °F) using a standard mild detergent. Prior to washing, undo the clip rivets (8) of the strap holders (7) and remove the entire support element (4) attached by hook and loop connection, as well as the bridging element (5), if attached, from the brace. Close the abdominal closure (3). We recommend using a washing bag for washing the brace. Do not expose it to the sun. Do not use an oven or radiator to dry the brace. Allow the brace to air dry. The support element (4) and bridging element (5) can be cleaned with a damp cloth and neutral soap.

After drying, re-establish the hook and loop connection of the support element (4) and, if necessary, the bridging element (5) on the brace, and attach the strap holders (7) to the abdominal closure parts (3) with the clip rivets (8).

## 11 Restriction for use

The 50R30 Lumbo TriStep is designed for use on **one** patient. This brace is only suitable for contact with healthy skin.

### Attention!

Parts to be fitted and those parts that come directly into contact with the skin can cause functional and hygienic risks if the brace is used by another person. Functional risks are for example:

- Formation of cracks in the material
- Material fatigue due to alternating stress

The medical indication determines the period of application and the duration of use.

### Attention!

The 50R30 Lumbo TriStep contains metal components and must be removed when the patient is examined in an MRI scanner.

## 12 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 12.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 12.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH