

50R54 Lumbo Direxa Stable

DE	Gebrauchsanweisung	4
EN	Instructions for use	8
FR	Instructions d'utilisation	13
IT	Istruzioni per l'uso	18
ES	Instrucciones de uso	23
NL	Gebruiksaanwijzing	28
SV	Bruksanvisning	33
DA	Brugsanvisning	37
NO	Bruksanvisning	42
FI	Käyttöohje	47
IS	Notkunarleiðbeiningar	51
PL	Instrukcja użytkowania	56
HU	Használati utasítás	61
CS	Návod k použití	66
SK	Návod na používanie	70
BG	Инструкция за употреба	75
TR	Kullanma talimatı	80
ET	Kasutusjuhhis	84
RU	Руководство по применению	89
ZH	使用说明书	94

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INFORMATION

Date of last update: 2020-04-06

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.

- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50R54 Lumbo Direxa Stable lumbar brace.

2 Intended use

2.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the back and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

2.2 Indications

The indication must be determined by the physician.

- Lower back pain
- Lumbar spine syndrome
- Sciatica
- Myostatic insufficiency/muscular imbalance of the lumbar trunk muscles
- Degenerative changes of the lumbar spine
- Intervertebral disc protrusion L4-S1
- Stabilization after discectomy L4-S1
- Stabilization after spinal fusion L3-S1
- Aftercare of traumatic or osteoporotic fractures L3-L5
- Stabilization of the lumbar spine in case of osteoporosis L3-S1
- Moderate irritation of the sacroiliac joint
- Minor to moderate spondylolisthesis L3-L5

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

The following conditions require consultation with a physician: skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the support is fitted, tumors, lymphatic congestion – including unclear soft tissue swelling distal to the support, and dysesthesia in the trunk.

2.4 Mechanism of Action

The brace promotes sensorimotor function and activates muscular stabilization of the lumbar spine. It supports pain relief and relieves the lumbar spine by increasing intra-abdominal pressure and by providing external stabilization.

3 Safety

3.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

3.2 General safety instructions

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the Size

- 1) Measure the waist circumference approx. **2 cm** above the iliac crest.
- 2) Determine the brace size (see size chart).

4.2 Fitting and Application

CAUTION

Direct skin contact with the product

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

Optional – inserting the pad:

- ▶ **Optional – inserting the pad:** Use the hook and loop closure to attach the pad to the inside of the brace so that it is at least 1 cm above the seat surface when seated.

Fitting and application

- 1) **Fitting the metal rods** (see fig. 1): Remove the 4 dorsal metal rods from the pockets and fit them to the lumbar spine individually. Put the metal rods back in the pockets.
- 2) **Application** (see fig. 2): Place the brace around the torso. Place the fingers in the hand straps of the abdominal closure. Apply equal force to both sides to pull the brace forward.
- 3) **Hook and loop material** (see fig. 3): Attach the right part to the top of the left part of the abdominal closure. Remove the fingers from the left and then from the right hand strap. Make sure that the hook and loop surfaces are closed evenly against each other and are not in contact with clothing or skin.
- 4) **Increasing compression:** Pull the ends of the straps forward (see fig. 4) and attach them to the closure (see fig. 5).

Final steps

- 1) **Final inspection:** Verify the correct fit and proper positioning of the brace during handover to the patient. The brace is positioned correctly when the lower edge of the brace extends just below the groin.
- 2) **Instruct the patient:** Explain to the patient how to apply the brace and ensure that it is positioned correctly. Practice the application and removal of the brace with the patient.

4.3 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the brace regularly:

- 1) Close all hook-and-loop closures.
- 2) Recommendation: Use a laundry bag or net.
- 3) Wash the brace in warm water at **40 °C** with standard mild detergent. Do not use fabric softener. Rinse well.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

6.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

6.2 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

6.3 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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