

## 2R38=5/4R21=5/4R52=5

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Date of the last update: 2013-02-11

- Please read this document carefully.
- Follow the safety instructions.

## 1 Scope of Delivery (Fig. 1)

### 1.1 Single Component ■

–

### 1.2 Single-Component Pack ●

#### 4D4 Single-Component Pack (4R21=5 / 4R52=5)

Consisting of: 2 straight pins (3 + 4), 1 cap screw (5).

### 1.3 Minimum Order Quantity ▲

506G3=M8×12-V Set Screws (2)

4X28=4 (6)

## 2 Description

### 2.1 Intended Use

The 2R38=5 Tube Adapter and the 4R21=5, 4R52=5 Clamp Adapters are to be used **exclusively** for the prosthetic fitting of the lower limbs.

### 2.2 Field of Application

Field of application according to the **Ottobock Mobility System MOBIS®**:

≤ 136 kg / 300 lbs

Approved for a patient with a **weight of up to 136 kg/300 lbs.**

The tube adapter together with the corresponding clamp adapters comprise structural components for modular lower limb prostheses. The adapter combinations allow for controlled angular and linear adjustments in the sagittal and frontal planes, and internal and external rotations in the transverse plane.

For maximum structural strength the Tube Clamp Adapter must be placed proximally. The clamping slot must always be positioned anteriorly.

### Attention!

To achieve maximum stability for the transtibial prostheses, we recommend using a 34 mm tube adapter (e.g. 2R57 or 2R76) together with a 34 mm clamp adapter (e.g. 4R82 or 4R91). This applies in particular if the prosthesis is to withstand increased pressure due to higher activity levels, longer foot or transtibial lever arms or similar factors.

### Attention!

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water and acids. Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

**Please inform your patients.**

## 3 Technical Data

	<b>4R21=5</b>	<b>2R38=5</b>	<b>4R52=5</b>
Weight	130 g	275 g	75 g
System height	33 mm	97 mm	33 mm
Material	INOX stainless steel	Titanium	Titanium
Diameter	30 mm	30 mm	30 mm
Maximum patient weight	<b>136 kg / 300 lbs</b>	<b>136 kg / 300 lbs</b>	<b>136 kg / 300 lbs</b>
Mobility level	–	–	–

## 4 Handling

### 4.1 Connection to the Clamp Adapter

The tube adapter must be shortened according to the patient's measurements. To avoid deforming the tube, care must be taken to carry out the steps accurately as described below. Given the low tolerances of Ottobock adapters, a deformed tube will not fit easily into the clamp adapter.

#### Note:

Do not clamp the tube into a vise. This could cause damage or deformation.

**The cut surface must be perpendicular to the longitudinal axis of the tube** (Fig. 2). Use of the 719R2 Tube Cutter or equivalent is recommended.

**All inner and outer burrs must be removed from the cut edge** (Fig. 3). Use of the 718R1 Reamer or equivalent is recommended.

**Note:**

To ensure structural solidity, the remaining pieces of the shortened tube adapter are not to be used as components in the prosthesis.

Loosen the Cap Screw (Fig. 4.1) of the Tube Clamp Adapter by two turns, check, whether the plastic ring has been inserted into the Tube Clamp Adapter completely and with the chamfer on the ring's outer diameter first, and then insert the tube 5 mm deep into the Tube Clamp Adapter while slightly rotating the tube. Next, place the Tube Clamp Adapter with the tube inserted vertically on the worktop and fully insert the Tube Adapter into the Adapter until the tube is in contact with the base rim inside the clamp.

**Attention!**

The interior surface of the Tube Clamp Adapter and the exterior surface of the Tube Adapter must be cleaned with 634A3 (Aceton) in the insertion area prior to final assembly (insertion of the Tube Adapter into the Tube Clamp Adapter).

**Warning**

For final use, the Cap Screw (Fig. 4.1) must be torqued to 5 Nm with the 710D1 Torque Wrench. Failing to comply with this instruction may lead to permanent over-extension or breaking of the screw.

## **4.2 Maintenance Instructions**

**Information**

As a basic principle, all Ottobock modular adapters are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

## **5 Adjustment Screws**

Adjustments may be made through adjusting the four 506G3 Adjustment Set Screws during fitting and at any time after the prosthesis has been

completed. The optimal position of the components may be preserved when exchanging parts or when servicing the prosthesis by removing only the two most deeply set adjacent screws.

Once positioning has been determined, replace any set screws that are too long (to prevent foam cover damage) or too short (lacking strength).

Screws are available in the following sizes:

<i>Article Nr.</i>	<i>Length</i>
<b>506G3=M8x10</b>	10 mm
<b>506G3=M8x12-V</b>	12 mm
<b>506G3=M8x14</b>	14 mm
<b>506G3=M8x16</b>	16 mm

### Use a 710D4 Torque Wrench to tighten screws

#### **Attention!**

(Applies to the application of the Tube Adapter 2R38=5 and Tube Clamp Adapters 4R21=5; 4R52=5)

**For testing:** with the torque wrench 710D4, screw in the set screws **506G3** of the Adapter clockwise at **10 Nm** (*7.4 lbf. ft / 88.5 lbf. in*).

**For finalizing the prosthesis: secure** the set screws 506G3 with Loctite 636K13, tighten the set screws clockwise at **10 Nm** (*7.4 lbf. ft / 88.5 lbf. in*) first and then tighten completely at **15 Nm** (*11.1 lbf. ft / 133 lbf. in*).

## 6 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components in accordance with the Ottobock Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

## 7 Warranty

Ottobock warrants this device / these devices from the date of purchase. The warranty covers defects that can be proven to be a direct result of

flaws in the material, production or construction. Because warranty conditions may differ by country, please contact your local Ottobock distribution company for more information.

## 8 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH