ottobock.

4R63=5, 4R63=5.1, 4R100=5

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INFORMATION

Date of the last update: 2013-12-12

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

Explanation of symbols

CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting/use.

1 Indications for use

The socket adapters 4R63=5/4R63=5.1 (stainless steel) and 4R100=5 (titanium) are to be **exclusively** used for the prosthetic fitting of amputations of the lower limb.

Field of application according to the **Mobility System MOBIS:**

4R63=5, 4R63=5.1

≤ 136 kg / 300 lbs

Approved for a patient **with a weight of up to 136 kg/300 lbs.**

4R100=5

≤ 102 kg / 225 lbs

Approved for a patient **with a weight of up to 102 kg/225 lbs.**

The socket adapter is recommended for below-knee prostheses and is integrated as lamination anchor into the laminate of the prosthetic socket.

⚠ CAUTION

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- Observe the safety information in these instructions for use.
- Make the patient aware of all safety instructions marked with "**In-form the patient**".

⚠ CAUTION

Excessive strain on the product

Fall due to breakage of load-bearing components

- Use the prosthesis components according to the MOBIS classification (see section "Area of Application").
- **Inform the patient.**

⚠ CAUTION

Unallowable combination of prosthesis components

Fall due to breakage or deformation of the product

- Only combine the product with prosthesis components that are approved according to the section "Combination Possibilities".
- Based on the instructions for use of the prosthesis components, verify that they may be combined with each other.

⚠ CAUTION

Exceeding the service life and reuse on another patient

Fall due to loss of functionality as well as damage to the product

- Ensure that the approved service life is not exceeded (see section "Service Life").
- Only use the product for a single patient.
- **Inform the patient.**

CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- Use caution when working with the product.
- Check a damaged product for proper functionality and readiness for use.
- In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- Take suitable measures if required (e. g. repair, replacement, inspection by the manufacturer's customer service etc.).
- **Inform the patient.**

CAUTION

Wrong processing of the anchor arms of the laminating anchor

Injuries from falling and breakage to anchor arms

- Avoid excessive or too frequent bending of the anchor arms
- Do not bend anchor arms made of titanium.
- Use the 711S4* bending iron for bending.

NOTICE

Use under environmental conditions which are not allowed

Damage to the product

- Do not use the product under prohibited environmental conditions (see the chapter "Environmental Conditions").
- If the product has been exposed to prohibited environmental conditions, check it for damage.
- If damage is apparent or in case of doubt, do not continue using the product.
- Take suitable measures if applicable (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).
- **Inform the patient.**

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthesis components relative to each other and by the development of noises.

2 Technical data

	4R63=5.1	4R63=5	4R100=5
System height	-7 mm	-7 mm	-7 mm
Weight	95 g	89 g	51 g
Material	INOX*	INOX*	Titan
Max. patient weight	136 kg / 300 lbs	136 kg / 300 lbs	102 kg / 225 lbs

*stainless steel

3 Handling

3.1 Recommendations for the layup

- Pull a soaked PVA bag 99B81 over the model.
- Pull over 2 layers of Perlon stockinette 623T3.
- Place offset 3 layers of carbon fibre cloth 616G12 (e.g. 15 x15 cm) at the distal limb end.
- Pull over 2 layers of Perlon stockinette 623T3.
- Insert the enclosed 4X52 Lamination Aid into the socket adapter and place the adapter on the distal end of the plaster model. Then align the socket adapter perpendicularly according to the marking of the 4X3 Lamination Protection.
- Thread 699B1 Fiberglass Roving through the bore holes of the anchor arms and let it hang in loops.
- Fixate (reinforce) the lamination anchor with a mixture of sealing resin and talcum.
- Place offset 3 layers of carbon fibre cloth 616G12 over the lamination anchor arms.
- Pull over 1 layer of Perlon stockinette 623T3.
- Pull over 2 layers of carbon stockinette 616G15.
- Pull over 2 layers of Perlon Stockinette 623T3 and tie off.

- Pull on a soaked PVA bag 99B81 and tie off proximally.
- Laminate with 617H55 C-Orthocryl®.

The enclosed 4X3 Lamination Protection covers the adjustment pyramid during the lamination procedure with 617H19 Orthocryl® Lamination Resin 80:20 or 617H55 C-Orthocryl® and provides an accurately fitting brim of the laminate.

3.2 Maintenance instructions

INFORMATION

As a basic principle, all Ottobock modular adapters are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

4 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

4.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

4.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

4.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

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Otto Bock HealthCare GmbH