## 2R36, 4R156, 4R156=1, 4R156=2


<b>DE</b> Gebrauchsanweisung .....	4
<b>EN</b> Instructions for Use .....	11
<b>FR</b> Instructions d'utilisation .....	18
<b>IT</b> Istruzioni per l'uso .....	25
<b>ES</b> Instrucciones de uso .....	32
<b>PT</b> Manual de utilização .....	40
<b>NL</b> Gebruiksaanwijzing .....	47
<b>SV</b> Bruksanvisning .....	54
<b>DA</b> Brugsanvisning .....	61
<b>NO</b> Bruksanvisning .....	68
<b>FI</b> Käyttöohje .....	75
<b>PL</b> Instrukcja użytkowania .....	82
<b>HU</b> Használati utasítás .....	89
<b>CS</b> Návod k použití .....	96
<b>RO</b> Instrucțiuni de utilizare .....	103
<b>HR</b> Upute za uporabu .....	110
<b>SL</b> Navodila za uporabo .....	117
<b>SK</b> Návod na používanie .....	124
<b>BG</b> Инструкция за употреба .....	131
<b>TR</b> Kullanma talimatı .....	138
<b>EL</b> Οδηγίες χρήσης .....	145
<b>RU</b> Руководство по применению .....	153
<b>JA</b> 取扱説明書 .....	161
<b>ZH</b> 使用说明书 .....	167
<b>KO</b> 사용 설명서 .....	173

**4R156**

**4R156=1**

**4R156=2**


**2R36**

1


2


3

**4R156, 4R156=1, 4R156=2**

**2R36**

# 1 Produktbeschreibung

## INFORMATION

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### INFORMATION

Date of the last update: 2013-01-29

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

### 1.1 Construction and Function

The 4R156\* Tube Clamp Adapters are specially designed for an adjustable connection between the pyramid adapter of the prosthetic hip joint and the pyramid adapter of the prosthetic knee joint or the rotator. The 2R36 Tube is used to connect two 4R156\* Tube Clamp Adapters.

### 1.2 Combination Possibilities

For combination possibilities, please see the catalogue 646K2\* or consult the manufacturer.

## 2 Application

### 2.1 Intended Use

The product is intended solely for lower limb prosthetic fittings.

### 2.2 Area of Application

Area of application according to the MOBIS® mobility system:

Approved for a body weight of up to **150 kg**.

The 2R36 Tube may be used only as a connecting element between two 4R156\* Tube Clamp Adapters.

## 2.3 Environmental Conditions

Allowable Environmental Conditions
Temperature range for use: -10 °C to +60 °C (14 °F to 140 °F)
Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e. g. talcum)

## 2.4 Service Life

This component has been manufacturer-tested with 3 million load cycles according to ISO 10328. Depending on the level of user activity, this corresponds to a service life of 3 to 5 years.

## 3 Safety

### 3.1 Explanation of warning symbols

 <b>CAUTION</b>	Warnings regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warnings regarding possible technical damage.

### 3.2 General Safety Instructions

 <b>CAUTION</b>
<b>Using the product without following the instructions for use</b> Deteriorating health condition and damage to the product because of failure to follow the safety information
<ul style="list-style-type: none"><li>▶ Observe the safety information in these instructions for use.</li><li>▶ Make the patient aware of all safety instructions marked with "<b>Inform the patient</b>".</li></ul>

 <b>CAUTION</b>
<b>Excessive strain on the product</b> Fall due to breakage of load-bearing components
<ul style="list-style-type: none"><li>▶ Use the prosthesis components according to the MOBIS classification (see section "Area of Application").</li><li>▶ <b>Inform the patient.</b></li></ul>

**⚠ CAUTION**

**Unallowable combination of prosthesis components**

Fall due to breakage or deformation of the product

- ▶ Only combine the product with prosthesis components that are approved according to the section "Combination Possibilities".
- ▶ Based on the instructions for use of the prosthesis components, verify that they may be combined with each other.

**⚠ CAUTION**

**Use under unallowable environmental conditions**

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if applicable (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).
- ▶ **Inform the patient.**

**⚠ CAUTION**

**Exceeding the service life**

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded (see section "Service Life").
- ▶ **Inform the patient.**

**⚠ CAUTION**

**Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check a damaged product for proper functionality and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).

- ▶ Take suitable measures if required (e. g. repair, replacement, inspection by the manufacturer's customer service etc.).
- ▶ **Inform the patient.**

### Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthesis components relative to each other and by the development of noises.

## 4 Scope of Delivery

The scope of delivery for the product is shown on page 2 (Figure 1).

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

2R36 Tube					
Fig.	Item No.		Quantity	Designation	Reference Number
-	-	■	1	Instructions for Use	647G748=1

4R156, 4R156=1, 4R156=2 Tube Clamp Adapter					
Fig.	Item No.		Quantity	Designation	Reference Number
-	-	■	1	Instructions for Use	647G748
4	1	▲	4	Set screw	506G3=M8x12-V
4	-	●	-	Single-components pack	4D28
	<b>Consisting of:</b>				
	2	-	1	Cap screw	-
	3	-	1	Straight pin	-
	4	-	1	Straight pin	-
4	5	■	1	Plastic ring	4X28=3

## 5 Preparation for Use

### CAUTION

#### Incorrect alignment or assembly

Injuries due to fall resulting from damaged prosthesis components

- ▶ Observe the alignment and assembly instructions.

## 5.1 Adapting and Mounting the Tube

### CAUTION

#### Incorrect processing of tube

Fall due to damage to the tube

- ▶ Do not clamp the tube into a vice.
- ▶ For shortening the tube, use only a tube cutter.

### CAUTION

#### Incorrect mounting of the tube

Fall due to breakage of load-bearing components

- ▶ Slide the tube all the way to the stop in the intended prosthesis component when mounting.

### CAUTION

#### Improper assembly of the screw connections

Fall due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the installation torques provided.
- ▶ Observe the instructions for securing the screw connections and the use of the correct length.

#### > **Recommended tools and materials:**

710D4 Torque Wrench, 719R3 Tube Cutter, 718R1 Tube Deburrer, degreasing cleaner (e.g. 634A3 Acetone)

- 1) Use the tube cutter to shorten the tube (see Fig. 5) to fit the user's measurements (see Fig. 2).
- 2) Carefully deburr the inside and outside of the tube with the 718R1 Deburring Knife (Fig. 3).
- 3) Loosen the tube clamp adapter cylinder head screw approx. **2 turns** (see Fig. 4, Item 2).
- 4) For definitive mounting, clean the insertion area of the connecting component on the inside and the tube on the outside using a degreasing cleaner.
- 5) Slide the tube **5 mm** into the connecting component while turning it slightly.
- 6) Take the connecting component with the inserted tube and place it on the work desk in a vertical position.
- 7) Slide the tube into the connecting component until the tube rests against the bottom of the clamp fitting.

- 8) Using the torque wrench, tighten the cap screw (for definitive mounting): **10 Nm** - see Fig. 6).

## 5.2 Mounting in the Modular Prosthesis

### CAUTION

#### **Improper assembly of the screw connections**

Fall due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the installation torques provided.
- ▶ Observe the instructions for securing the screw connections and the use of the correct length.

### Mounting

#### > **Recommended tools and materials:**

710D4 Torque Wrench, 636K13 Loctite®

#### 1) **For definitive mounting:**

Use Loctite® to secure the set screws (see Fig. 4, Item 1).

#### 2) Turn in the set screws.

#### 3) **For trial fitting:**

Use the torque wrench to tighten the set screws (**10 Nm**).

#### **For definitive mounting:**

Pre-tighten the set screws with the torque wrench (**10 Nm**) and then tighten (**15 Nm** - see Fig. 6).

#### 4) Use suitable set screws to replace any set screws that extend out too far or screw in too deeply (see selection table).

**Selection Table for Set Screws**

Reference Number	Length (mm)
506G3=M8x10	10
506G3=M8x12-V	12
506G3=M8x14	14
506G3=M8x16	16

### Adjustment

The 4 set screws can be used to make static adjustments at any time during alignment, trial fittings, and even after the prosthesis is finished.

## **Replacement and Disassembly**

When replacing or disassembling the prosthesis component, the previously adjusted position can be maintained if the two deepest screwed-in set screws, which are NOT located opposite to each other, are removed.

## **6 Maintenance**

- The prosthesis components should be inspected after the first 30 days of use.
- Inspect the entire prosthesis for wear during normal consultations.
- Conduct annual safety inspections.

## **7 Disposal**

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return and collection processes.

## **8 Legal Information**

### **8.1 Liability**

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components in accordance with the MOBIS Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorised by the manufacturer.

### **8.2 CE Conformity**

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

### **8.3 Warranty**

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction. Because warranty conditions may differ by country, please contact your local manufacturer distribution company for more information.

## 9 Technical Data

Reference Number	4R156	4R156=1	4R156=2	2R36
Weight [g]	140	165	175	210
System height [mm]	36	37	38	–
Min. system height [mm]	–	–	–	73
Max. system height [mm]	–	–	–	380
Material	Titanium			Aluminium
Diameter [mm]	34			
Tilt angle	10°	20°	30°	–
Max. body weight [kg]	150			

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Otto Bock HealthCare GmbH