4F34

DE	Gebrauchsanweisung	2
EN	Instructions for use	4
FR	Instructions d'utilisation	6
IT	Istruzioni per l'uso	8
ES	Instrucciones de uso	10
PT	Manual de utilização.....	12
SV	Bruksanvisning.....	14
NL	Gebruiksaanwijzing.....	16
DA	Brugsanvisning.....	18
PL	Instrukcja użytkowania.....	20
HU	Használati utasítás	22
CS	Návod k použití.....	24
TR	Kullanma talimatı.....	26
EL	Οδηγίες χρήσης	28

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1 Application

Last update: 2020-08-03


- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document in a safe place.

The 4F34 Lock Control Unit is to be used **exclusively** in lower limb components.

2 Description and function

A safe, easily operated manual knee lock control is very important for the geriatric patient.

The 4F34 Lock Control Unit may be used with cable controlled knee lock mechanisms that utilize a cable travel of no more than 11 mm maximum.

Lock Control Unit without cover,
left side, lock lever anterior, pressure activated

In the above illustration the knee lock is released when the lock lever (2) is pushed distally. Re-routing the cable (5) through the other entrance hole (6) will cause release of the lock when the lever is pulled proximally.

To change the position of the lock lever (2) from anterior to posterior, first pull off the cover (not shown). Remove the oval head screw (3) and pull the lock lever out of the housing (1). Rotate the lock lever 180° and reassemble in the reverse order.

The cable is secured by two set screws in the bushing (4) of the lock lever. A hexagon wrench for use in securing and adjusting cable length is enclosed for your convenience.

3 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

4 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

4.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

4.2 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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