1S101 SACH+


| | | | |
|---|----|---|-----|
| DE Gebrauchsanweisung | 4 | CS Návod k použití..... | 70 |
| EN Instructions for use | 9 | RO Instrucțiuni de utilizare..... | 75 |
| FR Instructions d'utilisation..... | 14 | HR Upute za uporabu | 81 |
| IT Istruzioni per l'uso | 19 | SL Navodila za uporabo | 86 |
| ES Instrucciones de uso | 25 | SK Návod na používanie | 91 |
| PT Manual de utilização..... | 30 | BG Инструкция за употреба | 96 |
| NL Gebruiksaanwijzing..... | 35 | TR Kullanma talimatı | 102 |
| SV Bruksanvisning | 40 | EL Οδηγίες χρήσης | 107 |
| DA Brugsanvisning | 45 | RU Руководство по применению | 112 |
| NO Bruksanvisning | 50 | JA 取扱説明書..... | 118 |
| FI Käyttöohje..... | 55 | ZH 使用说明书..... | 123 |
| PL Instrukcja użytkowania..... | 60 | KO 사용 설명서..... | 127 |
| HU Használati utasítás | 65 | | |

1

2

3

4

1 Produktbeschreibung

Deutsch

INFORMATION

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INFORMATION

Date of last update: 2021-08-02

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 1S101 SACH+ prosthetic foot (Solid Ankle Cushion Heel) is suitable for use in both modular and exoskeletal prostheses. The functional properties are achieved through the combination of a contoured core and functional foam.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

| Allowable ankle blocks and screw connections | |
|--|---|
| Foot size [cm] | Reference number (ankle block and screw connection) |
| 24 to 30 | 2K34=30 and 2Z22=M10 |

| Allowable foot adapters | | |
|-------------------------|----------------|---------------------------------|
| Body weight [kg] | Foot size [cm] | Reference number (foot adapter) |
| ≤100 | 22 to 30 | 2R8=M10, 2R31=M10, 2R54=M10 |
| 101 to 125 | 26 to 27 | 2R8=M10, 2R31=M10 |
| | 28 to 30 | 2R31=M10 |

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 1 (indoor walker) and mobility grade 2 (restricted outdoor walker).

- The maximum approved body weight is specified in the technical data (see page 14).

2.3 Environmental conditions

Storage and transport

Temperature range -20°C to $+60^{\circ}\text{C}$ (-4°F to $+140^{\circ}\text{F}$), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

Allowable environmental conditions

Temperature range: -10°C to $+45^{\circ}\text{C}$ (14°F to 113°F)

Chemicals/liquids: fresh water, soapsuds, chlorine water

Moisture: Submersion: max. 1 h in 2 m depth, relative humidity: no restrictions

Solids: Dust, occasional contact with sand

Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage (see page 13).

Prohibited environmental conditions

Solids: highly hygroscopic particles (e.g. talcum), dust in high concentrations (e.g. construction site), intensive contact with sand

Chemicals/liquids: salt water, perspiration, urine, acids, continuous use in liquid media

2.4 Lifetime

The product was tested by the manufacturer with 2 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 3 years.

3 Safety

3.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION!

Risk of injury and risk of product damage

- ▶ To avoid the risk of injury and product damage, do not use the product beyond the tested lifetime.
- ▶ To avoid the risk of injury and product damage, only use the product for a single patient.
- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 9).
- ▶ Note the combination possibilities/comboination exclusions in the instructions for use of the products.

NOTICE!

Risk of product damage and limited functionality

- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ Clean the product promptly after contact with salt water, water containing chlorine/soap or abrasive substances (e.g. sand).

- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

Signs of changes in or loss of functionality during use

Decreased forefoot resistance or changes in roll-over behaviour are noticeable indications of loss of functionality.

4 Scope of delivery

| Quantity | Designation |
|----------|----------------------|
| 1 | Instructions for use |
| 1 | Prosthetic foot |

Spare parts/accessories (not included in the scope of delivery)

| Fig. | Item | Designation | Reference number |
|------|------|---|-----------------------------------|
| 2 | 1 | Foot adapter with screw connection | 2R8*, 2R31*, 2R54* |
| 2 | 2 | Connection plate | 2R14 |
| 3 | – | Ankle block (exoskeletal design) | 2K34=* |
| 4 | – | Screw connection for ankle block (exoskeletal design) | 2Z22=M10 |
| – | – | Ottobock special adhesive with hardener | 636W28 or 636W18 and 636W19 |

5 Preparing the product for use

⚠ CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

INFORMATION

- ▶ Not all of the materials may be available in your country. In this case, please contact the local branch of the manufacturer to obtain information on alternative materials.

5.1 Installing the foot adapter

- ▶ Install the foot adapter on the prosthetic foot according to its instructions for use.

5.2 Installing the ankle block (exoskeletal design)

> **Required materials:** Screw connection (see page 11), 710D20 torque wrench, sandpaper (grain size 40), 636W28 Ottobock special adhesive and hardener, suitable solvent (e.g. 634A28)

- 1) Insert the threaded socket of the screw connection in the proximal side of the hole in the ankle block.
- 2) Insert the screw on the foot sole side in the prosthetic foot and screw the ankle block to the prosthetic foot (torque: **30 Nm**).
- 3) Conduct the fitting of the ankle block.
- 4) Loosen the screw connection between the prosthetic foot and the ankle block.
- 5) Sand the ankle block according to the contour of the prosthetic foot.
- 6) Roughen the connection surface of the prosthetic foot in an a – p direction using the sandpaper.
- 7) Remove the dust from the sanded surfaces.
- 8) Mix the adhesive for the connection between the ankle block and the prosthetic foot.

- 9) **INFORMATION: Do not put adhesive in the hole for the screw connection.**
Apply adhesive evenly to the connection surfaces of the ankle block and prosthetic foot.
- 10) Align the ankle block and prosthetic foot with each other.
- 11) Insert the screw on the foot sole side in the prosthetic foot and screw the ankle block to the prosthetic foot (torque: **30 Nm**).
- 12) Wipe off adhesive that comes out of the joint using a cloth.
- 13) Remove any remnants of the adhesive using a suitable solvent.
- 14) Allow the adhesive to cure (curing time: **16 h**).

5.3 Alignment

NOTICE

Grinding the prosthetic foot

Premature wear resulting from damage to the prosthetic foot

- ▶ Do not grind the prosthetic foot.

5.3.1 Bench Alignment

TT bench alignment

| Bench alignment process | |
|---|---|
| Required materials: 662M4 goniometer, 743S12 heel height measuring apparatus, 743A80 50:50 gauge, alignment apparatus (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly) | |
| Perform the assembly and alignment of the prosthetic components in the alignment apparatus according to the following specifications: | |
| Sagittal plane | |
| ① | Heel height: Effective heel height (shoe heel height - sole thickness of forefoot) + 5 mm |
| ② | Exterior foot rotation: approx. 5° |
| ③ | Anterior placement, middle of the prosthetic foot to the alignment reference line: 30 mm |

| Bench alignment process | |
|-------------------------|--|
| ④ | Connect the prosthetic foot and prosthetic socket using the chosen adapters. The instructions for use of the adapters must be observed. |
| ⑤ | Determine the centre of the prosthetic socket with the 50:50 gauge. Align the prosthetic socket centrally to the alignment reference line. Socket flexion: Individual residual limb flexion + 5° |
| Frontal plane | |
| ⑥ | Alignment reference line of prosthetic foot: between big toe and second toe Alignment reference line of prosthetic socket: along the lateral patella edge |
| ⑦ | Observe the abduction or adduction position. |

TF bench alignment

- ▶ Observe the information in the prosthetic knee joint instructions for use.

5.3.2 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.
- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) may be requested from Ottobock.

5.3.3 Dynamic Trial Fitting

- Adapt the alignment of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) to ensure an optimum gait pattern.

- **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal plane is achieved when the leg begins to bear weight after the heel strike. Avoid medial movement of the knee joint. If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction. If the medial movement occurs in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.

5.4 Optional: Installing the foam cover

The foam cover sits between the prosthetic socket and prosthetic foot. It is cut longer in order to compensate for the movements of the prosthetic foot and prosthetic knee joint. During flexion of the prosthetic knee joint, the foam cover undergoes posterior compression and anterior elongation. The foam cover should be stretched as little as possible in order to increase its service life. There is a connecting element (such as a connection plate, connection cap or connection cover) on the prosthetic foot.

> **Required materials:** degreasing cleaner (e.g. 634A58 isopropyl alcohol), 636N9 contact adhesive or 636W17 plastic adhesive

- 1) Measure the length of the foam cover on the prosthesis and add the length allowance.

TT prostheses: Distal allowance for movement of the prosthetic foot.

TF prostheses: Allowance proximal of the knee rotation point for flexion of the prosthetic knee joint and distal allowance for movement of the prosthetic foot.

- 2) Cut the pre-shaped foam cover to length and fit it in the proximal area on the prosthetic socket.
- 3) Pull the foam cover over the prosthesis.
- 4) Set the connecting element onto the footshell or prosthetic foot. Depending on the version, the connecting element engages in the edge or rests on the foot adapter.
- 5) Install the prosthetic foot on the prosthesis.
- 6) Mark the outer contour of the connecting element on the distal face of the foam cover.

- 7) Disassemble the prosthetic foot and remove the connecting element.
- 8) Clean the connecting element using a degreasing cleaner.
- 9) Glue the connecting element onto the distal face of the foam cover according to the marked outer contour.
- 10) Let the glue dry (approx. **10 minutes**).
- 11) Install the prosthetic foot and adapt the exterior cosmetic shape. Take into account compression caused by cosmetic stockings or SuperSkin.

6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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