


3R95

DE Gebrauchsanweisung	4
EN Instructions for use	10
FR Instructions d'utilisation	16
IT Istruzioni per l'uso	23
ES Instrucciones de uso	29
PT Manual de utilização	36
NL Gebruiksaanwijzing	42
SV Bruksanvisning	49
DA Brugsanvisning	73
NO Bruksanvisning	61
PL Instrukcja użytkowania	67
HU Használati utasítás	73
TR Kullanma talimatı	79
EL Οδηγίες χρήσης	85
RU Руководство по применению	92
JA 取扱説明書	99
ZH 使用说明书	104

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English

Last update: 2021-01-15

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document in a safe place.

1 Single components

1.1 Scope of delivery

1 pc. of 3R95 Modular Lightweight Knee Joint

1 pc. of Instructions for Use

1.2 4D17 Single-component pack

For 3R95 Modular Lightweight Knee Joints when using the 4R57 Rotation Adapter

1 pc. of Oval countersunk head screw

1 pc. of Rubber stop

2 Description

2.1 Indications for use


The 3R95 Modular Lightweight Single-Axis Knee Joint with integrated miniaturised hydraulic swing phase control are to be used **exclusively** for the exoprosthetic fitting of the lower limb.

Amputation level: knee disarticulation, transfemoral amputation.

2.2 Field of application

Application of 3R95

Field of application according to the **Ottobock MOBIS Mobility System**:

Recommended for **mobility grades 3 and 4**
(unrestricted outdoor walkers and unrestricted outdoor walkers with especially rigorous demands).

Approved for a patient weight of up to 150 kg/330 lbs.

2.2.1 Application conditions

Primarily recommended for transfemoral amputees who can actively contribute to stance phase stability and wish to have a harmonious swing phase control. The light weight and low structural height are special advantages when fitting women and adolescents.

The 3R95 Modular Lightweight Joints have been developed for everyday use and must **not** be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

2.3 Environmental conditions

CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- Do not expose the product to unallowable environmental conditions (see the table „Unallowable environmental conditions“ in this section).
- If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Allowable environmental conditions

Temperature range:

Use: -10 °C – +45 °C

Storage, Transportation: -20 °C – +60 °C

Relative humidity: 20 % – 90 % (Non-condensing)

Commercially available, solvent-free cleaning agents

Unallowable environmental conditions

Storage/Transportation: Mechanical vibrations, Impacts

Hygroscopic particles (e.g. Talcum powder), Dust, Sand, Fresh water, Salt water, Acids, Sweat, Urine

Cleaning agents containing solvents

2.3.1 Reuse and lifetime

CAUTION

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- Only use the product for a single patient.

CAUTION

Exceeding the lifetime

Fall due to change in or loss of functionality and damage to the product

- Ensure that the maximum lifetime defined in this section is not exceeded.

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

2.4 Safety instructions

2.4.1 General safety instructions

- The 3R95 Modular Lightweight Knee Joint may not be used for bathing prostheses.
- Please follow all indicated torque values.
- Use the Patient Information section of this pamphlet to inform your patient about the function and use of the modular lightweight knee joint.

2.4.2 Patient information

- The modular lightweight knee joint should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- Do not allow foreign particles or liquids to get into the modular lightweight knee joint. Avoid exposing prosthetic components to corrosive elements such as fresh water, salt water, acids and other liquids. Using this medical product in such environmental conditions will render all claims against Ottobock null and void.
- The modular lightweight knee joint as well as damaged components may be opened or repaired only by certified Ottobock Service technicians.
- The environment and usage conditions may compromise the functionality of the knee joint. To reduce the risk for the patient, the knee joint must not be used after noticeable functional changes occur. Noticeable functional changes may include poor response, incomplete extension, diminished swing phase control or stance phase instability, abnormal noises, etc. Please contact a specialist workshop to have the prosthesis inspected if you discover damage or noticeable functional changes.

2.5 Prosthetist qualifications

The fitting of patients with the 3R95 Modular Lightweight Knee Joints may only be carried out by a prosthetist.

2.6 Construction and function (Fig. 1)

Upper and lower joint sections are constructed of light metal and are connected to one another by a specially designed axis. For swing phase control, a miniaturised hydraulic unit is integrated into the joint body. It serves to create dynamic extension resistances that optimise the gait pattern by preventing excessive heel rise during flexion and terminal impact during extension.

The resistance to movement (damping behaviour) can be easily changed by adjusting the cross section of flow. Knee flexion and knee extension resistances are independently adjustable.

- The arrangement of the linkage geometry allows powerful flexion resistance even with a high degree of knee flexion.
- At high walking speed, a special feature of the hydraulic unit reduces the flexion angle.
- Extension resistance may be set to a low value, since it automatically increases gradually as the knee approaches full extension.

Unlike mechanical gait controls, the hydraulic unit automatically adjusts to different walking speeds thus ensuring a harmonious pattern of movements. The flexion angle of approx. 135° with the 3R95 is limited when the modular adapter contacts the lower joint section (Fig. 1 also applies to the 3R95).

3 Technical data

Article number	3R95
Proximal connection	Pyramid adapter
Distal connection	Pyramid adapter
Knee flexion angle	135°
Weight	360 g
System height	62 mm
Proximal system height up to alignment reference point	6 mm
Distal system height up to alignment reference point	56 mm
Max. user weight	150 kg / 330 lbs
Mobility Grade	3 + 4

4 Handling

4.1 Alignment recommendations and assembly

The three-dimensional arrangement of the prosthetic socket and the modular components affects the static and dynamic functions of the prosthesis. Function of the joint is influenced by the position of the axes.

The advantages of the 3R95 can only be made optimal use of in case of a correct alignment. **The optimal residual limb position must be anticipated when positioning the socket connector.** Plumb lines in the frontal and sagittal planes (drawn from the hip joint's centre of rotation and marked during plaster cast taking and trial fitting of the test socket) will facilitate correct positioning of the lamination anchor or socket adapter.

To align the prosthesis please proceed in two steps:

1. First make the bench alignment using an alignment tool such as 743L200 L.A.S.A.R. Assembly.
2. For subsequent optimisation of the static alignment, please use 743L100 L.A.S.A.R. Posture.

4.1.1 Bench Alignment with Alignment Tool (Fig. 4)

- ① Position the middle of the foot approx. **30 mm anterior** to the alignment reference line.
- ② **Add 5 mm** to the required heel height of the foot. Set correct outward rotation of the foot.
- ③ Clamp the 3R95 knee joint using the appropriate adapter inserts. For bench alignment, **place the rotation axis of the knee joint (alignment reference point) 15 mm posterior to the alignment reference line.** At that point the pyramid base should be horizontal. Pay attention to the knee-ground distance and outward rotation of the knee (adapter inserts provide for a rotation of approx. 5°)! Recommended positioning of the alignment reference point: 20 mm above the medial tibial plateau.
- ④ Connect the foot to the modular knee joint using a tube adapter.
- ⑤ Mark the centre of the socket proximally and distally on the lateral side. Draw a line through both marks from socket brim to the distal end of the socket.
- ⑥ Position the socket such that the alignment reference line passes through the proximal centre mark of the socket. Set the socket flexion to somewhere between 3° and 5°; however, the

individual situation (e.g. hip joint contractures) must be taken into account and, if necessary, more flexion should be provided. Also pay attention to the ischial tuberosity to ground distance!

Attention!

If the residual limb flexion is not taken into account, the joint will be positioned too far to the front. This will lead to malfunction and premature wear. If necessary, use the 4R118 Adapter Plate to optimise the alignment corrections.

- ⑦ Connect the socket and modular knee joint using a corresponding adapter (e.g. 4R111, 4R41, 4R55, 4R51 Socket Adapter).

4.1.2 Static Alignment Optimisation ONLY using 743L100 L.A.S.A.R. Posture (Fig.5)

The bench alignment can be substantially improved using L.A.S.A.R. Posture. In order to ensure appropriate stability combined with easy swing phase initiation, please proceed as follows:

- ① To make the load line visible, the transfemoral amputee stands on the L.A.S.A.R. Posture with the prosthetic side on the force plate and with the other leg on the height compensation panel. The prosthesis side should be sufficiently loaded (> 35 % of the body weight).
- ② Now adapt the alignment by only **adjusting the plantar flexion** of the foot. The **load line** (laser line) should be approx. **45 mm anterior to the knee axis** (alignment reference point).
- ③ After step 2, dynamic optimisation can take place between parallel bars.

4.2 Combination possibilities

CAUTION
<p>Failure to observe the manufacturer requirements regarding combination possibilities</p> <p>Injuries, malfunctions or product damage due to unallowable combination of prosthesis components</p> <ul style="list-style-type: none"> • Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's field of application. • Contact the manufacturer with questions.

INFORMATION
<p>In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.</p>

Unallowable combinations	
Prosthetic hip joints	7E10*

4.3 Adjustment of swing phase control

To avoid damages, the prosthetist must use the 4X50 Adjustment Key that is provided attached to the knee joint to adjust swing phase control! Once the knee joint has been adjusted, the prosthetist should keep this key for future use. Resistance to flexion and extension are independently adjusted by moving the pins on each side of the piston rod. The adjustment pins are **only accessible via the slots in the housing when the knee joint is flexed.**

4.3.1 Adjustment of flexion (Fig. 2)

3R95: Flex the knee joint up to the stop (Fig.4).

Insert the provided adjustment key into the slot on the side marked with „**Flexion**“ and adjust according to the scale.

Turn the key **posteriorly** = **low** flexion resistance = **0**

Turn the key **anteriorly** = **highest** flexion resistance = **8**

Note:

Remove the adjustment key from the slot while the knee joint is still flexed.

4.3.2 Adjustment of extension (Fig. 3)

3R95: Flex the knee joint up to the stop (Fig. 3)

Insert the provided adjustment key into the slot on the side marked with „**Extension**“ and adjust according to the scale.

Turn the key **anteriorly** = **low** extension resistance = **0**

Turn the key **posteriorly** = **highest** extension resistance = **8**

Attention!

When making extension damping adjustments, make sure that the leg still completely extends with each step.

4.4 Assembly instructions

Once alignment is complete, make sure that the 506G3 Set Screws of the adapters used have the appropriate length. Set screws which are too short may not provide adequate stability for the prosthesis. Screws which are too long can damage the cosmetic foam cover.

Use the 710D4 Torque Wrench to tighten the set screws. **For the tightening torques please refer to the Instructions for Use of the corresponding adapters.**

Use a 3R24 or 3S124 Foam Cover to complete the prosthesis.

Attention!

Do not use talcum powder to eliminate noise in the cosmetic foam cover. Talcum powder reduces the lubrication of the mechanical parts, which may lead to malfunctions and thus increase the risk of failure. Using this medical product after application of talcum powder will render all claims null and void.

Notice:

To reduce friction and to eliminate noise, please apply 519L5 Silicone Spray directly on the contact surfaces of the cosmetic foam cover.

4.5 Maintenance instructions

Attention!

Do not disassemble the knee joint! If you have a problem, please send the knee joint in for service.

Ottobock recommends readjusting the knee joint's settings once the patient has gotten used to the prosthesis. This period of time varies depending on individual patient characteristics.

Please check the knee joint for wear and proper function once a year and make adjustments, if necessary. Special attention should be paid to the flexion and extension resistances and abnormal noises.

Full flexion and extension must be guaranteed.

Whenever servicing the knee joint, replace the rubber stop.

5 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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