



3R40

DE Gebrauchsanweisung.....	5
EN Instructions for use.....	13
FR Notice d'utilisation.....	20
IT Istruzioni d'uso	28
ES Instrucciones de uso	36
PT Manual de utilização	45
NL Gebruiksaanwijzing	53
SV Bruksanvisning	61
DA Brugsanvisning	68
NO Bruksanvisning	76
FI Käyttöohje	84
HU Használati utasítás.....	92
CS Návod k použití	100
RO Instrucțiuni de utilizare	108
HR Upute za uporabu	116
TR Kullanım kılavuzu.....	124
EL Οδηγίες χρήσης.....	131
ZH 使用说明书	140

● 4D16	▲	■
(D) Einzelteile-Pack	Mindestmenge	Einzelteile
(GB) Single-Component Pack	Minimum order quantity	Single components
(F) Kit de pièces de rechange	Quantité minimum	Pièces à l'unité
(I) Imballo dei singoli componenti	Quantità minima	Singoli componenti
(E) Kit componentes	Cantidad mínima	Componentes
(P) Conjunto de peças de reposo	Quantidade mínima	Componentes
(NL) Onderdelenpakket	Minimikvantitet	Onderdelen
(S) Service-Set	Minimum aantal	Separata delar
(DK) Komponentpakke	Mindste mængde	Komponenter
(PL) Zestaw naprawczy	Ilość minimalna	Podzespoły
(H) Alkatrész-csomag	Minimum mennyiség	Alkatrészek
(CZ) Sada jednotlivých dílů	Minimální množství pro objednání	Jednotlivé díly
(RO) Set componente individuale	Cantitate minimă	Componente individuale
(HR) Pakiranje pojedinačnih dijelova	Minimalna količina	Pojedinačni dijelovi
(TR) Müfritler parçaları	Minimum miktar	Müfrit parçaları
(GR) Σετ μεμονωμένων εξαρτημάτων	Ελάχιστη ποσότητα	Μεμονωμένα εξαρτήματα
(RU) Комплект деталей	Минимальное количество	Отдельные детали
(J) 交換パーセツト	最低注文個数あり	個別販売
(CN) 维修组件	最低起订量	单个部件

INFORMATION

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INFORMATION

Last update: 2020-07-08

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

1.1 Construction and function

The 3R40 prosthetic knee joint is a monocentric modular lightweight knee joint. Thanks to the design, its weight has been reduced to 290 g. The upper and lower joint sections are manufactured from a high-strength light metal alloy.

The surface is protected against corrosion and wear with a special treatment.

The locking mechanism integrated in the lower joint section secures the joint in the extended position; the lock latch engages the stop screw. The locking mechanism is unlocked with the lock bale, which is attached to the switch axis, allowing knee flexion. A moveable plastic cap covers the front joint gap across the entire flexion range (Fig. 4).

1.2 Combination possibilities

⚠ CAUTION

Failure to observe the manufacturer requirements regarding combination possibilities

Injuries, malfunctions or product damage due to unallowable combination of prosthesis components

- ▶ Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's field of application.
- ▶ Contact the manufacturer with questions.

INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

Unallowable combinations

Prosthetic hip joints	7E10*
-----------------------	-------

1.3 Service life

⚠ CAUTION

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.
- ▶ Inform the patient.

As a basic principle, all Ottobock knee joints are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

2 Application

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of application

Area of application according to the Ottobock Mobility System (MOBIS)

Recommended for amputees with **mobility grade 1** (indoor walker)
Approved for a patient weight of up to **100 kg**.

2.3 Environmental conditions

⚠ CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table „Unallowable environmental conditions“ in this section).
- ▶ If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Allowable environmental conditions

Temperature range:

Use: -10 °C – +45 °C

Storage, Transportation: -20 °C – +60 °C

Relative humidity: 20 % – 90 % (Nicht kondensierend)

Commercially available, solvent-free cleaning agents

Unallowable environmental conditions

Storage/Transportation: Mechanical vibrations, Impacts

Hygroscopic particles (e.g. Talcum powder), Dust, Sand, Fresh water, Salt water, Acids, Sweat, Urine

Cleaning agents containing solvents

3 Safety

3.1 Explanation of warning symbols

⚠ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

3.2 General safety instructions

CAUTION

Reaching into the area of the joint mechanism

Pinching of fingers, pinching of skin

- ▶ Do not reach into the joint mechanism!

3.3 Fitting safety instructions

CAUTION

Excessive strain on the product

Injuries caused by falling and damage to the product

- ▶ In order to avoid risks to the safety of the user and increase the life expectancy of the product, handle the product and its components with care.
- ▶ Only use the product if it is approved for the user's required mobility grade (see the section "Area of Application").
- ▶ Replace the product immediately if the patient's mobility grade changes.
- ▶ Inspect the product for damage if it has been subjected to extreme strain (e.g. due to falling).
- ▶ Send the product to Ottobock Service if there is any doubt about its suitability for further use.

CAUTION

Excessive mechanical strain

Injuries caused by falling and damage to the product

- ▶ Do not subject the system components to mechanical vibrations or impacts.

CAUTION

Use under environmental conditions that are not allowable

Injuries caused by falling and damage to the product

- ▶ Do not expose prosthetic components to substances which could cause metal parts to corrode, such as fresh water, salt water and acids.
- ▶ Penetration of liquids, dust, sand, etc. and use of the prosthesis in environments with excessively high relative humidity must be avoided.

4 Delivery

4.1 Scope of delivery

The 3R40 Modular Lightweight Knee Joint has been set at the factory and inspected ready for installation. The stop screw (6) is secured with Loctite®.

- | | |
|---------------------------------------|---|
| (1) 4G118 Upper Joint Section | (12) 506G3=M3x3 Set Screw (2x) |
| (2) 4G120 Lower Joint Section | (13) 513T4=18.8x10.2x0.35 Plate Spring (2x) |
| (3) 4G123 Plastic Cap | (14) 4G125 Lock Bale with Cable Guide |
| (4) 4A58 Axis Pin (2x) | (15) 21A18=2x600 Perlon Cable |
| (5) 4A56 Joint Axis | (16) 21A12 Threaded Fitting, short |
| (6) 4Z50 Stop Screw | (17) 4F17=N Lock Slide, complete |
| (7) 4G127 Lock Latch | (18) 21A8 Cable Clamp |
| (8) 4A57 Switch Axis | (19) 4G124 Plastic Ring |
| (9) 504R1=2.5x20 Semicircular Rivet | (20) 501Z2=M6x20 Cap Screw |
| (10) 513D18=4.7x31 Compression Spring | (21) 4Z59=4x6 Extension Stop (2x) |
| (11) 4Z51 Spring Bolt | |

5 Preparation for use

5.1 Alignment Instructions

CAUTION

Incorrect alignment

Injuries due to falling and damage to the prosthesis components

- ▶ Observe the alignment instructions.

CAUTION

Unallowable combination of prosthesis components

Injuries due to falling and damage to the prosthesis components

- ▶ Pay attention to the information about the possible combinations when ordering and assembling.

CAUTION

Usage of highly hygroscopic particles (e.g. talcum) and materials that remove fat

Injuries due to falling, functional disturbances and damage to the product

- ▶ Do not allow the product to come into contact with the aforementioned materials.

CAUTION

Incorrect assembly of the tube

Injuries due to falling and damage to the prosthesis components

- ▶ During assembly, slide the tube fully into the intended prosthesis component until it has reached the bottom of the clamp fitting.

The three-dimensional alignment of the prosthesis socket and the modular components affects the static and dynamic function of the prosthesis. The function of the joint is influenced by the position of the axis. The advantages of the 3R40 Modular Lightweight Knee Joint can only be fully realised with correct alignment.

The optimum residual limb position must be taken into account for the socket connection.

Plumb lines in the frontal and sagittal planes, which are marked starting at the hip pivot point during plaster casting and trial fitting, simplify correct positioning of the lamination anchor or socket adapter.

The pyramid adapter connects the knee joint proximally (to the socket). The distal connection is realised by the clamp of the lower joint section, into which the tube of the tube adapter is inserted.

Perform the alignment in 2 steps:

Bench alignment is performed first using the 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly alignment apparatus. This is followed by static alignment optimisation using the 743L100 L.A.S.A.R. Posture.

5.1.1 Bench alignment with alignment apparatus (the following steps refer to Fig. 7)

1. Move the centre of the foot forward 30 mm relative to the alignment reference line.
2. Set the effective heel height of the foot and add 5 mm.
3. Set the outward rotation of the foot.
4. Clamp the knee joint, positioning the knee centre 15 mm posterior to the alignment reference line.
5. Pay attention to the knee-ground distance and outward rotation of the knee (approx. 5° – this is provided for by an adapter insert).
6. Connect the foot to the modular knee joint using the supplied tube adapter.
 - Cut the tube to length with a 719R2 Pipe Cutter.
 - Deburr the edge of the cut with a 718R1 Deburring Knife.
 - Insert the tube into the clamp fitting all the way to the stop.
 - Screw on the cap screw using the 710D1 Torque Wrench. (Tightening torque **13 Nm**, see Fig. 1)
 - Screw on the 506G3 Tube Adapter set screws.
 - For trial fitting:** Screw in the set screws located next to each other with a 710D4 Torque Wrench. (Tightening torque **10 Nm**)
 - For final assembly:** Pre-tighten the set screws located next to each other with a 710D4 Torque Wrench. (Tightening torque **10 Nm**).
Screw in all 4 set screws with a 710D4 Torque Wrench. (Tightening torque **15 Nm**)
7. Mark the lateral centre of the socket with a centred, proximal dot and a centred, distal dot. Draw a line through both marks from socket brim to the distal end of the socket.
8. Position the socket so that the proximal centre of the socket coincides with the alignment reference line. Adjust the socket flexion to 3° – 5°, but take the individual situation (e.g. hip joint contractures) and the ischial tuberosity-to-ground distance into account.

⚠ CAUTION

Alignment error (the joint is located too far anterior)

Injuries due to falling caused by malfunctions and/or premature wear

► Note the flexion of the residual limb.

9. Connect the socket and modular knee joint using a suitable adapter (e.g. 4R111, 4R41, 4R55, 4R51 Socket Adapter).

5.1.2 Static alignment optimisation (the following steps refer to Fig. 8)

Significant optimisation of the bench alignment is possible using the L.A.S.A.R. Posture. In order to achieve sufficient safety while providing easy swing phase initiation, please proceed as follows:

1. To measure the load line, the transfemoral amputee steps on the L.A.S.A.R. Posture force plate with the prosthetic fitting side and places the other leg on the height compensation plate. Sufficient weight should be placed on the prosthesis side (> 35 % body weight).
2. Now the alignment should only be adjusted by **changing plantar flexion** so that the **load line** (laser line) runs approx. **45 mm in front of the axis** (see Fig. 8).
3. Perform dynamic optimisation during trial walking.

6 Adjustment

CAUTION

Danger of falling for the user due to incorrect alignment/forward displacement

Knee flexion moment in the stance phase

- ▶ For safety during alignment, provide enough rear displacement!

The 3R40 Modular Lightweight Knee Joint has been set at the factory and inspected ready for installation. The stop screw (6) is secured with Loctite®.

A small amount of play has no functional disadvantages!

The locking mechanism can be readjusted with the stop screw (Fig. 5). Proceed as follows:

1. Remove the plastic ring from the plastic cap (Fig. 3).
2. Unlock the mechanism and flex the knee joint (Fig. 4).
3. Lock play with 6 mm pin wrench (Fig. 5):

The stop screw may be unscrewed by up to 1/2 turn.

Turning the stop screw to the **right** = more play

Turning the stop screw to the **left** = less play

CAUTION

Alignment error (the joint is located too far anterior)

Injuries due to falling caused by malfunctions and/or premature wear

- ▶ Note the flexion of the residual limb.

6.1 Final assembly

CAUTION

Incorrect assembly of the pull cable or incorrect lock engagement

Injuries due to falling and damage to the prosthesis components

- ▶ Check the locking function after making adjustments – the lock must not open unintentionally under knee flexion moment, and must engage when fully extended (unintentional flexion of the knee joint).

The locking mechanism engages with an audible click and can be unlocked easily with the locking mechanism cable. The lock slide can be exchanged for the 4F34 Lock Lever.

The locking mechanism cable is attached to the cable guide on the lock bale and should run along the lateral prosthesis side. The cable guide can be moved on the lock bale (Fig. 2) with a suitable tool (e.g. screwdriver).

6.1.1 Foam cover

The following foam covers are suitable for the 3R40 knee joint:

- 3R24
- 3S124

7 Maintenance/repair

INFORMATION

All Ottobock knee joints are generally tested for 3 million load cycles.

Depending on the amputee's activity level, this corresponds to a service life of 3 to 5 years.

In general, we recommend regular annual safety inspections.

7.1 Maintenance

Ottobock recommends adjusting the knee joint to the user's requirements again after the user has taken some time to get used to the prosthesis.

Please check the knee joint for wear and correct functionality at least once a year and make adjustments if necessary. Special attention should be paid to movement resistance and abnormal noises. Full flexion and extension must be ensured.

7.1.1 Spare parts

4D16 Single Component Pack for 3R40 (Fig. 1, •)

consisting of:

- 4G123 Plastic Cap
- 506G3=M3×3 Set Screw (2x)
- 4G125 Lock Bale with Cable Guide
- 21A18=2×600 Perlon Cable
- 21A12 Threaded Fitting, short
- 4F17=N Lock Slide, complete
- 21A8 Cable Clamp
- 4G124 Plastic Ring
- 4Z59=4x6 Extension Stop (2x)

7.2 Disposal/re-use

7.2.1 Disposal Information

INFORMATION

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

7.2.2 Information on re-use

The products are intended for exclusive use by one user. The manufacturer does not authorise use of the products by another person.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

9 Appendices

9.1 Technical data

Article Number	3R40
Proximal connection	Pyramid adapter
Distal connection	Tube clamp
Weight	290 g
Knee flexion angle	155°
System height	23 mm
Proximal system height up to the alignment reference point	1 mm
Distal system height up to the alignment reference point	22 mm
Max. user weight	100 kg
Mobility grade	1

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