


## 5R2=C

<b>DE</b>	Gebrauchsanweisung .....	5
<b>EN</b>	Instructions for use .....	9
<b>FR</b>	Instructions d'utilisation.....	13
<b>IT</b>	Istruzioni per l'uso .....	17
<b>ES</b>	Instrucciones de uso .....	21
<b>PT</b>	Manual de utilização.....	25
<b>NL</b>	Gebruiksaanwijzing.....	29
<b>SV</b>	Bruksanvisning .....	33
<b>DA</b>	Brugsanvisning .....	37
<b>NO</b>	Bruksanvisning .....	41
<b>FI</b>	Käyttöohje.....	44
<b>PL</b>	Instrukcja użytkownika.....	48
<b>HU</b>	Használati utasítás .....	52

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### INFORMATION

Date of last update: 2021-08-06

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

### 1.1 Construction and Function

The 5R2=C Carbon Socket Attachment Block is integrated into the prosthetic socket of a transfemoral or transtibial prosthesis using prepreg technology. It serves to connect the prosthetic socket to a socket adapter.

### 1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

## 2 Intended use

### 2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

### 2.2 Area of application

Approved for a body weight of **up to 150 kg**.

### 2.3 Environmental conditions

#### Storage and transport

Temperature range  $-20\text{ °C}$  to  $+60\text{ °C}$  ( $-4\text{ °F}$  to  $+140\text{ °F}$ ), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

#### Allowable environmental conditions

**Temperature range:**  $-10\text{ °C}$  to  $+45\text{ °C}$  ( $14\text{ °F}$  to  $113\text{ °F}$ )

**Chemicals/liquids:** fresh water, soapsuds, chlorine water

**Moisture:** Submersion: max. 1 h in 2 m depth, relative humidity: no restrictions

**Solids:** Dust, occasional contact with sand

**Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage** (see page 12).

#### Prohibited environmental conditions

**Solids:** highly hygroscopic particles (e.g. talcum), dust in high concentrations (e.g. construction site), intensive contact with sand

**Chemicals/liquids:** salt water, perspiration, urine, acids, continuous use in liquid media

### 2.4 Lifetime

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

## 3 Safety

### 3.1 Explanation of warning symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

#### ⚠ CAUTION!

#### Risk of injury and risk of product damage

- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 9).
- ▶ Note the combination possibilities/combination exclusions in the instructions for use of the products.
- ▶ Observe the maximum lifetime of the product.
- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

#### NOTICE!

#### Risk of product damage and limited functionality

- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

#### Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

### 4 Scope of delivery

Fig.	Item	Quantity	Designation	Reference number
–	–	1	Instructions for use	–
1	1	1	Socket attachment block, carbon	5R2=C
1	2	1	Lamination dummy	4X301
<b>The following parts are only for:</b>			<b>Connection to socket adapter: 4R22, 4R23, 4R54, 4R55, 4R73=A, 4R73=D, 4R74, 4R95</b>	
1	3	4	Countersunk head screw	501S41=M6x12- _12.9
<b>The following parts are only for:</b>			<b>Connection to socket adapter: 4R37, 4R51, 4R77</b>	
1	4	4	Countersunk head screw	501S41=M6x16

### 5 Preparing the product for use

#### ⚠ CAUTION!

#### Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

#### INFORMATION

The layup described in this document was approved for the maximum product user body weight. The prosthetist assumes full responsibility for any change to the layup.

## 5.1 Fabricating the prosthetic socket

### > Required tools and materials:

99B81 PVA Bag, 616B10 Carbon Fibre Woven Prepreg, 616B16 Peel Ply, Woven Absorbent Fabric, 636K8 Plastaband, 634A58 Isopropyl Alcohol, 633W8 Wax, 616B11 Carbon Fibre Prepreg, 4X301 Lamination Dummy

- 1) Soak a PVA bag and pull it over the model.
- 2) Apply a layer of carbon fibre woven prepreg to the entire model.
- 3) Apply three layers of carbon fibre woven prepreg (e.g. 15 cm x 15 cm) offset to each other on the distal socket end.
- 4) Apply a layer of peel ply to the entire model.
- 5) Apply a layer of woven absorbent fabric to the entire model.
- 6) Soak a PVA bag and pull it over the model.
- 7) Bake the model in the oven under vacuum.
- 8) Remove the PVA bag, absorbent fabric and peel ply.
- 9) Seal the threads on the underside of the carbon socket attachment block with Plastaband (see fig. 2).
- 10) Position the socket attachment block on the socket end as required for proper alignment and secure it with a mix of sealing resin and talcum.
- 11) Sand the hardened mix of sealing resin and talcum to shape the socket (see fig. 3).
- 12) Clean the prosthetic socket with a degreasing cleaning agent and remove dust.
- 13) Apply isolation wax to the threads at the top of the socket attachment block (see fig. 4).
- 14) Set the lamination dummy on the socket attachment block (see fig. 5).
- 15) Fill the edges between the socket attachment block and the lamination dummy with carbon fibre prepreg (see fig. 6).
- 16) Apply a layer of carbon fibre woven prepreg to the entire model and the socket attachment block (see fig. 7).
- 17) Apply a thread of carbon fibre prepreg in a circle around the base of the socket attachment block and pull it tight (see fig. 8).

- 18) Apply a layer of carbon fibre woven prepreg (width: **5 cm**) in a circle around the base of the socket attachment block so that the ends overlap by **4 cm** (see fig. 9).
- 19) Apply medial and lateral reinforcement of carbon fibre prepreg to the prosthetic socket according to the weight and mobility grade of the patient.
- 20) For TF prostheses, reinforce the prosthetic socket with carbon fibre prepreg on the ischial tuberosity according to the weight and mobility grade of the patient.
- 21) Apply two layers of carbon fibre woven prepreg to the entire model and the socket attachment block (see fig. 10).
- 22) Apply a layer of peel ply to the entire model (see fig. 11).
- 23) Apply a layer of woven absorbent fabric to the entire model.
- 24) Soak a PVA bag and pull it over the model.
- 25) Bake the model in the oven under vacuum.

## 5.2 Final assembly

### CAUTION

#### Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

#### Selection chart for countersunk head screws

Chosen socket adapter	Screws to be used
4R22, 4R23, 4R54, 4R55, 4R73=A, 4R73=D, 4R74, 4R95	501S41=M6x12_12.9
4R37, 4R51, 4R77	501S41=M6x16

> **Required materials:** sanding tool, 636K13 Loctite, torque wrench (e.g. 710D20)

1) **CAUTION! Do not sand the laminate around the edge of the carbon adapter!**

Sand the laminate on the distal end of the prosthetic socket down to the surface of the lamination dummy.

2) Remove the lamination dummy (see fig. 12).

3) Select the appropriate countersunk head screws to connect the socket adapter to the socket attachment block (see selection chart).

4) Apply Loctite to the threads of the socket attachment block.

5) **CAUTION! Observe the minimum screw insertion depth of 8 mm!**

Install the socket adapter according to its instructions for use.

## 6 Cleaning

1) Rinse the product with clear fresh water.

2) Dry the product with a soft cloth.

3) Allow to air dry in order to remove residual moisture.

## 7 Maintenance

► A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.

► Inspect the entire prosthesis for wear during normal consultations.

► Conduct annual safety inspections.

## 8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

### 9.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

## 10 Technical data

Reference number	5R2=C
Weight [g]	50
System height [mm]	10
Build height [mm]	10
Material	Carbon
Max. body weight [kg]	150

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