








## 4R152\* Harmony E2

<b>DE</b> Gebrauchsanweisung .....	6
<b>EN</b> Instructions for use .....	15
<b>FR</b> Instructions d'utilisation .....	24
<b>IT</b> Istruzioni per l'uso .....	34
<b>ES</b> Instrucciones de uso .....	44
<b>PT</b> Manual de utilização .....	53
<b>NL</b> Gebruiksaanwijzing .....	63
<b>SV</b> Bruksanvisning .....	72
<b>DA</b> Brugsanvisning .....	81
<b>NO</b> Bruksanvisning .....	91
<b>FI</b> Käyttöohje .....	99
<b>TR</b> Kullanma talimatı .....	108
<b>RU</b> Руководство по применению .....	117
<b>JA</b> 取扱説明書 .....	128

## INFORMATION

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### INFORMATION

Date of last update: 2020-07-16

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

### 1.1 Construction and Function

The 4R152 Harmony E2 electronic vacuum pump is a component of the Harmony system for the generation of a high vacuum in the prosthetic socket. The vacuum pump and a prosthetic socket fabricated using Harmony technology form a functional unit for transfemoral and transtibial prostheses. Adhesion between the prosthesis and the residual limb is improved by the high vacuum.

The vacuum pump is installed on the 4R153 adapter plate or the 4R153=1 adapter plate with offset. The adapter plate with offset can be used when a greater lateral clearance to the prosthetic knee joint or prosthetic foot is required.

## 1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

## 2 Intended use

### 2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

### 2.2 Area of application

- The maximum approved body weight is specified in the technical data (see Page 24).

### 2.3 Environmental conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +60°C
Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e.g. talcum)

## 2.4 Service Life

### Vacuum Pump

The Harmony E2 vacuum pump was tested by the manufacturer with 400 hours of pump operation. The vacuum pump starts up at irregular intervals that depend on the selected mode, patient activity level and leak tightness of the overall system, including the prosthetic socket. Accordingly the service life of the device is 3 to 5 years.

### Adapter Plate

In general, all modular adapters are tested by the manufacturer for 3 million load cycles. Depending on the patient's level of activity, this corresponds to a service life of 3 to 5 years.

## 3 Safety

### 3.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 3.2 General safety instructions

 <b>WARNING</b>
<b>Use of the vacuum pump in explosion-prone rooms</b>
Explosion hazard due to sparks
▶ Never use the vacuum pump in explosion-prone rooms.
▶ Remove the vacuum pump from the prosthesis prior to entering explosion-prone rooms. Do not take the vacuum pump along with you into the explosion-prone rooms.

**⚠ CAUTION**

**Excessive strain on the product**

Risk of injury due to breakage of load-bearing components

- ▶ Use the product according to the specified area of application (see Page 16).

**⚠ CAUTION**

**Unallowable combination of prosthetic components**

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

**⚠ CAUTION**

**Use under unallowable environmental conditions**

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

**⚠ CAUTION**

**Exceeding the service life and reuse on another patient**

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded.
- ▶ Only use the product on a single patient.

**⚠ CAUTION**

**Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

**Signs of changes in or loss of functionality during use**

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

**4 Scope of delivery**

4R152* Harmony E2 (see fig. 1)				
Fig.	Item	Quantity	Designation	Reference number
1	1	1	Harmony E2 vacuum pump	–
2	–	1	Adapter plate set	4R153 or 4R153=1

4R152* Harmony E2 (see fig. 1)				
Fig.	Item	Quantity	Designation	Reference number
1	2	1	Hose piece for socket bore	–
1	3	1	Hose piece for outlet	4X220
1	4	1	Adapter plug (US)	757S7=US
1	5	1	Adapter plug (EU)	757S7=EU
1	6	1	Power supply	757L2
1	7	1	USB charging plug	4X269
1	8	1	O-ring set for suction connection of the vacuum pump	4X315
1	9	1	Torx T8 screwdriver	–
1	10	2	Battery, rechargeable, size AA	625B13

4R153* Adapter plate (see fig. 2)			
Item	Quantity	Designation	Reference number
1	1	Adapter plate (with or without offset)	–
2	1	Torx offset screwdriver	709S51=T8
<b>Single-component pack, consisting of:</b>			4X305
3	1	PE solid rod (PE sinter filter)	–
4	1	TPU insert for 5R2 (thin)	–
		TPU insert for TEC plate (thick)	–
5	1	Alternative connection for valve plate	–
6	1	O-ring	–
7	1	Locking screw	–
8	4	M6x38 mm screw	–
9	4	M6x35 mm screw	–

Spare part	Compatible product	Reference number
<b>O-ring set (3 pieces) for suction connection of the vacuum pump</b>	Vacuum pump from serial number: SN201714001	4X315
	Vacuum pump with older serial number	4X316

## 5 Preparing the product for use

### CAUTION

#### **Incorrect alignment or assembly**

Risk of injury due to damaged prosthetic components

- Observe the alignment and assembly instructions.

### CAUTION

#### **Use of powder (e.g. baby powder or talcum)**

Skin irritation on the residual limb and loss of function of prosthetic components due to clogging with particles or deprivation of lubricant

- Keep the product away from powders.

### 5.1 Preparing the Vacuum Pump

### CAUTION

#### **Contact with defective electrical components**

Risk of injury due to electric shock

- ▶ Conduct a visual inspection before touching electrical components.
- ▶ Do not use damaged electrical components.

### 5.1.1 Inserting or Changing Batteries

The vacuum pump is operated with 2 rechargeable NiMH batteries, size AA. Operation with the supplied batteries is recommended by the manufacturer. Alternatively the product can be operated with non-rechargeable batteries.

- > **Required materials:** Torx T8 screwdriver, 2 batteries (size AA), torque wrench adjustable to 0.5 Nm
- > **Prerequisite:** Battery charger is disconnected
- 1) Pull out the offset screwdriver on the side of the vacuum pump (see fig. 3).
- 2) Pull the vacuum pump off the adapter plate and loosen the outlet hose (see fig. 4, see fig. 5).
- 3) Loosen the four screws on the battery compartment (see fig. 6) and remove the battery compartment cover.
- 4) Remove the old batteries if necessary.
- 5) Insert two batteries into the compartment, making sure to observe the correct polarity (see fig. 7).
- 6) Reinstall the battery compartment cover with the screws (tightening torque: **0.5 Nm**).
- 7) Connect the outlet hose and attach the vacuum pump to the adapter plate.
- 8) Insert the offset screwdriver into the opening on the side of the vacuum pump.

### 5.1.2 Charging the Batteries

#### NOTICE

#### Charging non-chargeable batteries

Damage to the product due to battery leakage

- ▶ Only connect the product to the battery charger when you are using rechargeable batteries.

#### INFORMATION

Rechargeable batteries lose capacity over the course of repeated charging cycles. This is a normal process.

- ▶ Replace the batteries when the operating time of the product decreases noticeably.

Depending on the performance of the batteries, the vacuum level and patient activity, the battery capacity of the vacuum pump is 2 to 3 days. The manufacturer recommends charging the batteries daily. A complete charging cycle takes 4 to 6 hours. Remove the vacuum pump from the product or take off the prosthesis during the charging process. The battery charge level is indicated by signals shown in the table.

#### Assembling the Battery Charger

- 1) Connect the chosen adapter plug to the power supply (see fig. 8).
  - 2) Ensure that the adapter plug engages in the power supply (see fig. 9).
  - 3) Connect the USB plug of the charging cable to the power supply (see fig. 10).
- The battery charger is assembled.

#### Charging the batteries

- > Only charge the batteries with the battery compartment cover in place.
- 1) Set the selector switch to standby mode.
- 2) **While wearing the prosthesis:** Pull the vacuum pump off the adapter plate. Remove the outlet hose from the vacuum pump if necessary.
- 3) Open the charging receptacle cover.
- 4) Insert the battery charger plug into the charging receptacle (see fig. 11); a signal indicates that the charging process has started.

- 5) After the charging process is complete, disconnect the plug from the charging receptacle and close the charging receptacle cover.
- 6) Disconnect the battery charger plug from the outlet.
- 7) **If necessary:** Connect the outlet hose and attach the vacuum pump to the adapter plate.

Light Signal	Audible Signal	Status
LEDs 1 and 2 flash alternating yellow (until automatic shutdown)	Signal tone ringing (1 x)	Batteries almost drained, remaining capacity $\leq 1$ h
LEDs 1 and 2 flash alternating red (until automatic shutdown)	Signal tone ringing (3 x)	Batteries drained, device shuts down
LEDs 1 and 2 flash alternating green	–	Charger plugged in, connection established
LEDs 1 and 2 flash alternating green, brightness dimmed	–	Batteries are being charged
LEDs 1 and 2 illuminated green, brightness dimmed	–	Batteries fully charged
LEDs 1 and 2 illuminated turquoise, brightness dimmed	–	Batteries fully charged, reduced battery performance; we recommend replacing the batteries

## 5.2 Installing the Adapter Plate

The adapter plate for connecting the vacuum pump is installed between the prosthetic socket and the distal components of the modular prosthesis. It is mounted together with the socket adapter. The adapter plate has two connectors for generating a vacuum in the socket and one connector for connecting to the vacuum pump. The adapter plate is additionally equipped with a PE sinter filter to protect the vacuum pump against dirt, and a non-return valve that maintains the vacuum when the pump is disconnected.

### CAUTION

#### Using a wrong prosthetic socket

Contusions and pressure sores due to wrong socket fit

- ▶ Use the vacuum pump only in combination with a prosthetic socket made with Harmony technology.

### 5.2.1 Preparing the Adapter Plate

A vacuum can be generated in the socket either via the proximal or the lateral connector of the adapter plate.

### INFORMATION

Using the proximal connector of the adapter plate is recommended only when using a 5R2 lamination disc or the 4A94=3 TEC plate (with TPU insert respectively). The lateral connector for the intake hose can be recognised by the internal thread for the hose connector.

#### Using the proximally located connector

### INFORMATION

To ensure an optimum sealed connection to the socket, the following approach is recommended:

- ▶ Drill the distal socket bore with a **6 mm** bit.
- ▶ Apply PU adhesive inside the bore and insert the included hose piece (see fig. 12).
- ▶ Remove excess adhesive and allow the adhesive to set.
- ▶ Trim the hose piece flush on both ends.

> **Required tools and materials:**

Torque wrench, M12x1.5 locking screw, TPU insert (thin) for 5R2 or TPU insert (thick) for TEC plate, PE solid rod (PE sinter filter)

- 1) Screw the locking screw into the lateral connector of the adapter plate and tighten (mounting torque: **3 Nm**) (see fig. 13).
- 2) Press the PE solid rod into the TPU insert to the stop (see fig. 14).
- 3) Slide the TPU insert with the PE solid rod into the connector on the proximal side of the adapter plate (see fig. 15).

**Using the lateral connector**

> **Required tools and materials:**

Torque wrench, 2 mm Allen wrench, M12x1.5 locking screw, PE solid rod (PE sinter filter), O-ring for hose connector, hose connector

- 1) Screw the locking screw into the proximal connector of the adapter plate and tighten (mounting torque: **3 Nm**) (see fig. 16).
- 2) Press the PE solid rod into the lateral connector of the adapter plate up to the stop (see fig. 16).
- 3) Insert the O-ring into the groove of the hose connector.
- 4) Use the Allen wrench to screw the hose connector into the lateral connector and tighten (installation torque: **1 Nm**) (see fig. 16).

**5.2.2 Installing the Adapter Plate**

**CAUTION**

**Improper assembly of the screw connections**

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

**INFORMATION**

Note the position of the Harmony E2 vacuum pump in the prosthesis. The vacuum pump must not bump against the prosthesis or project beyond the flexed prosthetic knee joint. Use the 4R153=1 adapter plate with offset if you require greater lateral clearance (see fig. 23, see fig. 24).

**Selection table for countersunk head screws**

Chosen socket adapter	Screws to be used
4R22, 4R23, 4R54, 4R55, 4R73=A, 4R73=D, 4R74, 4R95	501S128=M6X35
4R37, 4R51, 4R77	501S128=M6X38

> **Required tools and materials:**

710D4 Torque Wrench, 636K13 Loctite® 241

- 1) Position the adapter plate on the lamination disc with the desired alignment (see fig. 17).
- 2) Set the chosen socket adapter onto the adapter plate (see fig. 18).
- 3) **INFORMATION: Secure the screws with Loctite®.**  
Insert the screws through the socket adapter and the adapter plate into the lamination disc and, in diagonal sequence, tighten them using a torque wrench (mounting torque: **12 Nm**) (see fig. 19).

## 5.3 Installing the Hoses

### 5.3.1 Optional: Installing the Intake Hose

An intake hose only needs to be installed if the lateral connector of the adapter plate will be used.

- 1) Connect the intake hose to the lateral connector of the adapter plate.
- 2) Run the intake hose to the intended connector of the prosthetic socket and connect it.

### 5.3.2 Installing the Outlet Hose

The suctioned air and fluids are drained through the outlet hose. Fluids must not be drained onto metal components of the modular prosthesis. Increasing the length of the outlet hose lowers the noise level of the vacuum pump.

#### CAUTION

##### Improper installation of the outlet hose

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Run the outlet hose in such a way that the hose opening is not directed towards parts of the body or into the cosmetic cover.

- 1) Connect the outlet hose to the corresponding connector on the vacuum pump (see fig. 5).
- 2) **If a cosmetic foam cover is used:** Make an opening and route the outlet hose through the cosmetic foam cover to the outside.
- 3) **If no cosmetic cover is used:** Route the outlet hose along the prosthesis to the prosthetic foot and secure it.
- 4) Trim any excess outlet hose.

## 6 Usage

### 6.1 Installing the Vacuum Pump

Slide the vacuum pump onto the adapter plate. The connector for generating the vacuum is plugged into the adapter plate. The vacuum pump is secured with the Torx offset screwdriver to prevent it from being inadvertently pulled out.

- 1) **If the outlet hose is attached to the prosthesis:** Connect the vacuum pump with the outlet hose.
- 2) Set the vacuum pump onto the guide rails of the adapter plate diagonally from above, and slide it on horizontally to the stop (see fig. 20).
- 3) Insert the offset screwdriver into the opening on the side of the vacuum pump (see fig. 21).

### 6.2 Operating the Vacuum Pump

The vacuum pump offers several modes for selection by the user. It is operated via a selector switch on the housing. There is an automatic mode, 4 constant vacuum modes and a reverse mode. In automatic mode, the activity level is determined based on the pressure fluctuations within the system. The vacuum pump adjusts the vacuum dynamically in order to ensure optimum adhesion and comfort. In reverse mode, the vacuum pump conveys air into the socket. This function serves as the flushing mode and can be used to assist removal. The various modes are listed in the table:

Symbol on select-or switch	Selected mode
	Automatic mode
	Constant vacuum, level 1 to 4, indicated by circles of increasing size
	
	
	
	Reverse mode

Symbol on select- or switch	Selected mode
	Standby mode

► **NOTICE! Only use the vacuum pump when it is installed on the adapter plate.**

Use the selector switch to choose the desired mode. The vacuum pump reports certain operating states with signals.

Light signal	Audible signal	Status
LEDs 1 and 2 flash white	Signal tone rising (1x)	Restart of the device
LEDs 1 and 2 illuminated red	Signal tone falling (1x)	Device switching to standby mode
–	Click sound	Selector switch operated, mode change
LEDs 1 and 2 flash blue	–	Pump running (forwards or in reverse)
LEDs 1 and 2 flash green, every 10 s		Device is powered up and working without leaks
LEDs 1 and 2 flash alternating blue and red, every 10 s	–	Leak warning (see information)
LEDs 1 and 2 flash red	Signal tone, like siren (1x)	System error

### Leak warning information

The leak warning indicates that the desired vacuum cannot be achieved or maintained. The vacuum pump continues running. The leak warning is output under the following conditions:

- The vacuum pump operates continuously for more than 3 minutes
- Less than 2 minutes elapse between several pump intervals

The leak warning deactivates itself automatically when there is no problem for 10 minutes. It can also be deactivated by briefly activating the vacuum pump standby mode.

### 6.3 Flushing the Vacuum Pump

A reduced vacuum indicates contamination. In case of contamination, the system has to be flushed to clean it. This flushes the vacuum pump and the filter in the adapter plate. The reverse mode of the vacuum pump is used for flushing. The flushing process can also help if the system is plugged.

#### CAUTION

#### Improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- Clean the product regularly.

#### > Required materials:

Vessel with 30 ml to 60 ml distilled water, one collecting vessel

#### Preparations required:

Vacuum pump installed on adapter plate and in standby mode, prosthesis taken off the leg

- 1) Dip the outlet hose into the vessel with distilled water.
- 2) Hold the socket in such a way that the water runs into the collecting vessel during flushing.
- 3) Set the selector switch to reverse mode.
- 4) Allow the vacuum pump to run until the water has been pumped through the entire system.
- 5) Allow the pump to dry run for 15 s in order to pump out any remaining water.
- 6) Set the selector switch to standby mode.
- 7) Clean the interior surface of the socket.

## 7 Maintenance

- A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.

- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

## 8 Disposal

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

## 10 Technical data

Reference number	4R152* Harmony E2	4R153 adapter plate	4R153=1 adapter plate with offset
Weight [g]	185	125	160
System height [mm]	95	22	22
Material	–	Aluminium	Aluminium
Max. body weight [kg]	–	150	150
Battery charger operating voltage [V]	100 to 240	–	–
Battery charger operating frequency [Hz]	50 to 60	–	–
Battery charging temperature [°C]	0 to 40	–	–
Battery charging temperature [°F]	32 to 104	–	–

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