

6A41

DE	Gebrauchsanweisung	3
EN	Instructions for use	7
FR	Instructions d'utilisation	12
IT	Istruzioni per l'uso	17
ES	Instrucciones de uso	21
PT	Manual de utilização	26
NL	Gebruiksaanwijzing	30
SV	Bruksanvisning	35
DA	Brugsanvisning	40
NO	Bruksanvisning	44
FI	Käyttöohje	49
PL	Instrukcja użytkowania	53
HU	Használati utasítás	58
CS	Návod k použití	62
RO	Instrucțiuni de utilizare	67
HR	Upute za uporabu	71
SL	Navodila za uporabo	76
SK	Návod na používanie	80
BG	Инструкция за употреба	85
TR	Kullanma talimatı	90
EL	Οδηγίες χρήσης	94
RU	Руководство по применению	99
JA	取扱説明書	104
ZH	使用说明书	108
KO	사용 설명서	112

INFORMATION

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INFORMATION

Date of last update: 2020-04-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.

- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 6A41 Slider Plate is mounted at the MagnoFlex Lock and serves to optimise the prosthetic alignment. It allows for an offset of the distal prosthetic components by 12 mm or 24 mm in an a - p or m - l direction.

1.2 Combination Possibilities

6A41 Slider Plate	
Designation	Reference Number
Shuttle Lock	6A40
Socket Adapter	4R54

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

Approved for a body weight of up to **125 kg**.

2.3 Environmental conditions

Storage and transport
Temperature range -20 °C to +60 °C (-4 °F to +140 °F), relative humidity 20 % to 90 %, no mechanical vibrations or impacts
Allowable environmental conditions
Temperature range: -10 °C to +45 °C
Moisture: relative humidity: 20% to 90%, non-condensing
Unacceptable environmental conditions
Chemicals/liquids: fresh water, salt water, perspiration, urine, acids, soapsuds, chlorine water
Solids: dust, sand, highly hygroscopic particles (e. g. talcum)

2.4 Service life

In general, all modular adapters are tested by the manufacturer for 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum service life of 5 years.

3 Safety

3.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the prosthetic components according to the MOBIS classification.

CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

CAUTION

Exceeding the service life and reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded.
- ▶ Only use the product on a single patient.

CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.

- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

NOTICE

Use under unallowable environmental conditions

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop, etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of Delivery

6A41 Slider Plate					
Fig.		Item No.	Quantity	Designation	Reference Number
-	■	-	1	Instructions for Use	647G755
1	-	1	1	Slider Plate	-
1	▲	2	4	Countersunk Head Screws	501S41=M6X50

5 Preparation for Use

CAUTION

Incorrect alignment or assembly

Injuries due to fall resulting from damaged prosthesis components

- ▶ Observe the alignment and assembly instructions.

CAUTION

Dirty thread in screw connection

Injuries due to falling when thread loosens, breaks or tears out

- ▶ Use a degreasing cleaner (e.g. Acetone) to clean the thread before each mounting.
- ▶ Observe the alignment and assembly instructions.

> **Required tools and materials:**

710D4 Torque Wrench, 636K13 Loctite® 241

- 1) **If a socket adapter is installed at the shuttle lock:** remove the socket adapter.
- 2) Apply Loctite® to the thread.
- 3) Position the slider plate at the shuttle lock and screw it in using the countersunk head screws (see fig. 2, see fig. 3) (mounting torque: **12 Nm**).
- 4) Apply Loctite® to the thread of the countersunk head screws.
- 5) Position the socket adapter at the slider plate and then use the countersunk head screws to fasten it (see fig. 4, see fig. 5) (mounting torque: **12 Nm**).

6 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

7.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

8 Technical Data

Reference Number	6A41
Weight [g]	105
System height [mm]	9
Material	Aluminium
Max. body weight [kg]	125

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