
## 13E100

<b>DE</b> Gebrauchsanweisung (Fachpersonal) .....	3
<b>EN</b> Instructions for use (qualified personnel) .....	8
<b>FR</b> Instructions d'utilisation (Personnel spécialisé) .....	14
<b>ES</b> Instrucciones de uso (Personal técnico especializado) .....	19
<b>PT</b> Manual de utilização (Pessoal técnico) .....	25
<b>NL</b> Gebruiksaanwijzing (Vakmensen) .....	30
<b>SV</b> Bruksanvisning (Fackpersonal) .....	36
<b>NO</b> Bruksanvisning (Fagpersonell) .....	41
<b>PL</b> Instrukcja użytkowania (Personel fachowy) .....	47
<b>CS</b> Návod k použití (Odborný personál) .....	52
<b>RO</b> Instrucțiuni de utilizare (Personal de specialitate) .....	58

**INFORMATION**

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### INFORMATION

Date of last update: 2021-03-29

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The AnalogAdapter will be referred to simply as the product/adapter below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Function

The adapter allows for the analogue Ottobock hand systems to be controlled in combination with the DynamicArm. The digital signal of the DynamicArm is converted to an analogue signal by the adapter.

### 2.2 Combination possibilities

The product is intended exclusively for use in MyoBock components.

Prosthetic components that are especially well suited for combination with the product are listed below.

Name	Reference number
DynamicArm/ DynamicArm Plus	12K100N=*/ 12K110N=*
SensorHand Speed	8E38=8-*/ 8E39=8-*/ 8E41=8-*
MyoHand VariPlus Speed	8E38=9-*/ 8E39=9-*/ 8E41=9-*
System Electric Hand DMC plus	8E38=6-*/ 8E39=6-*/ 8E41=6-*
System Electric Hand Digital Twin	8E38=7-*/ 8E39=7-*/ 8E41=7-*
DMC VariPlus System Electric Greifer	8E33=9/ 8E34=9
Digital Twin System Electric Greifer	8E33=7/ 8E34=7
Electric Rotator	10S17

## 3 Intended use

### 3.1 Indications for use

The product is to be used solely for upper limb prosthetic fittings.

### 3.2 Conditions of use

The prosthesis is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

The product has been developed for use in everyday activities. It must not be used for activities where it would be exposed to unallowable environmental conditions or strains (e.g. extreme sports).

### 3.3 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### 3.4 Qualification

The fitting of a patient with the product may only be carried out by O&P professionals who have been authorised with the corresponding Ottobock training.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

## 4.2 Structure of the safety instructions

### WARNING

#### The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

## 4.3 General safety instructions

### WARNING

#### Non-observance of safety notices

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

### CAUTION

#### Damage to connecting cable due to kinking or small radii

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ Replace the damaged connecting cable immediately.

### CAUTION

#### Use of unsuitable prosthetic components

Injury due to unexpected product behaviour.

- ▶ Use the product only in combination with components listed in the section "Combination Possibilities" (see page 9).

### NOTICE

#### Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).

### CAUTION

#### Improper use

Injury due to faulty operation or malfunction of the product.

- ▶ Instruct the patient in the proper use of the product.

## 4.4 Patient information

### INFORMATION

Advise your patient of the information in this section.

### WARNING

#### Operating the product near active implanted systems

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the product.

- ▶ When operating the product in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

**⚠ CAUTION**

**Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm to HF communication devices is recommended.

**⚠ CAUTION**

**Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

**⚠ CAUTION**

**Remaining in areas outside the allowable temperature range**

Injury due to faulty control or malfunction of the product.

- ▶ Avoid remaining in areas outside the allowable temperature range (see page 13).

**⚠ CAUTION**

**Penetration of product with dirt and humidity**

Injury due to unexpected product behaviour or malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

**⚠ CAUTION**

**Independent manipulation of the product**

Injury due to product damage or malfunction.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

## 5 Scope of delivery

- 1 pc 13E100 AnalogAdapter
- 2 pc hook and loop strips
- 1 pc. Instructions for use (qualified personnel)

## The following components are also required:

- 633F11 Silicone Grease

## 6 Preparing the product for use

### 6.1 Preparation of the Adapter

- 1) Attach a hook and loop strip to the domed side of the adapter (see fig. 2).
- 2) Attach the second hook and loop strip to the inside of the DynamicArm.

**INFORMATION: The position of the hook and loop strip inside the DynamicArm can be selected freely.**

### 6.2 Wiring the Adapter

#### INFORMATION

Plugs marked with a colour dot may be connected only to ports with a colour dot!

#### 6.2.1 Connecting Wiring between the DynamicArm and Adapter

Use the two cables with the three-pin plugs for the connection between the DynamicArm and the adapter.

- 1) Lubricate the cable connectors with silicone grease (see fig. 3).
- 2) Connect the three-pin plug to the port on the adapter.

The colour-coded plug must be connected to the colour-coded port (see fig. 4).

#### 6.2.2 Connecting the Wiring between the Adapter and Electric Rotator

- 1) Attach the motor cable of the DynamicArm to one of the two no. 3 contacts (see fig. 6).
- 2) Connect the three-pin adapter plug to the no. 2 contacts.

The colour-coded plug must be connected to the colour-coded contact.

#### 6.2.3 Connecting the Wiring between the Adapter and Coaxial Plug

- 1) Connect the three-pin adapter plug to the no. 2 contacts of the coaxial plug (see fig. 7).
- 2) The motor cable must not be attached to the coaxial plug.

### 6.3 Finishing the Adapter

- 1) Press the adapter with the hook and loop strip onto the hook and loop strip inside the DynamicArm (see fig. 8).
- 2) Thread the connecting cable into the DynamicArm.

## 7 Cleaning and Care

- 1) Clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock Derma Clean) when needed.

Make sure that no liquids get into the product and product components.

- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

## 8 Legal information

### 8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 8.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 8.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

## 9 Technical data

<b>Ambient conditions</b>	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing

<b>General Information</b>	
Reference number	13E100
Product service life	5 years
Weight	approx. 7 g
Voltage	7.6 VDC
Tightness	IP22 (EN 60529)

## 10 Appendices

### 10.1 Symbols Used

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

Lot number (PPPP YYYY WW)  
PPPP – plant  
YYYY – year of manufacture  
WW – week of manufacture

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