
4R101

DE	Gebrauchsanweisung	3
EN	Instructions for Use.....	6
FR	Instructions d'utilisation	9
IT	Istruzioni per l'uso	13
ES	Instrucciones de uso	16
PT	Manual de utilização.....	19
NL	Gebruiksaanwijzing.....	23
SV	Bruksanvisning.....	26
DA	Brugsanvisning.....	29
PL	Instrukcja użytkowania.....	33
HU	Használati utasítás	36
EL	Οδηγίες χρήσης	40
RU	Руководство по применению	43

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English

Date of the last update: 2014-01-15

- Please read this document carefully.
- Follow the safety instructions.

1 Application

The 4R101 Sliding Adapter is to be used **exclusively** for the prosthetic fitting of lower limbs.

Field of application according to the **Otto Bock Mobility System MOBIS®**:

Approved for a patient with a **weight of up to 100 kg / 220 lbs.**

The Otto Bock MOBIS® is not intended to be used as a guide to obtaining reimbursement for prosthetic components in the USA.

The 4R101 Sliding Adapter is designed for permanent integration in the transfemoral prosthesis.

Attention!

In a transtibial prosthesis the adapter is to be used only in the form of an interim solution.

2 Description

The 4R101 Sliding Adapter increases the adjustment capabilities for modular lower limb prostheses. The adapter's aluminum upper and lower sec-

tions may be moved independently in two planes via a cross support. The Sliding Adapter is installed between 5R1 Socket Attachment Block or 5R6 Thermoplastic Adapter and any socket adapter (fig. 1). The socket material between the shells of the 5R6 does not affect adjustability when using the 4R101 Modular Adapter.

The 4R101's low structural height facilitates its use for relatively long residual limbs. A 4R37/4R51 Rotatable Socket Adapter may be used to add rotation adjustments.

Attention!

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water, acids and other liquids. Using this medical device under such conditions will render all claims against Otto Bock HealthCare null and void.

Please also inform your patients!

3 Technical Data

Article number	4R101
Weight	205 g
System height	25 mm
Material	Aluminum
Max. patient weight	100 kg
Displacement / m – l and a – p directions	± 11 mm
Mobility level	–

4 Handling

4.1 Assembly and Adjustment

Use the 501S41=M6×25 Countersunk Head Screws supplied with the adapter for its connection to the socket attachment block. The **501S41=M6×16 Countersunk Head Screws** are intended for the attachment of a socket adapter to the lower section of the adapter.

Attention!

Only the 501S41=M6x16 Countersunk Head Screws provided may be used for screwing the 4R101 Sliding Adapter to the 4R37/4R51/4R77 Socket Adapters.

When mounting the 4 countersunk head screws, tighten the screws in the following order:

1. Tighten both countersunk head screws on the displacement side with controlled torque to **12 Nm**.
2. Afterwards, tighten the two countersunk head screws on the opposite side with controlled torque **12 Nm**.

Make sure that the screws are screwed to minimum thread depth of **8 mm** to avoid stress fractures!

Using a 710D1 Torque Wrench.

For adjustment loosen the center screw by **maximal turning of 1–1½**, using a 4 mm wrench. Sliding is independently adjustable in the frontal plane (medial/lateral) and in the sagittal plane (anterior/posterior). The slide range of ± 11 mm is indicated on the scale (fig. 2). Using a 2 mm wrench, secure the desired adjustment by tightening the two clamping screws (fig. 3). The four clamping screws and the center screw must also be tightened (fig. 4) to secure the upper and lower sections on the dovetail guides of the cross support. For extended periods of trial walking or if the adapter is to remain in the definitive prosthesis, use a 710D1 Torque Wrench or equal to tighten the center screw to a **torque moment of 7 Nm (5.2 ft./lb. or 62 in./lb.)**.

Attention!

Do not use the Sliding Adapter without center screw.

4.2 Maintenance Instructions

Note:

As a basic principle, all Otto Bock modular adapters are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

5 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components in accordance with the Otto Bock Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

6 Warranty

Otto Bock warrants this product/these products from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction. Because warranty conditions may differ by country, please contact your local Otto Bock distribution company for more information.

7 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH