


# ottobock.

4R108=3L, 4R109=3L, 4R108=3T, 4R109=3T

|           |                                  |    |
|-----------|----------------------------------|----|
| <b>DE</b> | Gebrauchsanweisung .....         | 4  |
| <b>EN</b> | Instructions for Use.....        | 10 |
| <b>FR</b> | Instructions d'utilisation ..... | 16 |
| <b>IT</b> | Istruzioni per l'uso .....       | 22 |
| <b>ES</b> | Instrucciones de uso .....       | 28 |
| <b>PT</b> | Manual de utilização.....        | 34 |
| <b>NL</b> | Gebruiksaanwijzing.....          | 40 |
| <b>SV</b> | Bruksanvisning .....             | 46 |
| <b>DA</b> | Brugsanvisning.....              | 51 |
| <b>PL</b> | Instrukcja użytkowania.....      | 57 |
| <b>CS</b> | Návod k použití.....             | 63 |
| <b>EL</b> | Οδηγίες χρήσης .....             | 68 |

**INFORMATION**

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**INFORMATION**

Last update: 2017-03-21

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

## 1 Description and function

### 1.1 Indications for use

The Ottobock 4R108 and 4R109 Endoskeletal Alignment Adapters are to be used **solely** for fitting lower limb prostheses, especially transtibial applications. Each socket adapter functions as an adjustable connection between the socket and distal component section.

### 1.2 Introduction

The 4R108 and 4R109 Endoskeletal Alignment Adapters offer a broad range of quick and easy linear, angular, and rotary adjustments for thermoplastic or laminated prosthetic sockets. The adapters are designed to simplify socket fabrication and alignment, particularly for initial prostheses and remote fittings. The 4R108 incorporates a 30 mm pylon clamp, while the 4R109 has a fixed pyramid adapter.

Approved for a patient weight of up to 100 kg/220 lbs.

### 1.3 Environmental conditions

#### **Allowable environmental conditions:**

- ▶ Temperature range for use: -10 °C to +60 °C
- ▶ Allowable relative humidity 0 % to 90 %, non-condensing

#### **Unallowable environmental conditions:**

- ▶ Mechanical vibrations or impacts
- ▶ Perspiration, urine, fresh water, salt water, acids
- ▶ Dust, sand, highly hygroscopic particles (e.g. talcum)

## 1.4 Service Life

In general, all modular adapters are tested by the manufacturer for 3 million load cycles. Depending on the patient's level of activity, this corresponds to a service life of 3 to 5 years.

## 1.5 Features (fig. 1)

### NOTICE

**Do not bend or scribe the components.**

### Components for 4R108=3L and 3T (*italicized*)

- (1) Adapter with 2Y39 Tube Connector
- (2) 2Y38=3L Lamination Shell
- (3) 2Y38=3T Thermoplastic Shell
- (4) 2Y33=3 Adapter Plate
- (5) 501S41=M10×20 Countersunk Head Screw
- (6) 501S41=M10×25 Countersunk Head Screw
- (7) 4X100 Lamination Dummy
- (8) 501Z2=M6×25 Cylinder Head Screw
- (9) 506S1=5×16 Clamping Sleeve (2 pcs.)

### Components for 4R109=3L and 3T (*italicized*)

- (10) 2Y41 Adapter with Pyramid
- (11) 2Y38=3L Lamination Shell
- (12) 2Y38=3T Thermoplastic Shell
- (13) 2Y33=3 Adapter Plate
- (14) 501Z10 Oval Countersunk Head Screw
- (15) 501Z21=M10×45 Eyebolt DIN 580
- (16) 4X100 Lamination Dummy
- (17) 506S1=5×16 Clamping Sleeve (2 pcs.)

The 4R108 and 4R109 Adapters are mounted on the exterior of the socket or frame using the single central screw provided. In most cases, simply moving the adapter and plate to the new location will improve the overall alignment and will not require compensatory alignment changes. The adapter allows spherical sliding adjustments, i.e. in the frontal plane (abduction/adduction) by  $\pm 8$  mm and in the sagittal plane (flexion/extension) by  $\pm 6$ mm.

For definitive use, all attachments must be sealed with one drop of threadlock

and torqued to specifications. If increased loading is expected, the socket must be reinforced appropriately.

## 2 Safety

### 2.1 Explanation of warning symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 2.2 General safety instructions

#### **CAUTION**

##### **Excessive strain on the product**

Risk of injury due to breakage of load-bearing components

- ▶ Use the product according to the specified area of application.

#### **CAUTION**

##### **Unallowable combination of prosthetic components**

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

#### **CAUTION**

##### **Use under unallowable environmental conditions**

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

### **⚠ CAUTION**

#### **Exceeding the service life**

Risk of injury due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded.

### **⚠ CAUTION**

#### **Incorrect alignment or assembly**

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

### **⚠ CAUTION**

#### **Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see „Signs of changes in or loss of functionality during use“ in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

### **Signs of changes in or loss of functionality during use**

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

## **3 Alignment and connection**

### **3.1 Lamination socket (use of 4R108=3L or 4R109=3L)**

The 2Y38=3L Lamination Shell is integrated in the laminate as a lamination anchor.

Fabricate the socket liner with Pedilin<sup>®</sup>, providing a distal area build-up of

80 mm Ø, then pull PVA sheeting material over the socket liner. Insert the lamination dummy into the lamination shell and cover the adjustment area with PE adhesive tape. Next, place the lamination shell so that the longest sliding path is in the medial-lateral direction. The shell's marks should face anterior. (fig. 2).

Reinforcement is provided with two layers of Nylglass 623T9, partial application of carbon layers (fig. 3), 1 ½ layers of Fiberglass Stockinette 616G13 and four layers of Nylglass 623T9. Each layer is tied into the groove of the adapter shell (fig. 4). Finally, remove the lamination dummy from the socket. After the model has been removed from the mold, assemble the adapter plate and adapter with the countersunk screw(s).

### **3.2 Thermoplastic socket** (use of 4R108=3T or 4R109=3T)

The 2Y38=3T thermoplastic shell is integrated in the thermoplastic socket so that the plastic material is located between the thermoplastic shell and the adapter. When using a socket liner, this is built up to 80 mm Ø in the distal area (refer to point 2.1).

Pull the Nylon Stockinette 99B25 over the socket, then add the thermoplastic shell with the lamination dummy inserted.

Use plastilin to blend away any uneven cuts or edges, so that disassembling the thermoplastic shell is possible (fig. 5).

After vacuum forming, drill through the indentation in the lamination dummy using a 4 mm bit. Remove the model and lamination dummy and widen the hole to 26 mm.

Prior to assembly, make sure the thermoplastic material is evenly distributed over the thermoplastic shell area (min. 8 mm). Smooth out any irregularities on the attachment surface by grinding the thermoplastic material. This also provides a large surface area for the adapter.

#### **NOTICE**

If the outer edge of the adapter is not properly supported, there may be a **risk of fracturing** the central countersunk screw.

## **4 Adjustment and final assembly**

Linear, angular and rotary changes are accomplished by loosening the single clamping screw in the center of the assembly. After loosening the

central countersunk screw with a 6 mm pin wrench, the extension / flexion, abduction / adduction and rotation can be adjusted (figs. 6 and 7). After the trial fitting has been finished, bore holes in the lower section of the adapter through the bottom of the socket using a 4 mm bit. Insert the rolled pins to prevent rotary movement.

The 4R109 Adapter with Pyramid allows additional angular changes at any time.

#### **4.1 Important torque information**

For definitive use, the screws must be secured with Loctite 636K13 or equivalent, and the **countersunk screw** should be **torqued to 20 Nm** (14.7 ft. /lb. 177 in. /lb.)

Use the 701D1 Torque Wrench or equivalent to tighten the **clamping screw** at the tube connector with a torque of **13 Nm**. (9.6 ft. /lb. 115 in. /lb.)

## **5 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **5.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### **5.2 CE conformity**

This product meets the requirements of the European Directive 93 / 42 / EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

### **5.3 Warranty**

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

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Otto Bock HealthCare GmbH