














# 743A23

<b>DE</b> Gebrauchsanweisung .....	6
<b>EN</b> Instructions for Use.....	8
<b>FR</b> Instructions d'utilisation.....	10
<b>IT</b> Istruzioni per l'uso .....	12
<b>ES</b> Instrucciones de uso.....	15

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English

**Last update: 2013-07-01**

- **Please read this document carefully.**
- **Follow the safety instructions and the stated precautions in this accompanying document.**
- **If assistance is required during the start-up, use or maintenance of the product, as well as in the event of unexpected operating behaviour or circumstances, contact the manufacturer (see manufacturer's address on the back of this accompanying document).**

## 1 Intended Use

The 743A23 Alignment Tool for ErgoArm Elbow Components is to be used **exclusively** in upper limb prosthetics for alignment and foaming as well as for support of functional training and for rehabilitation purposes.

## 2 Safety Instructions

**Failure to follow the safety instructions mentioned below can lead to faulty control or malfunction of the pressure switch and result in risk of injuries for the patient.**

- **Before processing the Pedilen Rigid Foam 200 (Art. no. 617H12), please refer to the corresponding material safety data sheet!**
- **Before processing the 617H19 Orthocryl or 617H55 C-Orthocryl Lamination Resin and 617P37 Hardener read the corresponding safety data sheet!**
- **For your own protection wear appropriate respirator, protective gloves and goggles when sanding the socket.**
- **Ensure sufficient ventilation of working rooms and spaces!**

## 3 Delivery Scope


Fig. 1 1 Pc. Alignment Tool

Fig. 2 1 Pc. Foaming Cover

Fig. 3 2 Pcs. Twist-lock tie

1 Pc. 647H447 Manual

---

## 4 Description and Function

The 743A23 Alignment Tool for ErgoArm Elbow Components is used for adjusting an ErgoArm to the interim socket of the patient and for further alignment of the prosthesis in preparation for the lamination of the definitive above-elbow socket. Wearing the finished interim socket during the rehabilitation phase helps the patient train his movement sequences and so allows for a seamless change to the "real" prosthesis.

## 5 Range of Application

The 743A23 Alignment Tool for ErgoArm Elbow Components can be used for the 12K41, 12K42, 12K44 and 12K50 ErgoArm product family.

## 6 Preparatory Operations

- ① Tear protective foil off double-sided adhesive tape only in the lower area.
- ② Fix alignment aid to lamination ring by using the tape and the twist-lock tie.
- ③ Tighten twist-lock tie with pliers. Grab twist-lock tie with pliers and roll it over the twist-lock tie fastener. Cut off protruding end.

## 7 Working with the Amputee

- ④ Fix the MYOBOCK electrodes with 13E135 electrode holder to interim socket.
- ⑤ Lift off the four protective foils.
- ⑥ Adjust ErgoArm to provisional socket. **(Attention: The socket must be free of grease!)**
- ⑦ Use Lasar measuring device to check symmetry on patient.
- ⑧ Additional fixation of ErgoArm with adhesive tape to achieve stability.
- ⑨ Mark alignment positions on interim socket.
- ⑩ Finish alignment with electronic components.
- ⑪ Check ergonomics and electronic functions.
- ⑫ Transfer resp. examination of position possible with 743A100=3 Alignment Apparatus.

## 8 Preparation and Fabrication of the Definitive Socket

- ⑬ Place lamination cover on lamination ring (Photo without electrode holes)
- ⑭ Mark hole positions on plaster cast.
- ⑮ Remove interim socket from cast
- ⑯ Place corrected model into vacuum forming device
- ⑰ Insulate plaster model with 99B25 Sock
- ⑱ Form definitive socket with 616T53 Thermolin.soft

## 9 Working on the Definitive Inner Socket

- ⑲ Grind off the material covering the 13E127 casting dummy
- ⑳ Remove dummy with screw and pliers
- ㉑ Insert 13E128 casting dummy and cover definitive socket with 99B71 PVC-Plastic sheeting
- ㉒ Coat lamination cover with 633F11 Silicone Grease
- ㉓ Place alignment aid on marked positions on the definitive socket

## 10 Fabricating the Foam Core for the Outer Socket

- ㉔ Foam cover with space for twist-lock tie.
- ㉕ Place foam cover on lamination ring. The foam cover must be flush with the lamination ring.
- ㉖ Attach 627B5=19 Foam Adhesive Tape on to definitive socket
- ㉗ Wrap 616T3=1 Trolen Sheeting conically around definitive socket and secure with adhesive tape. Shape Trolen cone so that you can fill in the Pedilen Rigid Foam 200.
- ㉘ Let Pedilen Rigid Foam climb up to 3/4. Then place model horizontally and leave it to cure to rule out cavities.
- ㉙ After curing pull off Pedilen sheeting
- ㉚ Pull foam from socket

- 
- ③1 Cut through foam until reaching foam cover.
  - ③2 Remove foam cover from lamination ring.
  - ③3 Grind foam cone down to alignment aid.
  - ③4 Cut through alignment aid fixation arms.
  - ③5 Cut through twist-lock tie with oscillating saw.
  - ③6 Remove alignment aid from lamination ring.
  - ③7 Lift off lamination ring.
  - ③8 Finish the foam core.
  - ③9 The fixation arms remain in the foam core. For lamination procedures please refer to the 647H437 Instructions for Use of the 12K42/12K44/12K50 ErgoArm family or to the 647H438 Instructions for Use of the 12K41 ErgoArm.

## 11 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

## 12 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise.

The product may only be repaired by authorised Ottobock technicians.

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