1G9, 1H32, 1H34, 1H38, 1H40

DE Gebrauchsanweisung	3
EN Instructions for use	10
FR Instructions d'utilisation	17
IT Istruzioni per l'uso	25
ES Instrucciones de uso	33
PT Manual de utilização	41
NL Gebruiksaanwijzing	48
SV Bruksanvisning	56
DA Brugsanvisning	63
NO Bruksanvisning	71
FI Käyttöohje	78
PL Instrukcja użytkowania	85
HU Használati utasítás	93
CS Návod k použití	100
RO Instrucțiuni de utilizare	108
HR Upute za uporabu	115
SL Navodila za uporabo	123
SK Návod na používanie	130
BG Инструкция за употреба	137
TR Kullanma talimatı	145
EL Οδηγίες χρήσης	152
RU Руководство по применению	160
JA 取扱説明書	169
ZH 使用说明书	176
KO 사용 설명서	182

INFORMATION

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INFORMATION

Last update: 2016-10-25

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The 1G9, 1H32, 1H34, 1H38 and 1H40 single axis prosthetic feet are suitable for use in both modular and exoskeletal prostheses. The functional properties are achieved through the combination of a contoured core and functional foam.

1.2 Combination possibilities

This prosthetic component is part of the Ottobock modular system and can be combined with other products of this modular system.

Modular foot adapter	
1G9	2R51*
1H32, 1H34, 1H38, 1H40	2R10*, 2R33*, 2R51*

2 Intended use

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

The product is suitable for fittings of geriatric patients.

2.2 Area of application

The product is recommended for mobility grade 1 (indoor walker).

The product must be used in TF prostheses only.

Prosthetic foot	Max. body weight [kg]
1G9	75
1H32, 1H34, 1H38, 1H40	100

2.3 Environmental conditions

Allowable environmental conditions	
Temperature range for use: -10 °C to +60°C	
Allowable relative humidity 0 % to 90 %, non-condensing	
Unallowable environmental conditions	
Mechanical vibrations or impacts	
Perspiration, urine, fresh water, salt water, acids	
Dust, sand, highly hygroscopic particles (e.g. talcum)	
Other environmental conditions	
UV light	Not resistant
Hydrolysis	Not resistant

2.4 Service Life

This product has been tested by the manufacturer with 2 million load cycles. Depending on the patient's activity level, this corresponds to a service life of 2 to 3 years.

3 Safety

3.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

3.2 General safety instructions

 CAUTION
Excessive strain on the product
Risk of injury due to breakage of load-bearing components

- ▶ Use the prosthetic components according to the MOBIS classification (see section "Area of application").

⚠ CAUTION

Exceeding the service life and reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded.
- ▶ Only use the product on a single patient.

⚠ CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

NOTICE

Use under unallowable environmental conditions

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions (see the section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

NOTICE

Mechanical damage to the product

Change in or loss of functionality due to damage

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).

- Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover characteristics) is an indication of loss of functionality.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Prosthetic foot	–
1	Foam connecting cap (included with 1G9)	2R63

1G9 Spare parts/accessories (not included in the scope of delivery)	
Designation	Reference number
Single component pack for single axis feet	2D5
Single axis foot adapter with screw connection (aluminium)	2R51*

1H32, 1H34, 1H38, 1H40 Spare parts/accessories (not included in the scope of delivery)	
Designation	Reference number
Single component pack for single axis feet	2D5
Dorsal stop set	2S88*
Connection cap	2R22*
Single axis foot adapter with screw connection (titanium)	2R33*
Single axis foot adapter with screw connection (aluminium)	2R51*
Single axis foot adapter with screw connection (steel)	2R10*
Pedilan sole	
for 1H32	2Z67*
for 1H34	2Z64*

5 Preparation for use

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- Observe the alignment and assembly instructions.

5.1 Installing the foot adapter

- Install the foot adapter on the prosthetic foot according to its instructions for use.

5.2 Alignment

NOTICE

Grinding the prosthetic foot

Premature wear resulting from damage to the prosthetic foot

- ▶ Do not grind the prosthetic foot.

5.2.1 Bench Alignment

TF bench alignment

- ▶ Observe the information in the prosthetic knee joint instructions for use.

5.2.2 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.
- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) may be requested from Ottobock.

5.2.3 Dynamic Trial Fitting

- Adapt the setup of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) so that correct heel contact and optimal roll-over are ensured.
- Remove the adapter cover from the prosthetic foot after completion of the dynamic trial fitting and gait training exercises.

5.3 Optional: Installing the foam cover

A connecting element (such as a connection plate, connection cap or foam connection cap) provides a detachable connection between the foam cover and the prosthetic foot.

- > **Required materials:** Degreasing cleaner (e.g. 634A58 isopropyl alcohol), 636N9 contact adhesive or 636W17 plastic adhesive
- 1) Measure the distance from the knee rotation point to the desired end of the foam cover and add the compression allowance: For PE foam **10 mm**, for PUR soft foam **30 mm**. For TF prostheses, add twice the compression allowance.
 - 2) Cut the pre-shaped foam block to length.
 - 3) Pull the foam cover over the prosthesis.
 - 4) Set the connecting element onto the footshell or prosthetic foot. Depending on the version, the connecting element engages in the edge or sits on the foot adapter.
 - 5) Mount the prosthetic foot on the prosthesis.

- 6) Mark the outer contour of the connecting element on the distal face of the foam cover.
- 7) Dismount the prosthetic foot and remove the connecting element.
- 8) Clean the connecting element using a degreasing cleaner.
- 9) Glue the connecting element onto the distal face of the foam cover according to the marked outer contour.
- 10) Let the glue dry (approx. **10 minutes**).
- 11) Mount the prosthetic foot and adapt the exterior cosmetic cover shape. Take into account compression due to cosmetic stockings or SuperSkin.

6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

10 Technical data

1G9					
Sizes [cm]	23	24	25	26	27
Heel height [mm]	10 ±5				
System height with 2R51 [mm]	63	64	65	67	69
Product weight without adapter [g]	270	280	295	315	325
Max. body weight [kg]	75				
Mobility grade	1				

1H32								
Sizes [cm]	22	23	24	25	26	27	28	29
Heel height [mm]	25 ±5							
System height with 2R33/2R10 [mm]	45		46	47	49	51	53	
System height with 2R51 [mm]	49		50	51	53	55	57	
Product weight without adapter [g]	255	325	360	400	410	440	470	
Max. body weight [kg]	100							
Mobility grade	1							

1H34							
Sizes [cm]	22	23	24	25	26	27	
Heel height [mm]	35 ±5						
System height with 2R33/2R10 [mm]	43			44	45	46	
System height with 2R51 [mm]	47			48	49	50	
Product weight without adapter [g]	300	315	335	355	380	420	
Max. body weight [kg]	100						
Mobility grade	1						

1H38								
Sizes [cm]	21	22	23	24	25	26	27	28
Heel height [mm]	10 ±5							

1H38								
Sizes [cm]	21	22	23	24	25	26	27	28
System height with 2R33/2R10 [mm]	42	43	44	45	46			47
System height with 2R51 [mm]	46	47	48	49	50			51
Product weight without adapter [g]	255	275	305	335	360	365	420	435
Max. body weight [kg]	100							
Mobility grade	1							

1H40								
Sizes [cm]	22	23	24	25	26	27	28	29
Heel height [mm]	25 ±5							
System height with 2R33/2R10 [mm]	45			46	47	49	51	53
System height with 2R51 [mm]	49			50	51	53	55	57
Product weight without adapter [g]	295	305	320	370	400	440	470	530
Max. body weight [kg]	100							
Mobility grade	1							

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Otto Bock HealthCare GmbH