



743R6

DE Gebrauchsanleitung (Benutzer / Techniker)	4
EN Instructions for Use (User / Technician)	8
FR Mode d'emploi (Utilisateur / Technicien)	12
IT Istruzioni d'uso (utilizzatore / tecnici)	16
ES Instrucciones de uso (usuario / técnico)	20
PT Manual de instruções (Usuário / Técnico)	24
NL Gebruiksaanwijzing (gebruiker / technicus)	28
SV Bruksanvisning (användare / tekniker)	32
DA Brugsvejledning (bruger / tekniker)	36
NO Bruksanvisning (bruker / tekniker)	40
PL Instrukcja użytkowania (Użytkownik / Technik)	44
CS Návod k použití (pro uživatele / techniky)	48
HU Használati utasítás (használóknak / műszerészeknek)	52
TR Kullanım kılavuzu (Kullanıcı / Teknisyen)	56
EL Οδηγίες χρήσης (χρήστης / τεχνικός)	60
RU Руководство по применению (для пользователей / технического персонала)	64

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Date of the last update: 2013-09-25

- Please read this document carefully.
- Follow the safety instructions.

1 Introduction

The 743R6 Joint Alignment Fixture was developed as a tool for aligning lower limb orthotic joints. It is suitable for all orthotic applications, including thermoset lamination, vacuumassisted thermoforming, and fabrication from a schema/delineation.

To accommodate changes in fabrication technique, this alignment fixture was designed to be as universally applicable as possible. The large variety of Ottobock orthotic joints requires several adapters. To minimize the number of adapters, certain joint types were intentionally omitted.

1.1 Components of the Orthotic Joint Alignment Fixture (ill. 1)

- (1) 743Y55 Plaster Impression Adjustment
- (2) 743Y70=80 Schema/Delineation Alignment Tube
743Y70=50 Schema/Delineation Alignment Tube (accessory)
- (3) 709S10=2.5 Allen wrench
- (4) 743Y56=1 Alignment Adapter, size 1, with 501T28=6x8 Cap Screw
- (5) 743Y56=2 Alignment Adapter, size 2, with 501A1=12x7xM4 Shoulder Screw
- (6) 743Y56=3 Alignment Adapter, size 3, with 501A1=14x9xM6 Shoulder Screw
- (7) 743Y56=4 Alignment Adapter, size 4, with threaded retainer M4 and M6 and 501A12=3 Shoulder Screw
- (8) 743Y72 Allen Head Screws with bushing for alignment adapter size 1
- (9) 501A1=12x6xM4 Shoulder Screws for alignment adapter size 2
- (10) 501A1=12x8xM6 Shoulder Screws for alignment adapter size 3

1.2 Using Various Joints with the Alignment Adapters

The following is a list of the various adapters for the leg and system joints which may be used with the alignment set.

Alignment Adapter, size 1

with 743Y61 Retainer Bolt and
501T28=M6x8 Cap Screw
for knee orthoses with polycentric joints
–stainless steel–

17K43

17K45

17K46

17K47

17K49=6

Alignment Adapter, size 1

with Allen Head Screw with bushing

for bilateral knee orthosis joints
–plastic–

7U45

7U49

7U56

Alignment Adapter, size 2

with bushing flanges 743Y64=12×7.6 and 501A1=12×6×M4 Shoulder Screw

for System Knee Joints:

17B95=16, =B16

17B96=16

for Leg Orthoses (ankle joint):

17F46=6, =5

17F34=6, =5

17F47=6, =5

17F63=6, =5, =4

17F53=6, =5

17F64=6, =5, =4

17F65=6, =5

for Leg Orthoses (knee joint):

17K32=6, =5

17K33=6, =5

17K34=6, =5

17K29=6, =5

17K42=6, =5

17K49=6

Alignment Adapter, size 2

with bushing flanges 743Y64=12×7.6 and 501A1=12×7×M4 Shoulder Screw

for System Ankle Joints:

17B62=16

17B57=16

17B54=16

17B63=16

17B59=16

17B53=16

for Leg Orthoses (ankle joint):

17F31=26x2.5

17F32=26x2.5

17F24=4, =2

Alignment Adapter, size 3

with bushing flanges 743Y64=14×9.6 and 501A1=14×9×M6 Shoulder Screw

for System Ankle Joints:

17B62=20

17B57=20

17B54=20

17B63=20

17B59=20

17B53=20

for System Knee Joints:

17B23=16, =20

17B95=20

17B45=16, =20

17B96=20

17B92=16, =20

17B33=16, =20

17B91=16, =20

17B20=16, =20

17B21=16, =20

for Leg Orthoses (ankle joints):

17F31=28x3, =30x3

17F32=28x3, =30x3

for Leg Orthoses (knee joints):

17K32=4

17K33=4

17K34=4

17K29=4

17K42=4

17K49=5, =4

Alignment Adapter, size 4

with 743Y65=M4 Threaded Retainer

for System Knee Joints:

17B26=16, =20

17B47=16, =20

17B3 =16, =20

17B43=16, =20

17B46=16, =20

17B71=16, =20

Alignment Adapter, size 4

with Threaded Retainer 743Y65=M6

for System Knee Joint:

17B44=16, =20

1.3 Accessories for 743R6 (ill. 7)

(1) Stainless Steel Locating Bolt

For 743Y56 Alignment Adapter, to be used with System Knee Joints

Article no.	For bar width	For
743Y47=1	16 mm	17B36, 17B48, 17B90, 17B17, 17B16
743Y47=2	20 mm	17B36, 17B48, 17B90, 17B17, 17B16

(2) 743Y46 Stainless Steel Positioning Aid

Positioning aid for aligning 743Y56 Alignment Adapters.

Threaded width adjustment

(3) 743Y48 Adapter Sleeve

Adapter sleeve of rigid foam for receiving the 743Y56 Alignment Adapters in the plaster positive. The adapter sleeves are slid onto the 743Y55 Plaster Impression Axis and positioned in the plaster negative

29 cm length.

1 package includes 4 pieces

2 Directions for Assembly

2.1 Alignment Adapter, size 1

To attach polycentric joints to the adapter, remove the joint screws and threaded bushings. Use the Allen head screws provided to attach the joint assemblies to the adapter.

2.2 Alignment Adapters, sizes 2 and 3

To attach the joints, remove the joint screws and threaded bushings. Use the shoulder screws provided to attach the joint to the adapter. For assembly of the Joints please use the following Shoulder Screws:

17B95/96=20	→	501A1=12×8×M6
17B105=20/17B106=20	→	501A1=12×8×M6
17B105=16/17B106=16	→	501A1=10×7×M4

Be certain the oval flange on the adapter fits properly in the medial recess of the joint.

2.3 Adjustment Adapter, size 4

To attach the joint with joint screw (17B44 System Knee Joint), screw the joint with the protruding threaded connector onto the M6 threaded retainer and then align it using the rotation adjustment (see Section 5, Rotation Adjustment). The 501A12=3 shoulder screw is used to install the system knee joints with bearing race screw (17B26, 17B47, 17B3, 17B43, 17B46 and 17B71). When positioning the 743Y55 alignment axis in the plaster negative, observe the posterior offset or pivot point of the double joints with toothed segment.

3 Orthosis Alignment using Plaster Model

This method requires the use of the alignment axis. Use a wire and the 743A6 orthosis alignment tool to define the axis of rotation (Figure 2) when installing the alignment axis in the plaster negative. To insert the axis into the plaster negative, make an X-shaped cut or 16 mm bore at the joint pivot point and push the greased alignment axis through the plaster negative (Figure 3). Fill the negative to create a plaster positive model. Once the plaster positive has hardened, pull

the alignment axis out of the model. This will leave a square channel for the alignment adapters. The joints attached to the alignment adapters can now be positioned on the plaster model with correct axis orientation and aligned correspondingly. Should a correction of the joint bar in the longitudinal direction take place, i.e. in the sagittal plane, the rotation adjustment of the adjustment adapters is to be used.

4 Orthosis Alignment using Schema/Delineations

For this procedure, the Alignment Tube is used to hold the Alignment Adapters in a congruent position (ill. 5). The clamping screws allow precise adjustment of the mediolateral width between the joint surfaces.

5 Rotation Adjustment

Every adapter has a rotation adjustment that allows precise positioning in the sagittal plane. Rotation is altered by loosening the Allen head screws (2.5 mm). This adjustment can be used to insure that the joint extension stops are congruent and that the uprights are optimally placed on the plaster positive model (ill. 4).

For Alignment Adapters sizes 1–3, the neutral position for rotation is marked by an arrow etched into the Adapter which should initially be in line with the Allen head screw. When Alignment Adapter size 4 is used, rotation must always be adjusted to compensate for the varying thread lengths which protrude medially through the joint assemblies.

6 Maintenance

Although the Alignment Adapters (ill. 6) are made from stainless steel, repeated contact with plaster may result in corrosion. Clean the adapters and oil them periodically to prevent rust. The rotation adjustment is delivered with a long-lasting lubricant. Should rotation become difficult after repeated use and exposure to moisture, the Adapter may be disassembled. Remove the 2.5 mm 506G3=M5×5 Allen head screw (3) and use a 3 mm pin punch to drive out the 743Y71 Attachment Peg (4) to allow removal of the attachment piece (2) from the square axis (1). After cleaning all parts thoroughly, lubricate and reassemble.

7 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

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Otto Bock HealthCare GmbH