
ottobock.

757Z103=* / 757Z149

DE	Gebrauchsanweisung	3
EN	Instructions for Use.....	8
FR	Instructions d'utilisation	13
IT	Istruzioni per l'uso	18
ES	Instrucciones de uso	23
PT	Manual de utilização.....	28
PL	Instrukcja użytkowania.....	34
CS	Návod k použití.....	39
EL	Οδηγίες χρήσης	44

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Date of the last update: 2014-02-26

- **Please read this document carefully.**
- **Follow the safety instructions.**

1 Intended Use

1.1 Medical Purpose

The Ottobock Battery Mounting Set 757Z103=* / 757Z149 is to be used **solely** for the exoprosthetic fitting and mounting into the arm prosthesis.

1.2 Application

The Ottobock Battery Mounting Set 757Z103=* / 757Z149 may be used only in combination with components that were explicitly authorised for the use in this combination. It is therefore indispensable in this connection to observe the instructions for use of the components or to ask your Ottobock Myo-Service.

1.3 Qualifications

The fitting of a patient with the Ottobock Battery Mounting Set 757Z103=* / 757Z149 may only be carried out by a prosthetist.

2 Safety Instructions

- Wear appropriate personal protective equipment while laminating and mounting.
- Use the Patient Information section of this pamphlet (Section 3) to inform your patient about the function and use of the battery receptacle.

3 Patient Information

- Attention should be paid that no foreign particles can get into the battery receptacle. The battery receptacle should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- For faultless functioning, the electric contacts of the battery connection cable in the battery receptacle must be clean and free of grease. The contacts can be cleaned with a dry cotton swabs.

4 Components

757Z103=1 Battery Mounting Set

for X-Change Pack 757B15

1 pc. 757Z104=1 Battery Receptacle

1 pc. 757Z105=1 Cable Retention Piece with
627F1=3.1×1.6 O-Ring

1 pc. 757Z106=1 Lamination Pattern

1 pc. 757Z107=1 Spacer Block

757Z103=2 Battery Mounting Set

for X-Change Pack 757B15

The battery is by 2.5 mm inserted into device.

1 pc. 757Z104=2 Battery Receptacle

1 pc. 757Z105=2 Cable Retention Piece with
627F1=3.1×1.6 O-Ring

1 pc. 757Z106=2 Lamination Pattern

1 pc. 757Z107=2 Spacer Block

757Z149 Battery Mounting Set

for Interchangeable Battery 757B13

1 pc. 757Z150 Battery Receptacle

1 pc. 757Z151 Cable Retention Piece with
627F1=3.1×1.6 O-Ring

1 pc. 757Z152 Lamination Pattern

1 pc. 757Z153 Spacer Block

5 Integration into the Lamination Socket

Manufacture and foam the inner socket using a comfort technology and grind the foam.

Insert the lamination pattern into the hard foam so that the profile edge protrudes 2 mm at least.

Attention!

The surface under the lamination pattern must be plane.

Use Plastilin for modelling the shape around the lamination pattern. Be careful not to damage the lamination pattern edges (Fig. 1).

Pull over the PVA bag and laminate the outer socket as customary.

Before processing the C-Orthocryl Lamination Resin 617H55 or Orthocryl Lamination Resin 617H19 and Hardener 617P37 read the corresponding safety data sheet and observe the specified processing instructions.

Once the lamination resin has hardened, grind flush the visible profile edge of the lamination pattern (Fig. 2).

Screw the wood screw or Fuse Puller 710Z2 into the lamination pattern and pull out the latter. Continue dismantling. Insert the battery receptacle so that the cable output is located distally from residual limb.

Insert the spacer block into the battery receptacle and mount into the lamination socket so that the cable retention piece is positioned proximally. Subsequently glue the battery receptacle using Orthocryl Sealing Resin 617H21.

Attention!

Be careful not to glue the cable retention piece when applying the sealing resin in the surrounding area.

Once the sealing resin has hardened, remove the spacer block and mount the Battery Connection Cable 13E51=2, 13E51=3, 13E51=4 or 13E132 and Cable Retention Piece 757Z105 / 757Z151 in the battery receptacle.

6 Liability

Otto Bock Healthcare GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer does not assume liability for damage caused by component combinations which it did not authorize.

The product may only be opened and repaired by authorized Ottobock technicians.

7 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines.

The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH