743A6

| | | |
|----|--|----|
| DE | Gebrauchsanleitung (Techniker/Benutzer) | 2 |
| EN | Instructions for use (technician/user) | 12 |
| FR | Mode d'emploi (Technicien/Utilisateur) | 22 |
| ES | Instrucciones de uso (técnicos/usuarios) | 32 |

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Table of contents EN

| | |
|--|-----------|
| 1 Introduction..... | 13 |
| 2 Area of application..... | 13 |
| 3 Components..... | 13 |
| 4 Mounting the 743A6 alignment aid..... | 14 |
| 5 Function and description of the components..... | 14 |
| 6 Application examples | 18 |
| 7 Liability..... | 21 |

1 Introduction

INFORMATION

Date of the last update: 2021-02-12

- Please read this document carefully.
- Follow the safety instructions.

Before you install and start operation of the device, please read these operating instructions carefully. Correct handling of the device will avoid possible hazards, lead to reproducible alignment results and make your job easier. If the device is not properly handled, the guarantee


expires and you are responsible for possible damages.

2 Area of application

The orthotic alignment aid is designed for the three-dimensional alignment of lower limb orthoses and the positioning of ankle and knee joint pivot points on a plaster negative or plaster positive. The orthotic alignment aid is secured in a vice for use.

3 Components

- 1 Base part
- 2 Double cone
- 3 Cylinder pin
- 4 Footplate
- 5 Clamping lever
- 6 Heel plate
- 7 Adjustment column
- 8 Guide column
- 9 Holding fork
- 10 Height adjustment unit
- 11 Star grip screw
- 12 Lower slider plate
- 13 Spindle
- 14 Rod attachments, 4
- 15 Square shafts, 2


4 Mounting the 743A6 alignment aid


The orthotic alignment aid can be secured in a vice using the base part (1).

The integrated cylinder pins (3) in the base part make horizontal alignment in the vice easier.


5 Function and description of the components


The footplate attached to the base part (4) can be adjusted by loosening the clamping lever (5).

Inward and outward rotation of the foot and length adjustment can be performed. On the back of the footplate, the rotation adjustment of the footplate can be read off a scale up to **15°** (arrow).

A heel height adjustment of up to 125 mm is made with the heel plate (6) and can be read off the adjustment column (7).

Height-adjustable holding forks (9) for the ankle joint and knee joint are attached to the guide columns (8). To adjust the height of the holding forks, simply loosen the clamping lever on the height adjustment unit (10) and then retighten it.


The minimum distance of the ankle holding fork to the footplate is **30 mm**; the maximum distance for the knee holding fork is **580 mm**.

The star grip screw (11) is loosened so both holding forks can be adjusted in the AP direction. The lower slider plate (12) can be displaced by **50 mm**.


For a different rotational adjustment of the horizontal knee and ankle axis, the corresponding holding fork can be rotated **30°** to either side by loosening the two star grip screws in each case.


By adjusting the spindle (13), a clearance width of up to **210 mm** can be achieved in both holding forks.

The spindle attachment can be removed from the holding fork by unscrewing the corresponding star grip screw.

If the plaster model is to be aligned with round steel, the rod attachments (14) provided for this are used.

If the adjustment set for 743R6 orthotic joints is used during orthosis alignment, the included square shafts (15) can be inserted into the holding forks for work on the plaster negative.

6 Application examples

Transfer of the relevant markings (e.g., knee and ankle pivot points) in the plaster negative.

Place the plaster negative on the footplate and clamp it between the holding forks at the knee and ankle by adjusting the spindles.

When the model is secured, the pivot points (e.g., knee pivot points according to Nietert) can be traced easily. During the AP measurement, the calliper can rest on the holding fork.

After the adjustment of the heel height and the foot position, the frontal and sagittal alignment are performed with the help of the 743L30=* LaserLine device.

If the pivot points are to be marked with round steel, the spindle attachments are replaced with the rod attachments (see sec. 5, note on (14)), and the plaster model is aligned in the holding forks.

Using the square shafts creates a square channel in the plaster positive after the model is poured out, which makes orthosis alignment possible with the proven adjustment adapters of the adjustment set for 743R6 orthotic joints.

Preparation of the plaster negative for pouring out: to preserve the specified pivot points in every position, the square shafts are secured after greasing with 640Z5=1/640Z5=5 plaster isolation cream with plaster languettes.

The forefoot can be built up with plaster while the plaster positive is clamped in. When doing this, the footplate must be isolated with plaster isolation cream at the contact points for protection.

Finished plaster positive with marked construction lines.

7 Liability

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according

to the instructions for use. The manufacturer is not responsible for damages caused by components and spare parts not approved by the manufacturer.

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