1C59 Taleo Adapt

DE	Gebrauchsanweisung	9
EN	Instructions for use	18
FR	Instructions d'utilisation	27
IT	Istruzioni per l'uso	37
ES	Instrucciones de uso	47
PT	Manual de utilização	57
NL	Gebruiksaanwijzing	66
SV	Bruksanvisning	76
DA	Brugsanvisning	85
NO	Bruksanvisning	94
FI	Käyttöohje	103
PL	Instrukcja użytkowania	112
HU	Használati utasítás	122
CS	Návod k použití	131
RO	Instrucțiuni de utilizare	140
HR	Upute za uporabu	150
SL	Navodila za uporabo	159
SK	Návod na používanie	168
BG	Инструкция за употреба	177
TR	Kullanma talimatı	188
EL	Οδηγίες χρήσης	197
RU	Руководство по применению	208
JA	取扱説明書	218
ZH	使用说明书	227
KO	사용 설명서	233

ISO 22675 - P n - m kg*)


*) Body mass limit not to be exceeded!

Stiffness	Weight (kg)	Label text
1	51	ISO 22675-P3*-R3-51-kg
2	58	ISO 22675-P3-R3-58-kg
3	67	ISO 22675-P3-R3-67-kg
4	77	ISO 22675-P4-R3-77-kg
5	88	ISO 22675-P4-R3-88-kg
6	100	ISO 22675-P5-R3-100-kg
7	115	ISO 22675-P5-R3-115-kg
8	130	ISO 22675-P6-R3-130-kg

*Actual test load reflect maximum body mass.

1

Bench alignment preparation

Heel wedge:

Try the firmest (black) heel wedge first.

Hydraulic ankle:


Set the ankle to the neutral position (solid line).

Alignment reference line:

On the alignment marking on the ankle, see zoom in above. Ignore the sagittal plane alignment marks on the footshell.

2

Bench alignment, transtibial

Transfemoral or knee disarticulation prosthesis: Position the alignment reference line on the alignment marking on the ankle. Check the instructions for use of the prosthetic knee joint.

3

Range of motion


Max. dorsiflexion (dashed line)

Max. plantarflexion (dashed line)

Neutral position (solid line)
2° from max. dorsiflexion
10° from max plantarflexion


Plantarflexion resistance (P): Determines how quickly the prosthetic foot plantarflexes during loading response.

Dorsiflexion resistance (D): Influences the movement of the ankle from loading response to the end of the stance phase.

5




Static alignment

1. Move foot in sagittal plane: Posterior, if user is falling backwards. Anterior, if user is falling forwards.

2. Adjust the plantarflexion/dorsiflexion angle in the sagittal plane:

Use the alignment adjustment screws on the adapter at the top of the foot. Adjust the screws until the solid line (neutral position) on the ankle is showing against the edge of the frame.

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INFORMATION

Date of last update: 2024-04-23

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The 1C59 Taleo Adapt prosthetic foot is equipped with a hydraulic ankle joint. The ankle joint has a range of motion of 12°. When the prosthesis is aligned in the neutral position, 10° of plantar flexion will support loading

response when walking on level ground as well as down ramps. The dorsiflexion of 2° facilitates an easy rollover and walking up ramps. In addition to the hydraulics, the carbon and polymer spring elements also enable perceptible plantar flexion. They also support natural rollover and offer a high energy return. An exchangeable heel wedge provides absorption at heel strike.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

The MOBIS classification describes the mobility grade and body weight, and makes it easy to identify compatible components.

The product is recommended for mobility grade 3 (unrestricted outdoor walker) and can be used for mobility grade 4 (unrestricted outdoor walker with particularly high demands).

The following table shows the appropriate spring stiffness for the level of load exerted by the user on the prosthetic foot.

TIP: Also try the stiffness of the higher load level/weight class so the user can experience the difference.

Spring stiffness relative to body weight and load level		
Body weight [kg] (lbs)	Medium load level	High load level
Up to 51 (112)	1	2
52 (115) to 58 (128)	2	3
59 (130) to 67 (148)	3	4
68 (150) to 77 (170)	4	5
78 (172) to 88 (194)	5	6
89 (196) to 100 (220)	6	7
101 (223) to 115 (254)	7	8
116 (256) to 130 (287)	8	

2.3 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2.4 Environmental conditions

Allowable environmental conditions
Temperature during use: -10 °C to +45 °C
Chemicals/liquids: Fresh water, salt water, perspiration, urine, soapsuds, chlorine water
Moisture: Submersion: max. 1 h in 3 m depth, relative humidity: no restrictions
Solids: Dust, occasional contact with sand
Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and tear and damage (see page 26).
Prohibited environmental conditions
Solids: Highly liquid-binding particles (e.g. talcum), continuous contact with sand
Chemicals/liquids: Acids, continuous use in liquid media
Storage and transport
Storage temperature: -20 °C to +60 °C, relative humidity: 20 % to 90 %, no mechanical vibrations or impacts

2.5 Lifetime

Prosthetic foot

Depending on the patient's activity level, the maximum lifetime of the product is 3 years.

Footshell, protective sock

The product is a wear part, which means it is subject to normal wear and tear.

3 General safety instructions

CAUTION! Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ Note the combination possibilities/comboination exclusions in the instructions for use of the products.
- ▶ Observe the maximum lifetime of the product.
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not expose the product to excessive strain (section: "Area of application" see page 19) or to prohibited environmental conditions (see page 20).
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

- ▶ The product may only be used by one person; reuse on another person is prohibited.
- ▶ To avoid the risk of pinching, do not reach into the joint mechanism.

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover) or delamination of the spring are indications of loss of functionality. Play in the hydraulic ankle joint, oil leakage and unusual noises may indicate a loss of functionality.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Prosthetic foot	–
1	Protective sock	SL=Spectra-Sock-7
1	Heel wedge set	2F50=*

Spare parts/accessories	
Designation	Reference number
Footshell	2C8=*
Connection cover	2C21=*

5 Preparing the product for use

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

NOTICE

Altering the prosthetic foot or footshell

Premature wear resulting from damage to the product

- ▶ Do not alter the prosthetic foot or footshell.

INFORMATION

A plastic part protects the connecting section of the product from scratches during the alignment and fitting.

- ▶ Remove the plastic part before the patient leaves the fitting area.

5.1 Applying/removing the footshell

INFORMATION

- ▶ Pull the protective sock over the prosthetic foot to prevent noises in the footshell.

▶ CAUTION! Always use the prosthetic foot with the footshell.

Apply or remove the footshell as described in the footshell instructions for use.

5.2 Bench Alignment

INFORMATION

Use the prosthetic foot only with an installed heel wedge.

Preparation

Three angle indicator lines are on the side of the ankle joint of the prosthetic foot. The lines indicate the neutral position (solid line), maximum dorsiflexion and maximum plantar flexion (both indicated by dashed lines). The neutral position is 2 degrees away from maximum dorsiflexion.

- 1) Insert the firmest (black) heel wedge into the prosthetic foot and start the alignment with it.
- 2) Move the prosthetic foot into the neutral position using the lines on the side of the ankle joint (see fig. 3):
- 3) Set the two valves to maximum resistance (clockwise) so that the neutral position is maintained during the bench alignment process (see fig. 4).

Bench alignment

- > **Required materials:** 662M4 goniometer, 743S12 heel height measuring device, 743A80 50:50 gauge, alignment apparatus (e.g. 743A200 PROS.A. Assembly)
 - > **Illustration to prepare for bench alignment (see fig. 1)**
 - > **Transfemoral or knee disarticulation prosthesis:** Observe the instructions for use for the prosthetic knee joint.
- 1) Position the prosthetic foot in the alignment apparatus.
 - **Heel height: Effective heel height of everyday shoe without additional elevation under the heel.**
 - **Ankle joint:** Neutral position, both valves at maximum resistance.
 - 2) **Sagittal plane:** Align the prosthetic foot with the alignment reference line.

→ Alignment reference line through the alignment marking on the ankle joint of the prosthetic foot.

INFORMATION: Ignore the markings on the footshell and known general alignment guidelines.

- 3) **Frontal plane:** Align the prosthetic foot with the alignment reference line.
 - Alignment reference line through the markings on the footshell (**external rotation approx. 5°**).
- 4) Install the structural components between the prosthetic foot and prosthetic socket.
- 5) **Sagittal plane:** Align the prosthetic socket with the alignment reference line.
 - Determine the centre of the prosthetic socket proximally and distally with the 50:50 gauge and mark the centre line.
 - Mark the socket reference point on the centre line: **For transtibial prostheses, at the height of the centre of the patella.**
 - Place the socket so that the alignment reference line runs vertically through the socket reference point.
 - Adjust the socket flexion by rotating around the socket reference point: **individual residual limb flexion/contracture + 5°**
- 6) **Frontal plane:** Align the prosthetic socket with the alignment reference line.
 - Take the abduction or adduction of the residual limb into account.
- 7) **Before the end of the bench alignment process, check whether the prosthetic foot is still in the neutral position (solid line) (see fig. 3).**
 - If the prosthetic foot is not still in the neutral position, move the ankle back to the neutral position and recheck the alignment.
 - If the prosthetic foot is still in the neutral position, continue on to the static alignment.

5.3 Static Alignment

During the static alignment process, the settings of the prosthesis (length, angle, displacements) are checked and adjusted on the user. The aim is to allow the user to stand in a relaxed manner while the prosthetic foot is in the neutral position.

> **Ottobock recommends checking the alignment of the prosthesis using the 3D L.A.S.A.R. Posture (see fig. 5).**

- 1) Set the two valves of the prosthetic foot to the medium resistance so the prosthetic foot can move.

- 2) Ask the user to put on the prosthesis, stand comfortably and distribute their weight evenly across both legs.
- 3) Check the length of the prosthesis and adjust it as necessary.
- 4) **Shift the prosthetic foot in the sagittal plane:** The user should stand comfortably while using minimal muscle strength to maintain the position.
 - If the user has the feeling that they are falling backwards, move the prosthetic foot in the posterior direction.
 - If the user has the feeling that they are falling forwards, move the prosthetic foot in the anterior direction.
- 5) **Adjusting the angle in the sagittal plane:** The dorsiflexion and plantar flexion ranges can be changed using the adjustment screws of the adapter on top of the prosthetic foot.
 - **If the solid neutral line is covered:** There is too much range for dorsiflexion movement and not enough plantar flexion range of motion. Add more plantarflexion above the ankle until the neutral line is visible.
 - **If space is visible between the dorsiflexion stop edge of the frame and the solid neutral line:** There is too much plantar flexion range of motion and not enough dorsiflexion range. Add more dorsiflexion above the ankle until the neutral line sits right along the frame.
 - It is possible to align the prosthetic foot at maximum dorsiflexion if the user prefers this or has a high need for stability (e.g. bilateral amputee).
- 6) **Checking the load line with the 3D L.A.S.A.R. Posture:** The projected load line allows you to control the prosthetic alignment precisely:
 - **Knee joint:** The load line should be **15 mm** in front of the compromise pivot point according to Nietert. Due to the kinematics of the prosthetic foot, the distance between the load line may be smaller.
 - **Prosthetic knee joint:** Observe the instructions for use for the prosthetic knee joint.
 - **Prosthetic foot:** The load line should run close to the alignment marking on the prosthetic ankle.

5.4 Dynamic Trial Fitting

Once optimal static alignment has been achieved, the alignment of the prosthesis should not be changed in the sagittal plane. The behaviour of the prosthetic foot, from loading response to rollover is adjusted using the res-

instances of the hydraulics. The alignment in the frontal plane is checked and adjusted to achieve an optimal gait pattern.

Adjusting the hydraulics

The resistances of the hydraulics can be set separately for plantar flexion and dorsiflexion. Start with the plantar flexion resistance. After the initial adjustment, further adjustments are made during the fitting process as the user becomes accustomed to walking with the prosthetic foot.

- 1) **Adjust the plantar flexion resistance.** This determines how quickly the prosthetic foot plantarflexes during loading response after heel strike.
 - **Resistance too low:** The forefoot drops down too quickly (audible clapping possible). Or the user complains that it is hard to move over the foot.
 - **Resistance too high:** The heel feels too stiff to the user. Knee flexion is stimulated to a greater extent and the foot has a more propulsive effect.
- 2) **Adjust the dorsiflexion resistance.** This influences the movement of the ankle after loading response up to the end of stance phase.
 - **Resistance too low:** The transition to the dorsiflexion stop feels quick and hard to the user.
 - **Resistance too high:** The rollover process feels tough and strenuous to the user (like walking uphill).

Tips for transtibial prostheses

- 1) **Pay attention to knee flexion during loading response.**
 - If the knee remains fully extended during loading response, increase the plantar flexion resistance.
 - If increasing the resistance is insufficient, try a firmer heel wedge.
- 2) **Avoid medial movement of the knee.**
 - If the knee moves in the medial direction in the first half of stance phase, move the prosthetic foot in the medial direction.
 - If the knee moves in the medial direction in the second half of stance phase, reduce the exterior rotation of the prosthetic foot.

Tips for transfemoral and knee disarticulation prostheses

- ▶ The hydraulic plantar flexion of the Taleo Adapt can influence functions of the prosthetic knee joint (e.g. initiating the yielding movement or triggering stance phase control by bearing weight).
- ▶ Observe the instructions for use for the prosthetic knee joint.

5.4.1 Replacing the heel wedge

In addition to changing the plantar flexion resistance, the behaviour of the prosthetic foot at heel strike and loading response can be adjusted by exchanging the heel wedge. Heel wedges of varying firmness are included in the scope of delivery.

Heel wedge firmness: The colour of the heel wedge indicates the firmness (see fig. 6). Ottobock recommends starting with the firmest heel wedge.

- 1) Bend the heel spring of the prosthetic foot slightly and remove the installed heel wedge.
- 2) Align the other heel wedge so the Ottobock lettering is right side up and the tip points in the anterior direction.
- 3) Insert the heel wedge into the prosthetic foot (see fig. 9). To do so, first insert the tip into the front part of the retainer ①, then push the heel wedge into the rear of the retainer ②.

6 Cleaning

> **Allowable cleaning agent:** pH-neutral soap (e.g. 453H10 Derma Clean)

- 1) **NOTICE! To avoid product damage, only use the allowable cleaning agents.**

Clean the product with clear fresh water and a pH neutral soap.

- 2) Remove dirt from water drainage contours using a toothpick and rinse.
- 3) Rinse the soap away with clear fresh water. In doing so, rinse the foot-shell until all dirt has been removed.
- 4) Dry the product with a soft cloth.
- 5) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ During regular inspection: Inspect the prosthesis for signs of wear and check its function.
- ▶ Conduct annual safety inspections.

8 Disposal

Do not dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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