Kenevo 3C60/3C60=ST 647G1415=all_INT
 C-Leg 3C98-3*/3C88-3* G647G1375=all_INT
 Genium 3B1-2* 647G1148/647G1389=all_INT
 Genium 3B1-3* 647G1380=all_INT
 Genium X3 3B5-2* 647G1149/647G1388=all_INT
 Genium X3 3B5-3* 647G1374=all_INT

Taleo 1C50
 Taleo Vertical Shock 1C51
 Taleo Harmony 1C52
 Taleo Low Profile 1C53
 Empower 1A1-2
 Maverick Comfort AT
 Maverick Xtreme
 Maverick Xtreme AT
 Promenade
 Kintrol
 Restore
 Thrive
 Freestyle Swim

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e.g. 31 / 7: first number = foot size [cm], second number = Stiffness category [1 to 9]

Kenevo 3C60 / 3C60=ST

- **Empfohlene Kombinationen / Recommended combinations:** Taleo 1C50, Taleo Vertical Shock 1C51, Taleo Harmony 1C52, Taleo Low Profile 1C53, Kintrol¹, Restore¹, Promenade¹

Körpergewicht/ Body weight	1C50	1C53
≤ 115 kg (253 lbs)	22-30 / 1-7	22-30 / 1-7
≤ 125 kg (275 lbs)	–	22-29 / 1-8

¹ Ottobock Systemhöhe beachten / observe Ottobock systemheight

C-Leg 3C98-3*/3C88-3*

- **Empfohlene Kombinationen / Recommended combinations:** Taleo Vertical Shock 1C51, Taleo Harmony 1C52, Taleo Low Profile 1C53, Empower 1A1-2, Maverick Comfort AT¹, Promenade¹, Thrive¹

Körpergewicht/ Body weight	1C51 / 1C52	1C53	1A1-2	Promenade ¹	Maverick Comfort AT ¹	Thrive ¹
≤ 100 kg (220 lbs)	22-30 / 1-7	29, 30 / 1-6	22-30 / 1-7	29-31 / 1-3	23-27 / 1-9 28 / 1-7	22-26 / 1-9
≤ 115 kg (253 lbs)		22-28 / 1-7		22-28 / 1-4	–	–
≤ 125 kg (275 lbs)	–	–	26-27 / 1-8	–	–	–
≤ 136 kg (299 lbs)	–	–	–	–	–	–

¹ Ottobock Systemhöhe beachten / observe Ottobock systemheight

Genium 3B1-2*, Genium 3B1-3*

- **Empfohlene Kombinationen / Recommended combinations:** Taleo Vertical Shock 1C51, Taleo Harmony 1C52, Taleo Low Profile 1C53, Empower 1A1-2, Maverick Comfort AT¹, Maverick Xtreme¹, Maverick Xtreme AT¹, Thrive¹

Körpergewicht/ Body weight	Maverick Xtreme AT ¹	Thrive ¹
≤ 125 kg (275 lbs)	23-30 / 1-9	22-31 / 1-9
≤ 150 kg (330 lbs)	23-27 / 1-9 28 / 1-7	22-26 / 1-9

¹ Ottobock Systemhöhe beachten / observe Ottobock systemheight

Genium X3 3B5-2*, Genium X3 3B5-3*

- **Empfohlene Kombinationen / Recommended combinations:** Taleo Vertical Shock 1C51, Taleo Harmony 1C52, Taleo Low Profile 1C53, Empower 1A1-2, Maverick Comfort AT¹, Maverick Xtreme¹, Maverick Xtreme AT¹, Thrive¹, Freestyle Swim¹

Körpergewicht / Body weight	Freestyle Swim ¹
≤ 100 kg (220 lbs)	22-31 / 1-6
≤ 125 kg (275 lbs)	–

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INFORMATION

Date of last update: 2021-04-09

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

The content of this additional documentation extends the permissible combination possibilities according to the instructions for use listed on the cover sheet.

This additional documentation is intended for qualified personnel that have been authorised by Ottobock through training.

The content of this additional documentation does not replace the instructions for use included with the respective products.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

INFORMATION

Calculating the Ottobock system height for the Maverick Comfort AT, Maverick Xtreme, Maverick Xtreme AT, Kintrol, Restore, Promenade, Thrive, Freestyle Swim prosthetic feet

To calculate the Ottobock system height (e.g. for input in the adjustment software), the build height according to the technical data for the listed prosthetic feet has to be reduced by about 18 mm.

Example: The build height of the "Maverick Comfort AT" prosthetic foot in size 26 is 117 mm.

Therefore, the system height is: $117 \text{ mm} - 18 \text{ mm} = 99 \text{ mm}$. This is only a reference value. Therefore, measure and verify the distances on the patient before shortening the tube adapter.

2 Safety instructions

⚠ CAUTION

Failure to observe the table on the inside cover

Falling due to breakage of load-bearing components of the knee joint.

- ▶ Depending on the knee joint used and the patient's body weight, the listed prosthetic feet may be combined only with the respective described foot sizes [cm] and carbon spring categories (stiffnesses) [1 through 9].
- ▶ If the carbon spring category (stiffness) is subsequently changed, ensure that the specified limit values are not exceeded.

⚠ CAUTION

Errors during prosthesis alignment

Falling due to unexpected behaviour of the product.

- ▶ Regardless of the alignment information in the instructions for use for the prosthetic feet, the alignment must be carried out according to the information in the instructions for use of the knee joints or the adjustment software.
- ▶ Note that the selection of the carbon spring category (stiffness) for the prosthetic feet influences the use of the knee joint stance function.
- ▶ Observe the instructions for use for all prosthesis components.

3 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

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Otto Bock Healthcare Products GmbH