



Vacuum Assisted Socket Systems (VASS) Evidence Essentials

	Mobility need or deficit of the patient	Evidence for benefits of VASS vs. standard sockets in transtibial amputees
Safety	<p>Patient stumbles and/or falls repeatedly</p> <p>Patient avoids activities due to fear of falling</p> <p>Patient sustained fall-related injuries</p>	<ul style="list-style-type: none"> - Significant reduction in risk of multiple falls of up to 75% (Rosenblatt et al., 2017a) - Significant improvements in balance and indicators for the risk of falling, such as Timed-up-and-go-test, Berg Balance Scale, Four Square Step Test (FSST), ABC scale, etc. (Ferraro et al., 2011; Rosenblatt et al., 2017b; Samitier et al., 2014)
Mobility	<p>Patient is limited in his/her mobility</p>	<ul style="list-style-type: none"> - Significant increase in walking speed, also in individuals with dysvascular amputation (Kuntze Ferreira et al., 2015; Rosenblatt et al., 2017b; Samitier et al., 2014) - Significant increase in walking capacity and distance in timed walk test in individuals with dysvascular amputation (Samitier et al., 2014)
Socket fit	<p>Patient has problems with socket fit, residual limb volume fluctuations, pistoning</p>	<ul style="list-style-type: none"> - Significant increase in Socket Comfort Score of up to 3 points (Rosenblatt et al., 2017b) - Significant prevention of residual limb volume loss during the day (Board et al., 2001; Gholizadeh et al., 2016; Goswami et al., 2003; Kahle et al., 2014; Sanders et al., 2011; Stevens et al., 2019; Youngblood et al., 2020) - Significantly reduced pistoning / relative movement between residual limb and socket with reduction of pressures and shear forces (Beil et al., 2022; Board et al., 2001; Darter et al., 2016; Gholizadeh et al., 2016; Kahle et al., 2014; Klute et al., 2011; Stevens et al., 2019) - Significantly reduced perspiration in the socket (Klute et al., 2016)

Residual limb health	Patient suffers from recurring residual limb pressure sores, ulcers, or wound healing delay after amputation	<ul style="list-style-type: none"> - Significantly earlier prosthesis fitting and mobilization after amputation without interference with wound healing (Traballesi et al., 2014) - Significantly improved mobility in the first 5 months after amputation despite wound healing delay (Traballesi et al., 2012) - Patients can maintain prosthetic mobility despite residual limb ulcers without interfering with wound healing (Hoskins et al., 2013) - Significantly improved perfusion of the residual limb (Rink et al., 2016)
Musculo-skeletal pain	Patient suffers from joint and back pain due to gait asymmetry and excessive loading	<ul style="list-style-type: none"> - Significant improvement in gait symmetry and, thus, loading of the locomotor system (Board et al., 2001; Kuntze Ferreira et al., 2015)

References

- Beil TL, Street GM, Covey SJ. Interface pressures during ambulation using suction and vacuum-assisted prosthetic sockets. *J Rehabil Res Dev* 2002; 39(6): 693-700.
- Board WJ, Street GM, Caspers C. A comparison of trans-tibial amputee suction and vacuum socket conditions. *Prosthet Orthot Int* 2001;25(3):202-209.
- Darter BJ, Sinitski K, Wilken JM. Axial bone–socket displacement for persons with a traumatic transtibial amputation: The effect of elevated vacuum suspension at progressive body-weight loads. *Prosthet Orthot Int* 2016;40(5):552-557. [Download](#)
- Ferraro C. Outcomes study of transtibial amputees using elevated vacuum suspension in comparison with pin suspension. *J Prosthet Orthot* 2011;23(2):78-81.
- Gholizadeh H, Lemaire ED, Eshragi A. The evidence base for elevated vacuum in lower limb prosthetics: Literature review and professional feedback. *Clin Biomech* 2016; 37: 108-116.
- Goswami J, Lynn R, Street GM, Harlander M. Walking in a vacuum-assisted socket shifts stump fluid balance. *Prosthet Orthot Int* 2003; 27(2): 107-113. [Download](#)
- Hoskins RD, Sutton EE, Kinor D, Schaeffer JM, Fatone S. Using vacuum-assisted suspension to manage residual limb wounds in persons with transtibial amputation: A case series. *Prosthet Orthot Int* 2013;38(1):68-74.
- Kahle JT, Orriola JJ, Johnston W, Highsmith MJ. The effects of vacuum-assisted suspension on residual limb physiology, wound healing, and function: a systematic review. *Technol Innov* 2014; 15: 333-341.
- Klute GK, Berge JS, Biggs W, Pongnumkul S, Popovic Z, Curless B. Vacuum-assisted socket suspension compared with pin suspension for lower extremity amputees: Effect on fit, activity, and limb volume. *Arch Phys Med Rehabil* 2011; 92(10):1570-1575.
- Klute GK, Bates KJ, Berge JS, Biggs W, King C. Prosthesis management of residual-limb perspiration with subatmospheric vacuum pressure. *J Rehabil Res Dev* 2016; 53(6): 721-728.
- Kuntze Ferreira AE, Neves EB. A comparison of vacuum and KBM prosthetic fitting for unilateral transtibial amputees using the Gait Profile Score. *Gait Posture* 2015;41(2):683-687.

Rink C, Wernke MM, Powell HM, et al. Elevated vacuum suspension preserves residual-limb skin health in people with lower-limb amputation: Randomized clinical trial. *J Rehabil Res Dev*. 2016;53(6):1121-1132.

Rosenblatt NJ, Ehrhard T. The effect of vacuum assisted socket suspension on prospective, community based falls by users of lower limb prostheses. *Gait Posture* 2017a;55:100-104.

Rosenblatt NJ, Ehrhardt T, Fergus R, Bauer A, Caldwell R. Effects of vacuum-assisted socket suspension on energetic costs of walking, functional mobility, and prosthesis-related quality of life. *J Prosthet Orthot* 2017b; 29(2): 65-72.

Sanders JE, Harrison DS, Myers TR, Allyn KJ. Effects of Elevated Vacuum on In-Socket Residual Limb Fluid Volume: Case Study Results Using Bioimpedance Analysis. *J Rehabil Res Dev* 2011;48(10):1231-48.

Samitier CB, Guirao L, Costea M, Camos JM, Pleguezuelos E. The benefits of using a vacuum-assisted socket system to improve balance and gait in elderly transtibial amputees. *Prosthet Orthot Int* 2014;40(1):83-88.

Stevens PM, DePalma RR, Wurdeman SR. Transtibial socket design, interface, and suspension: a clinical practice guideline. *J Prosthet Orthot* 2019; 31(3): 172-178.

Traballesi M, Delussu AS, Fusco A, Iosa M, Averna T, Pellegrini T, Brunelli S. Residual limb wounds or ulcers heal in transtibial amputees using an active suction socket system. A randomized controlled study. *Eur J Phys Rehabil* 2012; 48(4):618-623.

Youngblood RT, Brzostowski JT, Hafner BJ, Czerniecki JM, Allyn KJ, Foster RL, Sanders JE. Effectiveness of elevated vacuum and suction prosthetic suspension systems in managing daily residual limb fluid volume change in people with transtibial amputation. *Prosthet Orthot Int* 2020 Mar 18:309364620909044. doi: 10.1177/0309364620909044. Online ahead of print.

Ottobock North America, Reimbursement
P 800 328 4058 F 800 230 3962
US <https://shop.ottobock.us>
CA <https://shop.ottobock.ca>
reimbursement911@ottobock.com