
Lamination Anchors

DE Gebrauchsanweisung	5	CS Návod k použití.....	107
EN Instructions for use	12	RO Instrucțiuni de utilizare.....	115
FR Instructions d'utilisation.....	20	HR Upute za uporabu	123
IT Istruzioni per l'uso	29	SL Navodila za uporabo	130
ES Instrucciones de uso	37	SK Návod na používanie	138
PT Manual de utilização.....	45	BG Инструкция за употреба	146
NL Gebruiksaanwijzing.....	53	TR Kullanma talimatı	154
SV Bruksanvisning	61	EL Οδηγίες χρήσης	162
DA Brugsanvisning	69	RU Руководство по применению	170
NO Bruksanvisning	76	JA 取扱説明書.....	179
FI Käyttöohje.....	84	ZH 使用说明书.....	187
PL Instrukcja użytkowania.....	91	KO 사용 설명서.....	194
HU Használati utasítás	99		

5

6

7

8

9

10

11

12

13

1 Produktbeschreibung

Deutsch

INFORMATION

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INFORMATION

Date of last update: 2020-05-04

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

1.1 Construction and Function

Lamination anchors are laminated into a prosthetic socket. They serve to connect the prosthetic socket to the distal prosthetic components. This document applies to the following lamination anchors:

Lamination anchor, four-arm

Reference number	Distal connection	Figure
4R111	Pyramid receiver, rotating	see fig. 2
4R111=N, 4R111=T	Threaded connector	see fig. 3
4R116, 4R116=T, 4WR95=2	Pyramid adapter, rotating	see fig. 4
4R117, 4R117=T	Pyramid adapter, rotating	see fig. 5
4R119=N, 4R119=NT	Threaded connector	see fig. 6
4R119, 4R119=T, 4WR95=1	Pyramid receiver, rotating	see fig. 7

Lamination anchor, three-arm

Reference number	Distal connection	Figure
4R41	Pyramid receiver, rotating	see fig. 8
4R42	Pyramid adapter	see fig. 9
4R42=1, 4R42=5	Pyramid adapter with 9.5 mm bore for shuttle lock	
4R43	Threaded connector	see fig. 10
4R89	Pyramid adapter, rotating	see fig. 11

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

Product approved for max. body weight

- The maximum approved body weight is specified in the technical data (see Page 19).

Additional constraints

- 4R41, 4R42, 4R42=1, 4R42=5, 4R43, 4R89, 4R111, 4R111=N, 4R116, 4R117, 4R117=T, 4R119=NT, 4R119=T, 4WR95=1:** Use of products permissible in TT prostheses and TF prostheses.
- 4R111=T, 4R116=T, 4R119, 4R119=N, 4WR95=2:** Use of product only permitted in TF prostheses.

2.3 Environmental conditions

Storage and transport

Temperature range -20 °C to +60 °C (-4 °F to +140 °F), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

4R41, 4R42=1, 4R42=5, 4R43, 4R89, 4R111, 4R111=N, 4R116, 4R117

Allowable environmental conditions

Temperature range: -10 °C to +45 °C

Moisture: relative humidity: 20% to 90%, non-condensing

Unacceptable environmental conditions

Chemicals/liquids: fresh water, salt water, perspiration, urine, acids, soapsuds, chlorine water

Solids: dust, sand, highly hygroscopic particles (e. g. talcum)

4R42, 4R119, 4R119=N

Allowable environmental conditions

Temperature range: -10 °C to +45 °C

Moisture: relative humidity: 20% to 90%, non-condensing

Chemicals/liquids: fresh water as dripping water, occasional contact with salty air (e.g. near the ocean)

Solids: dust

Prohibited environmental conditions

Chemicals/moisture: salt water, perspiration, urine, acids, soapsuds, chlorine water

Solids: dust in high concentrations (e.g. construction site), sand, highly hygroscopic particles (e.g. talcum)

4R111=T, 4R116=T, 4R117=T, 4R119=NT, 4R119=T, 4WR95=1, 4WR95=2

Allowable environmental conditions

Temperature range: -10 °C to +45 °C

Chemicals/liquids: fresh water, soapsuds, chlorine water

Moisture: Submersion: max. 1 h in 2 m depth, relative humidity: no restrictions

Solids: Dust, occasional contact with sand

Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage (see Page 18).

Prohibited environmental conditions

Solids: highly hygroscopic particles (e.g. talcum), dust in high concentrations (e.g. construction site), intensive contact with sand

Prohibited environmental conditions

Chemicals/liquids: salt water, perspiration, urine, acids, continuous use in liquid media

2.4 Service life

In general, all modular adapters are tested by the manufacturer for 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum service life of 5 years.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

3.2 General safety instructions

 **CAUTION**

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the product according to the specified area of application (see Page 13).

 **CAUTION**

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

 **CAUTION**

Use under unallowable environmental conditions

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

⚠ CAUTION

Exceeding the service life

Risk of injury due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded.

⚠ CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

Lamination anchor, 3-arm		
Quantity	Designation	Reference number
1	Instructions for use	–
1	Lamination anchor	–
1	For 4R41, 4R89: Cap screw (steel)	501Z2=M5x30
1	For 4R42*: Lamination dummy	4X3
1	For 4R41: Screwed insert, pyramid receiver	4R44=N
4	For 4R41: Set screw	506G3=M8x12-V
1	For 4R89: Screwed insert, pyramid adapter	4R87

Lamination anchor, 4-arm		
Quantity	Designation	Reference number
1	Instructions for use	–
1	Lamination anchor	–
1	Cap screw (steel)	501Z2=M5x30
1	Rounded washer	507U16=5.2-NIR-O
1	Screwed insert, pyramid receiver	
	For 4R111, 4R119:	4R44=N
	For 4R119=T:	4R44=T
4	Set screw	
	For 4R111, 4R119, 4R119=T, 4WR95=1:	506G3=M8x12-V
1	Screwed insert, pyramid adapter	
	For 4R116, 4R117:	4R87
	For 4R116=T, 4R117=T:	4R87=T

Lamination anchor, 4-arm		
Quantity	Designation	Reference number
1	Cap screw, titanium	
	For 4R111=T, 4R116=T, 4R117=T, 4R119=NT, 4R119=T, 4WR95*:	501T24=M5x25

Spare parts/accessories (not included in the scope of delivery)	
Designation	Reference number
Lamination dummy for lamination anchor with threaded connector (flat)	4X46
Lamination dummy for lamination anchor with threaded connector (high)	4X46=ST

5 Preparing the product for use

⚠ CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

⚠ CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

⚠ CAUTION

Improper processing of the anchor arms of the lamination anchor

Risk of injury due to breakage of the anchor arms

- ▶ Only bend anchor arms of lamination anchors that are made of steel.
- ▶ Avoid bending the anchor arms too strongly or too frequently.
- ▶ Use the 711S4* bending irons for bending.

INFORMATION

With waterproof lamination anchors, the scope of delivery includes a titanium screw for improved corrosion protection.

- ▶ After trial fitting the prosthesis, replace the pre-installed steel cap screw with the titanium cap screw. Install the cap screw with Loctite®.

5.1 Laminating a transfemoral socket

> **Required materials:** torque wrench (e.g. 710D20), 711S4* bending iron, 99B81 PVA bag, 623T3* perlon stockinette, 616B1 carbon fibre webbing, 616G6 Dacron felt, 616G12 woven carbon fibre, 616G15 woven carbon fibre stockinette, 617H119 Orthocryl lamination resin 80:20 PRO

- 1) **If the lamination anchor has a threaded connector:** Fully screw in the screwed insert.
- 2) Tighten the cap screw of the threaded clamping using the torque wrench (**10 Nm**).
- 3) Pull a soaked PVA bag over the model.
- 4) Apply one layer of Dacron felt and thin out the fibres along the edge.
- 5) Pull two layers of Perlon stockinette over the model.
- 6) Place three layers of woven carbon fibre (e.g. **15 cm x 15 cm**) with offset fibre alignment on the distal end of the model.
- 7) Pull two layers of Perlon stockinette over the model.

- 8) Pull a soaked PVA bag over the model.
- 9) Complete the lamination process with Orthocryl lamination resin.
- 10) Once the lamination resin has cured, remove the PVA bag.
- 11) **INFORMATION: Align an anchor arm that has been angled at the factory in the posterior direction.**
Move the lamination anchor arms to the correct **a-p** and **m-l** position. Never align the threaded clamping in the anterior or posterior direction.
- 12) Place the lamination anchor at the socket end as required for proper alignment and bring it into position.
- 13) **Optional:** Adapt the anchor arms of the lamination anchor to the model with a bending iron.
- 14) Install the lamination dummy (see Page 18).
- 15) Stabilise and fix the lamination anchor with a sealing resin/talcum powder mixture.
- 16) Pull one layer of perlon stockinette over the model.
- 17) Place two layers of woven carbon fibre (e.g. **15 cm x 15 cm**) offset over the lamination anchor arms.
- 18) Pull one layer of perlon stockinette over the model.
- 19) Pull two layers of carbon woven stockinette over the model.
- 20) Pull two layers of perlon stockinette over the model and tie off below the cap screw.
- 21) Pull a soaked PVA bag over the model.
- 22) Complete the lamination process with Orthocryl lamination resin.
- 23) After the lamination resin has cured, remove the lamination dummy.

5.2 Laminating the transtibial socket

> **Required materials:** torque wrench (e.g. 710D20), 711S4* bending iron, 99B81 PVA bag, 623T3 perlon stockinette, 616B1 carbon fibre webbing, 699B1 fibreglass roving, 616G12 woven carbon fibre, 616G15 woven carbon fibre stockinette, 617H119 Orthocryl lamination resin 80:20 PRO

- 1) **If the lamination anchor has a threaded connector:** Fully screw in the screwed insert.

- 2) Tighten the cap screw of the threaded clamping using the torque wrench (**10 Nm**).
- 3) Pull a soaked PVA bag over the model. The PVA bag should be sealed by heat welding at the distal end.
- 4) Pull two layers of Perlon stockinette over the model.
- 5) Place three layers of woven carbon fibre (e.g. **15 cm x 15 cm**) with offset fibre alignment on the distal end of the model.
- 6) Pull two layers of Perlon stockinette over the model.
- 7) **INFORMATION: Align an anchor arm that has been angled at the factory in the posterior direction.**
Move the lamination anchor arms to the correct **a-p** and **m-l** position. Never align the threaded clamping in the anterior or posterior direction.
- 8) Place the lamination anchor at the socket end as required for proper alignment and bring it into position.
- 9) **Optional:** Adapt the anchor arms of the lamination anchor to the model with a bending iron.
- 10) Install the lamination dummy (see Page 18).
- 11) Stabilise and fix the lamination anchor with a sealing resin/talcum powder mixture.
- 12) Pull the fibreglass roving through the holes in the anchor arms and allow it to hang in a loop. If there are no holes, allow the fibreglass roving to hang around the anchor arms in a loop.
- 13) Place three layers of carbon fibre cloth (e.g. **15 cm x 15 cm**) in an offset arrangement over the lamination anchor arms.
- 14) Pull one layer of Perlon stockinette over the model.
- 15) Pull two layers of carbon woven stockinette over the model.
- 16) Pull two layers of perlon stockinette over the model and tie off below the cap screw.
- 17) Pull a soaked PVA bag over the model.
- 18) Complete the lamination process with Orthocryl lamination resin.
- 19) After the lamination resin has cured, remove the lamination dummy.

5.3 Install the lamination dummy

INFORMATION

- ▶ Use the **4X46** lamination dummy when the lamination anchor is to be used with a screwed insert.
- ▶ Use the **4X46=ST** lamination dummy when the next prosthetic component is to be screwed directly into the lamination anchor (e.g. 4R57=ST). The lamination dummy is somewhat higher proximally so that the prosthetic component does not contact the laminate when fully screwed in.

Lamination anchor without threaded connector

- ▶ Set the **4X3** lamination dummy on the pyramid adapter of the lamination anchor (see fig. 12).

Lamination anchor with threaded connector

- 1) Tighten the cap screw of the threaded clamping to a point where the screwed insert can still be screwed in and out.
- 2) Unscrew the screwed insert from the lamination anchor.
- 3) Screw in the **4X46*** lamination dummy (see fig. 13).

5.4 Assembling the distal prosthetic components

Installing the screwed insert

The threaded connector is used for installing a connecting component with an external thread (e.g. pyramid adapter or prosthetic knee joint).

- > **Required materials:** 710D4 torque wrench, 636K13 Loctite®
- 1) Screw in the screwed insert (e.g. pyramid adapter) completely into the thread of the lamination anchor.
- 2) The screwed insert may be unscrewed by a **maximum of 1/4 turn** out of the thread in order to align the insert.
- 3) **If a titanium screw is on hand:** After trial fitting, replace the steel cap screw with the titanium cap screw.
- 4) Secure the cap screw of the threaded clamping using Loctite® and tighten using the torque wrench (**10 Nm**).

Connecting the pyramid adapter and pyramid receiver

The pyramid adapter is fixed with the set screws of the pyramid receiver.

- > **Required materials:** torque wrench (e.g. 710D20), 636K13 Loctite 241
- 1) **Fitting:**
Screw in the set screws.
Use the torque wrench to tighten the set screws (**10 Nm**).
- 2) **Definitive assembly:**
Secure the set screws with Loctite.
Screw in the set screws.
Pre-tighten the set screws with the torque wrench (**10 Nm**) and then tighten them (**15 Nm**).
- 3) Replace any set screws that are protruding or are recessed too much with appropriate set screws (see selection chart).

Selection table for set screws

Reference number	Length (mm)
506G3=M8X12-V	12
506G3=M8X14	14
506G3=M8X16	16

Alignment

The set screws in the pyramid receiver can be used to make static adjustments during alignment, trial fittings and after the prosthesis is finished.

Replacement and disassembly

The set position of the prosthetic component can be maintained during replacement or disassembly. In order to do this, unscrew the two set screws that are screwed in the furthest and located next to each other.

6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.

- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

⚠ CAUTION

Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

▶ Observe the following maintenance instructions.

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

9.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

10 Technical data

Lamination anchor, three-arm

Reference number	4R41	4R42	4R42=1	4R42=5
System height [mm]	39	-5	-5	-5
Weight [g]	170	130	125	125
Material	Stainless steel			
Max. body weight [kg]	125	150	136	136

Reference number	4R43	4R89
System height [mm]	8	-3
Weight [g]	95	180
Material	Stainless steel	
Max. body weight [kg]	125	125

Lamination anchor, four-arm

Reference number	4R111	4R111=-N	4R111=T	4R116
Weight [g]	155	80	165	165
System height [mm]	44	13	2	13

Reference number	4R111	4R111=- N	4R111=T	4R116
Material	Stainless steel			
Max. body weight [kg]	150	150	125	150

Reference number	4R117	4R119	4R119=N	4R119=- NT
Weight [g]	145	165	165	95
System height [mm]	2	44	44	13
Material	Stainless steel			
Max. body weight [kg]	150	150		

Reference number	4R116=T	4R117=T	4R119=T
Weight [g]	155	145	135
System height [mm]	44	2	44
Material (lamination anchor)	Stainless steel		
Material (screwed insert)	Titanium		
Max. body weight [kg]	125	150	150

Reference number	4WR95=1	4WR95=2
Weight [g]	165	165
System height [mm]	44	2
Material	Stainless steel	
Max. body weight [kg]	150	125

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