

Prosthetic gloves for children

Gloves for the Electric Hand 2000 come in 9 new colours

Our youngest users and their parents set high expectations when it comes to the quality of prosthetic gloves. That is why we focus on robustness and flexibility as two main characteristics in prosthetic gloves for the Electric Hand 2000. Additionally we added nine eye-catching colours to our glove portfolio for your customers to choose from.

The following new colour versions are available

Colour	RAL	Colour code
Black	9005	20
Red	3028	25
Orange	2004	30
Yellow	1023	35
Green	6018	40
Blue	5015	45
Dark blue	5002	48
Purple	4001	50
Pink	4003	55



Benefits

- Suitable for boys and girls
- Good mechanical stability stands for durability, robustness, tear-resistant, minimum abrasion
- To be used in combination with Electric Hand 2000
- Cleaning with Special Cleaner (640F12) and spray Pump Bottle (640F13)

Technical data

Article number	Side	Size	Sleeve length	Sleeve end circumference
8S20=136X41L*	left (L)	5	210 mm	200 mm
8S20=136X41R*	right (R)	5	210 mm	200 mm
8S20=147X45L*	left (L)	5 ½	215 mm	200 mm
8S20=147X45R*	right (R)	5 ½	215 mm	200 mm
8S20=162X56L*	left (L)	6	220 mm	210 mm
8S20=162X56R*	right (R)	6	220 mm	210 mm
8S20=177X64L*	left (L)	6 ½	240 mm	220 mm
8S20=177X64R*	right (R)	6 ½	240 mm	220 mm

* When ordering, please add the colour code to the end of the order number, e.g., 8S20=136x41L25