

**2/4Rx - Titan Articles**

Material number 002/4Rx

Page: 1 of 7

**1. Product and company identification****Product identifier**

Trade name: 2/4Rx - Titan Articles

**Relevant identified uses of the substance or mixture and uses advised against**General use: Titanium-Article for orthopedic procedures  
Reserved for industrial and professional use.**Details of the supplier of the safety data sheet**Company name: Otto Bock Health Care  
Street/POB-No.: 3820 W. Great Lakes Drive  
Postal Code, city: Salt Lake City, UT 84120  
USAWWW: [www.ottobockus.com](http://www.ottobockus.com)

Telephone: +1 (801) 956-2400

Telefax: +1 (801) 956-2401

Dept. responsible for information:

Quality Department,  
Telephone: +1 (801) 954-2304 (7 AM – 3 PM, Mountain Time), Email:  
[USRegulatory@ottobock.com](mailto:USRegulatory@ottobock.com)

Additional information: This safety data sheet pertains to the following products:

2R38 - Tube Adaptor  
2R38=10 - Tube Adaptor  
2R57 - Internal Tube Adaptor  
2R58 - Internal Tube Adaptor  
2R216 - SACH Shin Kit  
2R217 - Single Axis Shin Kit  
2R219 - Universal Tube Kit  
4R52 - Tube Clamp Adaptor  
7E7-T - Modular Single Axis Hip Joint (Titan)**Emergency phone number****CHEMTREC, Telephone: +1 (800) 424-9300****2. Hazards identification****Emergency overview**Appearance: Form: solid  
Color: gray

Odor: odorless

Classification: This substance is classified as not hazardous.

**Regulatory status**

This material is not considered hazardous by the U.S. OSHA Hazard Communication Standard (29 CFR 1910.1200) and SIMDUT in Canada.

**Hazards not otherwise classified**

Processing, e.g. by cutting, sawing or grinding, can produce particles and dust. For risks which have to be observed thereby, see chapter 7: Handling, chapter 8: Exposure controls / personal protection and chapter 11: Toxicology.

Titanium-dust: Danger of dust explosion.

see section 11: Toxicological information

**3. Composition / Information on ingredients**

Chemical characterization: Article of Titanium

**4. First aid measures**

In case of inhalation: In case of troubles after inhalation of dust:  
Move victim to fresh air. Seek medical attention.

Following skin contact: Titanium-dust:  
Wash with plenty of water. In case of troubles:  
Take off immediately all contaminated clothing. Seek medical attention.

After eye contact: Titanium-dust:  
With eyelids open, wash out eyes for several minutes under flowing water. In case of troubles or persistent symptoms, consult an ophthalmologist.

After swallowing: Swallowing is not regarded as a possible way of exposition.  
Titanium-dust:  
If person is clearly conscious, have them drink two glasses of water to dilute ingested material.  
Do not induce vomiting. Seek medical attention.

**Most important symptoms/effects, acute and delayed**

Dust:  
In case of eye contact / in case of inhalation: May cause irritations.  
Lung damage is possible in a chronic situation.

**Information to physician**

Treat symptomatically.

**5. Fire fighting measures**

Flash point/flash point range:  
no data available

Auto-ignition temperature: Titanium-dust: approx. 480.2 °F  
Titanium-in pieces: approx. 2199.2 °F

Suitable extinguishing media:  
Product is non-combustible. Extinguishing materials should therefore be selected according to surroundings.

Titanium-dust:  
Extinguishing powder on the basis of NaCl or pulverized limestone.  
In case of fire, use dry sand or fire extinguisher of fire class D. Never use water.

Extinguishing media which must not be used for safety reasons:

Titanium-dust:  
 Never extinguish with a halon or carbon dioxide extinguisher or water.  
 Not a foam extinguisher.

**Specific hazards arising from the chemical**

Titanium-dust: Danger of dust explosion.  
 Titanium-dust, burning:  
 After contact with water: Danger of explosion!

**6. Accidental release measures**

- Personal precautions: Avoid generation of dust. Do not breathe dust.  
 Provide adequate ventilation.  
 In the case of the formation of dust: Wear protective equipment. Avoid contact with skin and eyes.
- Environmental precautions: Do not allow to penetrate into soil, waterbodies or drains.
- Methods for clean-up: Metal parts/dust:  
 Take up mechanically, placing in appropriate containers for disposal.  
 Avoid generation of dust.
- Additional information: Titanium-dust: Eliminate all ignition sources if safe to do so.

**7. Handling and storage**

**Handling**

- Advices on safe handling: For mechanical processing:  
 Avoid respiration of swarf. Wear protective equipment.  
 Provide adequate ventilation. Keep workplace dry.
- Precautions against fire and explosion:  
 For mechanical processing:  
 Avoid generation of dust. Danger of dust explosion.  
 Keep away from combustible material. Keep away from sources of ignition.

**Storage**

- Requirements for storerooms and containers:  
 Store in a dry place.

**8. Exposure controls / personal protection**

**Exposure guidelines**

Occupational exposure limit values:

Type	Limit value
USA: ACGIH: TWA	10 mg/m <sup>3</sup>
USA: ACGIH: TWA	3 mg/m <sup>3</sup>
USA: OSHA: TWA	15 mg/m <sup>3</sup>
USA: OSHA: TWA	5 mg/m <sup>3</sup>

**Engineering controls**

Provide good ventilation and/or an exhaust system in the work area.  
See also information in chapter 7, section storage.

**Personal protection equipment (PPE)**

- Eye/face protection For mechanical processing:  
Tightly sealed safety glasses according to OSHA Standard - 29 CFR: 1910.133 or ANSI Z87.1-2003. or face protection shield.
- Skin protection For mechanical processing: Wear suitable protective clothing.  
For mechanical processing:  
Protective gloves according to OSHA Standard - 29 CFR: 1910.138.  
Observe glove manufacturer's instructions concerning penetrability and breakthrough time.
- Respiratory protection: For mechanical processing:  
Respiratory protection must be worn whenever the TLV (WEL) levels have been exceeded. Particulates filter P2 according to OSHA Standard - 29 CFR: 1910.134 or ANSI Z88.2.
- General hygiene considerations:  
For mechanical processing:  
Do not breathe dust.  
Keep workplace dry.  
Wash hands before breaks and after work.  
Provide a conveniently located eye rinse station.

**9. Physical and chemical properties**

**Information on basic physical and chemical properties**

- Appearance: Form: solid  
Color: gray
- Odor: odorless
- Odor threshold: no data available
- pH value: no data available
- Melting point/freezing point: 3034.4 - 3050.6 °F
- Initial boiling point and boiling range: 5900 - 6332 °F
- Flash point/flash point range: no data available
- Evaporation rate: no data available
- Flammability: no data available
- Explosion limits: LEL (Lower Explosion Limit): Titanium-dust, dry: approx. 50g/m<sup>3</sup>
- Vapor pressure: no data available
- Vapor density: no data available
- Density: 4.5 g/cm<sup>3</sup>
- Water solubility: insoluble
- Partition coefficient: n-octanol/water: no data available
- Auto-ignition temperature: Titanium-dust: approx. 480.2 °F  
Titanium-in pieces: approx. 2199.2 °F
- Thermal decomposition: no data available
- Additional information: no data available

**2/4Rx - Titan Articles**

Material number 002/4Rx

Page: 5 of 7

**10. Stability and reactivity**

- Reactivity: no data available
- Chemical stability: Product is stable under normal storage conditions.
- Possibility of hazardous reactions  
Titanium-dust: Danger of dust explosion.
- Conditions to avoid: For mechanical processing:  
Avoid formation of dust/air mixtures because of explosion hazard.  
Keep away from sources of ignition. Protect from excessive heat.  
Protect from moisture contamination.
- Incompatible materials: For mechanical processing:  
Strong oxidizing agents, strong reducing agents
- Thermal decomposition: no data available

**11. Toxicological information****Toxicological tests**

- Toxicological effects: Acute toxicity (oral): Lack of data.  
Acute toxicity (dermal): Lack of data.  
Acute toxicity (inhalative): Lack of data.  
Skin corrosion/irritation: Lack of data.  
Eye damage/irritation: Lack of data.  
Sensitisation to the respiratory tract: Lack of data.  
Skin sensitisation: Lack of data.  
Germ cell mutagenicity/Genotoxicity: Lack of data.  
Carcinogenicity: Lack of data.  
Reproductive toxicity: Lack of data.  
Effects on or via lactation: Lack of data.  
Specific target organ toxicity (single exposure): Lack of data.  
Specific target organ toxicity (repeated exposure): Lack of data.  
Aspiration hazard: Lack of data.

**Symptoms**

- In case of inhalation: Dust: May cause irritations.  
Lung damage is possible in a chronic situation.  
After eye contact: Dust: May cause irritations.

**12. Ecological information****Ecotoxicity**

- Further details: no data available

**2/4Rx - Titan Articles**

Material number 002/4Rx

Page: 6 of 7

**Mobility in soil**

no data available

**Persistence and degradability**

Further details: Methods for the determination of biodegradability are not applicable to inorganic substances.

**Additional ecological information**

General information: Do not allow to penetrate into soil, waterbodies or drains.

**13. Disposal considerations****Product**

Recommendation: Special waste. Dispose of waste according to applicable legislation.

**14. Transport information****USA: Department of Transportation (DOT)**

Proper shipping name: Not restricted

**Sea transport (IMDG)**

Proper shipping name: Not restricted

Marine pollutant: no

**Air transport (IATA)**

Proper shipping name: Not restricted

**Further information**

No dangerous good in sense of these transport regulations.

**15. Regulatory information****National regulations - U.S. Federal Regulations**

Refer to 29 CFR Part 1910 subpart q - Welding, Cutting, Brazing 1910.252.


Titanium: TSCA listed

**National regulations - Great Britain**

Hazchem-Code: -

### 16. Other information

Hazard rating systems:



NFPA Hazard Rating:

Health: 0 (Minimal)

Fire: 0 (Minimal)

Reactivity: 0 (Minimal)

HMIS Version III Rating:

Health: 0 (Minimal)

Flammability: 0 (Minimal)

Physical Hazard: 0 (Minimal)

Personal Protection: X = Consult your supervisor

HEALTH	0
FLAMMABILITY	0
PHYSICAL HAZARD	0
	X

Reason of change: Changes in section 1: Address

Date of first version: 8/6/2008

#### Department issuing data sheet

Contact person: see section 1: Dept. responsible for information

The information in this data sheet has been established to our best knowledge and was up-to-date at time of revision. It does not represent a guarantee for the properties of the product described in terms of the legal warranty regulations.