



# Silicone Footshell

Unique, custom shape and color

Ottobock's Custom Silicone Services team offers a custom footshell to match your client's unique color and shape. Featuring 6-to-8 skin tones, the footshell complements the sound side, with multi-colored silicone nails, and the option of a split-toe design. The footshell is fabricated using high-temperature vulcanized silicone and a dyneema toe cap for added strength and durability at the distal end of the carbon foot.

# Silicone Footshell

A custom solution as unique as your client

Because skin color and foot shape are unique to each individual, Ottobock will consult with you to match your color and determine the right shape of the footshell, as much as the prosthetic componentry allows.



## Information on the ordering process

- Start by sending us photos of the sound side foot, based on our color determination guidelines.
- To create the custom shape, a detailed cast of the sound side foot is required.
- Ottobock will fabricate the silicone footshell, according to your specifications, and ship it within 20 working days.
- Within a short period of time, you can fit your client with the new, functional and aesthetically appealing custom silicone footshell.
- Note: The original footshell of the selected foot will be used and discarded during the fabrication process, and the silicone footshell will be fabricated to become the equivalent size of the sound side foot. For optimal results, we encourage you to speak with our custom silicone specialists.

## Warranty

Ottobock offers a one-year warranty on manufacturer defects under normal usage. The life expectancy of the prosthesis will vary depending on the user's daily activity level. Wearing a shoe or sandal on the silicone footshell will prolong the life of the prosthesis.

To obtain a quote or more information regarding the Silicone Footshell (article # 88A35=1), or any other custom silicone product, please call us at 800 665 3327, email us at [SiliconeNA@ottobock.com](mailto:SiliconeNA@ottobock.com), or talk to your Ottobock sales representative today.