

AFO

Order Form



Account Information		Ship To Information	
Date	<input type="text"/>	Name	<input type="text"/>
Account Number	<input type="text"/>	Address	<input type="text"/>
Bill to	<input type="text"/>	City/State/Zip	<input type="text"/>
Phone Number	<input type="text"/>	Phone Number	<input type="text"/>
Email	<input type="text"/>	Email	<input type="text"/>
Buyer	<input type="text"/>	Desired Delivery Date	<input type="text"/>
PO Number	<input type="text"/>		

Please mail your completed Ottobock order form and a negative impression of the patient's limb to the address below. An Ottobock Fabrication Coordinator will contact you.

Shipping Options:

UPS Next Day

UPS Ground

UPS 2-Day

Other

Once this form is complete, please send to Ottobock via fax or email.

Order Quote Only

For clinical questions, call 800 328 4058.

Special Instructions / Comments

Patient Information:

Patient Name _____ Patient Weight _____ lbs. Right _____
 Left Pathology _____
 Activity Level: Low (household ambulator) Moderate (community ambulator) High (unrestricted ambulator)

Select Type of AFO

Standard AFO

Trim: Rigid Standard (trim 1/2" anterior to malleolar apex)
 Semi Rigid (trim at malleolar apex)
 Custom (indicate trimline on negative wrap/cast)

Articulated AFO

Free Motion
 Limited Motion
 Indicate _____

SMO

Floor Reaction

PTB Feature AFO

Bivalve/Clam Shell (Must provide modified cast or model)

Specify material for anterior shell:

Double Upright
 Unilateral
 Check Orthosis

Correction to Cast

Ankle

Leave as casted to 90° (standard)
 to _____° plantarflexion
 to _____° dorsiflexion

Forefoot/Heel

Tone reducing feature
 Valgus/varus to neutral

Material Selection

Laminated*
 Specify color _____
 * Check orthosis required prior to fabrication
 Fred's Legs Pattern _____

Material Selection

Bocklite
 Space-Tech*
 Microcell Puff 3MM or 4 MM
 Bocklite 3MM or 4 MM
 Airflex 3MM

Other _____
 * Only available in 1/8"

Other Options

Liner options (mark all that apply):

- Footplate
- Calf
- Arch Pad
- Full Liner

Malleolus Pad

- Medial
- Lateral

Components

(Ottobock Fabrication Services uses Ottobock components unless otherwise specified):

Ankle Joints _____
 Manufacturer _____
 (Please indicate joint axis and clearances on cast)
 Part # _____

Special Instructions / Comments

Standard AFO supplied with laterally attached 2" hook-and-loop strap, and chafe fastened on the medial side. Extra straps are available for an additional fee.

