

connectgo.pro

connectgo.pro is the new set-up app for O&P professionals to quickly and conveniently program Ottobock microprocessor knees.

Benefits

- Easy and direct connection to the prosthesis for complete set-up, no Bluetooth dongle required
- Modern and intuitive user interface
- Quick access to training videos
- Select and customize MyMode settings to user needs for sports and other activities
- Automatic storage of prosthesis settings
- Export stored settings when needed
- Existing product certifications remain valid and can be transferred for app access
- Designed for use on either a smartphone or tablet
- iOS and Android compatible
- Frequent and convenient updates directly via the App Store and Google Play

MyOttobock

- A MyOttobock account in which certifications are stored makes the app convenient to use on different mobile devices.
- The MyOttobock account is the key to a growing number of new apps, digital features and related services from Ottobock.
- The my.ottobock.com website is the central location for managing the account, personal settings, privacy and available connections.

Availability

- Currently available for *C-Leg* and *Genium X4* microprocessor knees. Additional products will be released through regular updates.
- Data Station PC software will remain available for existing products.

Operating system

- The minimum requirements for *connectgo.pro* are Android 12 and higher and iOS 15 and higher.

