

Transfemoral unique custom liners.

Customer

Contact person _____

Customer no. _____ Date _____

PO # _____

Company _____

Address _____

Postal/zip code _____ City/State/Province _____

Email _____

Phone _____

Shipping address (if different from customer address)

Company _____

Address _____

Postal/zip code _____ City/State/Province _____

Email _____

Shipping Options

UPS Next Day UPS 2-Day
 UPS Ground Other

Patient Information


ID imprinted on custom liner: _____

Full name and/or ID: _____

Affected side: Left Right

Measurements

Extend dimensions as needed. Mark all problem areas on the scan, graphic or plaster.



Important

Regardless of the limb length, the circumference of 50 mm from the distal end in the proximal direction is necessary.

Product Selection

6Y430 Uneo Unique TF (Polyurethane)

6Y730 Skeo Unique TF (Silicone)

New Liner

Replacement Liner (Serial No. _____)

Wall thickness tolerances of ± 1 mm may occur

Locking Mechanism

Yes*: Pin KISS DVS

No

*Locking Liners must have exterior textile selected

Textile

Uncovered

Partial Cover: Bonding distal to proximal ____cm

Full Textile Cover

Textile Color: Anthracite Grey Tan

For liners with textile cover, the wall thickness increases by 1 mm

Wall Thickness

Tapered*

Uniform: 4 mm 5 mm 6 mm

*Minimal residual limb length of >20 cm

*Uneo Unique: From 6 mm distal to 3 mm proximal

*Skeo Unique: From 5 mm distal to 3 mm proximal

Manufacturing tolerance is ± 1 mm

Distal Cushion

Standard 13 mm

7 mm 9 mm 11 mm 15 mm

17 mm 19 mm 21 mm

For liners with distal connection, the distal wall thickness increases by 10 mm

Additives (Uneo polyurethane liners only - Select one or none)

With Skinguard - antibacterial additive

With Fresh scent

Questions or looking to order a replacement liner? Contact Liner Customer Service at 512-806-2596 or USSLLiners@Ottobock.com

Send scan files to USSLCCad@Ottobock.com or submit via iCC.

Please ship casts to our Fabrication Facility: 3820 West Great Lakes Drive, Salt Lake City, UT 84120

Ottobock maintains all rights, title, and ownership to the custom liner mold and will keep the mold on file for 2 years from the last order date.