

617Sx1 - PU-flexible foam

Material number 617Sx1

Page: 1 of 7

1. Product and company identification**Product identifier**

Trade name: 617Sx1 - PU-flexible foam

Relevant identified uses of the substance or mixture and uses advised against

General use: Article for orthopedic procedures
Pads:
617S119, 617S67, 617S68
Foam cover:
3R24, 3R48, 3R6, 3S106, 3S107, 3S124, 3S26, 6R6=B, 6R7

Details of the supplier of the safety data sheet

Company name: Otto Bock Health Care
Street/POB-No.: 3820 W. Great Lakes Drive
Postal Code, city: Salt Lake City, UT 84120
USA
WWW: www.ottobockus.com
Telephone: +1 (801) 956-2400
Telefax: +1 (801) 956-2401
Dept. responsible for information:
Quality Department,
Telephone: +1 (801) 954-2304 (7 AM – 3 PM, Mountain Time), Email:
USRegulatory@ottobock.com

Additional information: This safety data sheet pertains to the following products:
3R24 - PU-Flexible Foam
3R48 - PU-Flexible Foam
3R6 - PU-Flexible Foam
3S106 - PU-Flexible Foam
3S107 - PU-Flexible Foam
3S124 - PU-Flexible Foam
3S26 - PU-Flexible Foam
6R6=B - PU-Flexible Foam
6R7 - PU-Flexible Foam
617S119 - Dyatec
617S67 - PPT, Perforated
617S68 - PPT

Emergency phone number**CHEMTREC, Telephone: +1 (800) 424-9300****2. Hazards identification****Emergency overview**

Appearance: Form: solid
Color: varying colors,
flesh-colored:
3R24, 3R48, 3R6, 3S106, 3S107, 3S124, 3S26, 6R6=B, 6R7, 617S67, 617S68
red:
617S119
Odor: weak characteristic

617Sx1 - PU-flexible foam

Material number 617Sx1

Page: 2 of 7

Classification: This substance is classified as not hazardous.

Regulatory status

This material is not considered hazardous by the U.S. OSHA Hazard Communication Standard (29 CFR 1910.1200) and SIMDUT in Canada.

Hazards not otherwise classified

Processing, e.g. by cutting, sawing or grinding, can produce particles and dust. For risks which have to be observed thereby, see chapter 7: Handling, chapter 8: Exposure controls / personal protection and chapter 11: Toxicology.

Handling and/or processing of this material may generate a dust which can cause mechanical irritation of the eyes, skin, nose and throat.

see section 11: Toxicological information

3. Composition / Information on ingredients

Chemical characterization: Article: polyurethane-flexible foam

CAS-Number: 9009-54-5

RTECS-Number: TR7875000

4. First aid measures

In case of inhalation: In case of troubles after inhalation of dust:
Move victim to fresh air. Seek medical attention.

After eye contact: Polyurethane-dust:
Immediately flush eyes with plenty of flowing water for 10 to 15 minutes holding eyelids apart. Subsequently consult an ophthalmologist.

After swallowing: Swallowing is not regarded as a possible way of exposition.
Polyurethane-dust:
If person is clearly conscious, have them drink two glasses of water to dilute ingested material.
Seek medical attention.

Most important symptoms/effects, acute and delayed

Handling and/or processing of this material may generate a dust which can cause mechanical irritation of the eyes, skin, nose and throat.

Information to physician

Treat symptomatically.

5. Fire fighting measures

Flash point/flash point range:

no data available

Auto-ignition temperature: not self-igniting

Suitable extinguishing media:

Water fog, Extinguishing powder, carbon dioxide.

Extinguishing media which must not be used for safety reasons:

High power water jet

617Sx1 - PU-flexible foam

Material number 617Sx1

Page: 3 of 7

Specific hazards arising from the chemical

In case of fire may be liberated:
Nitrous fumes, carbon monoxide and carbon dioxide.

Protective equipment and precautions for firefighters:

Wear self-contained positive pressure breathing apparatus and full firefighting protective clothing.

Additional information: NFPA Sprinkler Classification: Extra Hazard

6. Accidental release measures

Personal precautions: Avoid generation of dust. Do not breathe dust.
In the case of the formation of dust: Wear protective equipment.

Environmental precautions: Do not allow to penetrate into soil, waterbodies or drains.

Methods for clean-up: Polyurethane-dust:
Take up mechanically, placing in appropriate containers for disposal.
Avoid generation of dust.

7. Handling and storage**Handling**

Advices on safe handling: For mechanical processing:
Avoid respiration of swarf. Wear protective equipment.
Provide adequate ventilation.

Precautions against fire and explosion:
Keep away from sources of ignition. - No smoking.
Take precautionary measures against static discharges.

Storage

Requirements for storerooms and containers:
Store at room temperature. Keep container dry.

8. Exposure controls / personal protection**Exposure guidelines**

Occupational exposure limit values:

Type	Limit value
USA: ACGIH: TWA	10 mg/m ³
USA: ACGIH: TWA	3 mg/m ³
USA: OSHA: TWA	15 mg/m ³
USA: OSHA: TWA	5 mg/m ³

Engineering controls

Provide good ventilation and/or an exhaust system in the work area.
See also information in chapter 7, section storage.

617Sx1 - PU-flexible foam

Material number 617Sx1

Personal protection equipment (PPE)

- Eye/face protection For mechanical processing:
Tightly sealed goggles according to OSHA Standard - 29 CFR: 1910.133 or ANSI Z87.1-2010.
- Skin protection Suitable protective clothing.
For mechanical processing:
Protective gloves according to OSHA Standard - 29 CFR: 1910.138.
Observe glove manufacturer's instructions concerning penetrability and breakthrough time.
- Respiratory protection: For mechanical processing: Respiratory protection must be worn whenever the TLV (WEL) levels have been exceeded. Dust mask filter P according to OSHA Standard - 29 CFR: 1910.134 or ANSI Z88.2
- General hygiene considerations:
For mechanical processing:
Do not breathe dust.
Wash hands before breaks and after work.
Provide a conveniently located eye rinse station.

9. Physical and chemical properties

Information on basic physical and chemical properties

- Appearance: Form: solid
Color: varying colors,
flesh-colored:
3R24, 3R48, 3R6, 3S106, 3S107, 3S124, 3S26, 6R6=B, 6R7, 617S67,
617S68
red:
617S119
- Odor: weak characteristic
- Odor threshold: no data available
- pH value: no data available
- Melting point/freezing point: no data available
- Initial boiling point and boiling range: no data available
- Flash point/flash point range: no data available
- Evaporation rate: no data available
- Flammability: no data available
- Explosion limits: no data available
- Vapor pressure: no data available
- Vapor density: no data available
- Density: no data available
- Water solubility: at 68 °F: insoluble
- Partition coefficient: n-octanol/water: no data available
- Auto-ignition temperature: not self-igniting
- Thermal decomposition: no data available
- Explosive properties: Product is not explosive.

617Sx1 - PU-flexible foam

Material number 617Sx1

Page: 5 of 7

10. Stability and reactivity

- Reactivity: no data available
- Chemical stability: Product is stable under normal storage conditions.
- Possibility of hazardous reactions
No hazardous reaction when handled and stored according to provisions.
- Conditions to avoid: Keep away from sources of ignition. - No smoking.
Take precautionary measures against static discharges.
- Incompatible materials: no data available
- Hazardous decomposition products:
In case of fire may be liberated:
Nitrous fumes, carbon monoxide and carbon dioxide.
Not readily flammable according to DIN 75200:
3R24, 3R48, 3R6, 3S106, 3S107, 3S124, 3S26, 6R6=B, 6R7
- Thermal decomposition: no data available

11. Toxicological information**Toxicological tests**

- Acute toxicity: LD50 Rat, oral: > 5000 mg/kg
- Toxicological effects: Acute toxicity (oral): Lack of data.
Acute toxicity (dermal): Lack of data.
Acute toxicity (inhalative): Lack of data.
Skin corrosion/irritation: Lack of data.
Eye damage/irritation: Lack of data.
Sensitisation to the respiratory tract: Lack of data.
Skin sensitisation: Lack of data.
Germ cell mutagenicity/Genotoxicity: Lack of data.
Carcinogenicity: Lack of data.
Reproductive toxicity: Lack of data.
Effects on or via lactation: Lack of data.
Specific target organ toxicity (single exposure): Lack of data.
Specific target organ toxicity (repeated exposure): Lack of data.
Aspiration hazard: Lack of data.
- In case of inhalation: Polyurethane-dust: May cause irritations.
- After swallowing: Polyurethane-dust: After intake of large amounts: nausea, vomiting.
- After eye contact: Polyurethane-dust: May cause irritations.

Symptoms

Handling and/or processing of this material may generate a dust which can cause mechanical irritation of the eyes, skin, nose and throat.

617Sx1 - PU-flexible foam

Material number 617Sx1

Page: 6 of 7

12. Ecological information**Ecotoxicity**

Further details: no data available

Mobility in soil

no data available

Persistence and degradability

Further details: no data available

Additional ecological information

General information: Do not allow to penetrate into soil, waterbodies or drains.

13. Disposal considerations**Product**

Recommendation: Incinerate according to applicable local, state and federal regulations.

Contaminated packaging

Recommendation: Dispose of waste according to applicable legislation.

14. Transport information**USA: Department of Transportation (DOT)**

Proper shipping name: Not restricted

Sea transport (IMDG)

Proper shipping name: Not restricted

Marine pollutant: no

Air transport (IATA)

Proper shipping name: Not restricted

Further information

No dangerous good in sense of these transport regulations.

15. Regulatory information**National regulations - U.S. Federal Regulations**

TSCA Inventory: not listed

Carcinogen Status:

IARC Rating: Group 3

OSHA Carcinogen: not listed

NTP Rating: not listed

National regulations - Great Britain

Hazchem-Code: -

617Sx1 - PU-flexible foam

Material number 617Sx1

Page: 7 of 7

16. Other information

Hazard rating systems:



NFPA Hazard Rating:

Health: 0 (Minimal)

Fire: 1 (Slight)

Reactivity: 0 (Minimal)

HMIS Version III Rating:

Health: 0 (Minimal)

Flammability: 1 (Slight)

Physical Hazard: 0 (Minimal)

Personal Protection: X = Consult your supervisor

HEALTH	0
FLAMMABILITY	1
PHYSICAL HAZARD	0
	X

Reason of change: Changes in section 1: Address

Date of first version: 9/24/2008

Department issuing data sheet

Contact person: see section 1: Dept. responsible for information

The information in this data sheet has been established to our best knowledge and was up-to-date at time of revision. It does not represent a guarantee for the properties of the product described in terms of the legal warranty regulations.