

ottobock.

I am unique.

The Unique Liner Family



Quality for life

Information for practitioners

The right solution is Unique. As unique as your patient's needs.

With Unique Liners from Ottobock you can offer highly customized concepts and thus achieve innovative fitting results.



- 1** Easy Scanning with iPad, Structure scanner and Custom4U App
- 2** Once scanned simply fill in the order form and e-mail via the App for image processing at iFab
- 3** Individualized mold is created and liner is fabricated and shipped back to you.

“My responsibility is to make my contribution in developing a quality product.”
Stephanie, Production Operator



Ordering process:
You can order the liner simply by using the Ottobock Custom4U-App. Furthermore you can generate a plaster cast and send it to the Ottobock production center.

Download on App Store.

Uniqueness optimized.

The next level of individual fitting.

In order to gain an optimal result, the liner may need to be cut precisely to the accurate size. As soon as the right fit and the right amount of protection is achieved, your patient is ready to go.



“The liner is created specific to the patient’s body. Thus I get the best possible fitting for my patient.”

Jens, CPO

1 Skeo Unique



2 Uneo Unique



Unique Liners are especially helpful for short residual limbs, proximal femoral focal deficiency, symes, knee disartic, and uneven distal ends. The Uneo Unique urethane liner has been well accepted for over 20 years and now the Skeo Unique silicone liner expands our custom portfolio.



1 Skeo Unique

Being a part of the tried and tested Skeo family, Skeo Unique is recommended when high stability and durability are required. The liner is appropriate for very conical or scarred residual limbs. Furthermore the Skeo Unique Locking Liners each have an integrated matrix which minimizes pistoning during the swing phase.

- Stabilizes residual limbs, even those with significant soft tissue
- Skin-friendly silicone enables a superior connection with the residual limb
- Easy handling, easy cleaning

2 Uneo Unique

Thanks to the extraordinary flow characteristics the Uneo Unique achieves an exact and comfortable fit. Due to the qualities of the polyurethane gel, Uneo Unique will cushion even highly sensitive or scarred residual limbs.

- Good pressure distribution within the socket
- Comfortable impact absorption
- Skinguard antibacterial or Fresh scent options

- Black
- Tan
- Silver
- Recommended solution
- Possible solution

Order options

	Skeo Unique 6Y700	Uneo Unique 6Y400
Wall thickness		
Uniform	4.5 mm uniform distal wall to proximal	6 mm uniform distal wall to proximal
Tapered	4.5 mm MPT to 2.5 mm proximal	6 mm MPT to 3 mm proximal
Distal thickness		
Locking	21 mm	18 mm
Cushion	13 mm	13 mm
Locking mechanisms		
KISS	TF	-
Shuttle Lock	TT+TF	TT+TF
Dynamic Vacuum System	TT	TT

Exterior

Uncovered / Partial Cover / Full Textile Cover

Lykra Fabric 0.6 mm	● X	X
Wearforce 1.6 mm	● X	X
Silver Fabric 1.0 mm	● X	X
Additives		
Fresh (scented)	-	X
Skinguard (antibacterial)	-	X

Fitting opportunities

	Skeo Unique 6Y700		Uneo Unique 6Y400	
	TF	TT	TF	TT
KISS Lanyard System	■	-	□	-
Shuttle Lock System	■	■	□	□
Valve System	■	□	□	■
Dynamic Vacuum System	-	□	-	■
Harmony (P3 & P4)	□	□	□	■
Harmony (E2)	■	□	□	■

Step by step to a new quality of life. Unique Liners prove their worth in everyday use.

“My experience with the Unique Liner has exceeded my expectations. I was not sure if the new prosthesis could naturally fit in my everyday life. But after using the Skeo Unique for a while now, it pays off in every detail.”

Doreen



One of the greatest benefits of the Unique Liner is the exceptional fit. This liner can help your patients focus on getting back to their everyday life.



A crucial part in enjoying an active lifestyle is the reliability and secure feeling from the prosthesis.



