

**ottobock.**



***Varos.***

Right from the start.



László owns a building company, so after he lost his leg in an accident, it was essential to return to his normal life quickly and with a prosthesis he could really trust. He faces new challenges at work every day. He says, “I walk a lot and have to manage many different obstacles at building sites.” László was lucky and was fitted with a prosthesis that gave him a great deal of safety and independence right from the start.



“

*My technician found a solution for me that allowed me to quickly return to my normal life.”*

The *Varos Adjustable Socket* makes his daily routine easier. Just a few twists of the adjustment dials and the socket fits perfect. This allows him to enjoy even longer walks with his dogs. László is relieved: “In the beginning, I could only walk for limited distances, but now I’m on the go from morning to evening again.”



### **Varos Liner**

- Magnetic connection to socket
- Comfortable gel material
- Easy to put on and take off



### **Varos Adjustable Socket**

- Top and bottom fit can be adjusted independently by the user
- Very comfortable



### **Rotation adapter**

- Easy to put on and remove shoes/clothes
- Relaxed sitting



### **Kenevo - Microprocessor Knee**

- Safety with every step
- Custom adaptation to the user's needs



### **Terion K2 foot**

- Safe heel strike
- Comfortable rollover



# Laszlo's Innovative Socket Solution.

## **Independent. Safe. Comfortable.**

The *Varos* is an adjustable socket solution that benefits transfemoral amputees starting in the rehabilitation phase as well as later during treatment with the final prosthesis. The socket circumference can be easily adjusted at any time using two adjustment wheels. This is important because the circumference of the residual limb changes, for example due to differences in temperature and also depending on how long a prosthesis is worn during the day. Various activities are also easier to manage if the socket is somewhat tighter or looser.

The flexible shells of the *Varos* socket and *Varos* liner work in combination for easy donning and doffing, even when sitting.

Compared with a conventional socket, the *Varos* socket can be quickly custom fitted with minimal visits to the technician.

For more information about *Varos* please visit:

**[www.ottobock.com](http://www.ottobock.com)**