

## F23 *Maverick* *Vertical Shock.*

Trust the next step.

The *Maverick Vertical Shock* was designed for highly active users who value flexibility and high durability when performing high-impact activities.

### Product features

- Vertical shock unit offers higher comfort and flexibility with its torsion range and shock absorption.
- Effectively absorbs torsion (+/- 10°) to relieve the residual limb and boost comfort in everyday life in combination with vertical shock absorption (up to 15 mm).
- Innovative heel/spring design delivers a high degree of durability and flexibility without limiting energy return
- Split forefoot and rearfoot springs for enhanced adaptation to terrain.
- A full length, unbolted sole plate ensures a smooth stance phase rollover.
- Weatherproof design allows for more reliable performance in harsh conditions.\*

\*Please reference the *Maverick Vertical Shock* IFU for weatherproof details.

More than a foot.  
A foundation.

### Technical data

K-Level / Mobility grade	3, 4
Max. body weight	324 lbs / 147 kg
Side	left (L), right (R)
Sizes	23-30 cm
Footshell color	beige, light brown
Weight with footshell**	37.81 oz / 1072 g
Heel height	0.393" / 10 mm
Build height**	8.543" / 217 mm

\*\* Reference size is 26 cm.



# F23 Order form.

## Customer

Customer no.		
Company		
Street		
Postal code	City	State
O&P professional		

## Shipping address (if different from customer address)

Customer no.		
Company		
Street		
Postal code	City	State
Customer ref.		

### Selection of the spring stiffness relative to 1. body weight and activity as well as

Body weight	Activity level		
	Low	Moderate	High
100 – 115 lbs / 44 – 52 kg	1	1	2
116 – 130 lbs / 53 – 59 kg	1	2	3
131 – 150 lbs / 60 – 68 kg	2	3	4
151 – 170 lbs / 69 – 77 kg	3	4	5
171 – 195 lbs / 78 – 88 kg	4	5	6
196 – 220 lbs / 89 – 100 kg	5	6	7
221 – 255 lbs / 101 – 116 kg	6	7	8
256 – 285 lbs / 117 – 130 kg	7	8	9
286 – 325 lbs / 131 – 147 kg	8	9	-

### 2. foot size

Foot size	Stiffness								
	1	2	3	4	5	6	7	8	9
23							-	-	-
24								-	-
25								-	-
26	-								-
27	-	-							
28	-	-							
29	-	-	-						
30	-	-	-						

### Example for ordering F23=L26-4-P/4

Quantity	Article no.	Side	Size	Stiffness	Color
	F23	=		-	P /
		Side	Size (cm)	Stiffness	Color
		R Right	23, 24, ..., 30	1-9	4 Beige
		L Left			15 Light brown

### Alternative ordering option - Article number foot module (without footshell)

Quantity	Article no.	Stiffness	Size	Side
	F23	- 00 - 0	A	- S
		Stiffness	Size (cm)	Side
		1-9	23, 24, ..., 30	R Right
				L Left

### Alternative ordering option - Article number footshell

Quantity	Article no.	Size	Color	Side
	FTC	- 2 A - 1		4 - S
		Size (cm)	Color	Seite
		23, 24, ..., 30	L Beige	R Right
			D Light brown	L Left

### Scope of delivery

The scope of delivery includes the F23 foot module with an attached functional ring and pre-compression kit, FTC-2A-1\* footshell, 2F11=1 heel wedge set, one standalone functional ring and a black Spectra Sock. Additionally, a tool for changing the footshell can be ordered (ACC-00-10300-00).

