

ottobock.

## VS5 Restore. Your reliable companion.

The *Restore* was designed for less active users who mainly navigate indoor environments and use a walking aid outdoors, and who place a high value on reliable support from their prosthetic foot.

### Product features

- The full length, unbolted spring design provides the required stability and a soft rollover
- The C-shape offers stability from heel contact through to toe-off
- The partly split forefoot spring and the flexibility of the fibreglass provide sufficient ground contact even on uneven terrain
- Custom heel stiffness for a better balance thanks to four different heel wedge options
- Protected against fresh, salt and chlorinated water

More than a foot.  
A foundation.

### Technical data

K-Level / Mobility grade	1, 2
Max. body weight	300 lbs / 136 kg
Side	left (L), right (R)
Sizes	22-30 cm
Footshell colour	beige, light brown
Weight with footshell*	23.77 oz / 674 g
Heel height	3/8" / 10 mm
Build height*	3 1/2" / 89 mm

\* Reference size is 26 cm.

[www.ottobock.com](http://www.ottobock.com)

US: P 800 328 4058 | F 800 962 2549  
CA: P 800 665 3327 | F 800 463 3659



## VS5 Restore. Order form.

### Customer

Customer no.	
Company	
Street	
Postal code	City
O&P professional	

### Shipping address (if different from customer address)

Customer no.	
Company	
Street	
Postal code	City
Customer ref.	

### Selection of the spring stiffness relative to

1. body weight and activity as well as

Body weight (kg)	Stiffness
45 – 59	1
60 – 81	2
82 – 109	3
110 – 136	4

2. foot size

Stiffness	1	2	3	4
	Foot size			
22			-	-
23				-
24				
25				
26				
27				
28	-			
29	-			
30	-	-		

### Example for ordering VS5=L26-3-P/4

Quantity	Article no.	Side	Size	Stiffness	Colour
	VS5	=		-	P /

  

Side	Size (cm)	Stiffness	Colour
R Right	22, 23, ..., 30	1-4	4 Beige
L Left			15 Light brown

### Alternative ordering option - Article number foot module (without footshell)

Quantity	Article no.	Stiffness	Size	Side
	VS5	- 00 - 0	A	- S

  

Stiffness	Size (cm)	Side
1-4	22, 23, ..., 30	R Right
		L Left

### Alternative ordering option - Article number footshell

Quantity	Article no.	Size	Colour	Side
	FTC	- 2 B - 1		4 - S

  

Size (cm)	Colour	Side
22, 23, ..., 30	L Beige	R Right
	D Light brown	L Left

### Scope of delivery

The scope of delivery includes the VS5\* foot module, the FTC-2B-1\* footshell, the 2F13=\* heel wedge set (comprising four different degrees of hardness) and a black Spectra sock. Additionally, a tool for changing the footshell can be ordered (ACC-00-10300-00).



### Compatible mechatronic knee joints



Kenevo

[www.ottobock.com](http://www.ottobock.com)

US: P 800 328 4058 | F 800 962 2549  
CA: P 800 665 3327 | F 800 463 3659