




BionicLink PC compatibility.

			
		current/blue	new/white
Product Name	Version	BionicLink PC 60X5	BionicLink PC LE 60X5=2
C-Leg	Purchased prior to August 2023*	✓	✓
	Purchased after August 2023*	✗	✓
Kenevo	Purchased prior to November 2023*	✓	✓
	Purchased after November 2023*	✗	✓
Genium & Genium X3	All versions	✓	✓
Meridium	All versions	✓	✓
Dynamic Arm	All versions	✓	✗
Michelangelo/AXON	All versions	✓	✓
System 2000	All versions	✓	✓

To identify if the new BionicLink PC adapter is compatible with your C-Leg or Kenevo, check:

- if the article number ends in “-3”
- if the inscription on the hydraulic unit reads “LE DUAL”



*C-Legs and Kenevos produced using the newer LE-DUAL (low energy) Bluetooth chip will only pair with the new 60X5=2 (white) BionicLink PC adapter.