

## 6Y112 Skeo Sealing TT.

The 6Y112 *Skeo Sealing TT* liner is the latest addition to our Skeo Sealing Liner family, offering transtibial amputees a high level of comfort and reliability. The sealing ring provides a consistent seal against the inner socket wall for secure suspension. Its silky-smooth surface with no textile cover simplifies cleaning as well as the donning and doffing process.

### Designed for optimal outcomes

- One sealing height provides a higher vacuum area for higher adhesion as well as security and peace-of-mind for the patient
- Integrated sealing ring promotes durability throughout the lifetime of the liner
- The inside of the liner is textured for a pleasant feeling on the skin
- Silky-smooth outer coating for easier donning and doffing
- Orientation aid on the inside of the liner for an optimal liner position on the residual limb
- An integrated matrix for reduced elongation and good tissue management
- The *Skeo Sealing TT* can be combined with the NEW 21Y2 *QuickValve* or the *Harmony* system plus 21Y2 *QuickValve*.
- Waterproof: Protected against fresh, salt and chlorinated water



# 6Y112

## Skeo Sealing TT.

Article number	Circumference	Ring height
6Y112=160X80	160 mm	80 mm
6Y112=180X80	180 mm	80 mm
6Y112=200X80	200 mm	80 mm
6Y112=210X80	210 mm	80 mm
6Y112=220X80	220 mm	80 mm
6Y112=235X80	235 mm	80 mm
6Y112=250X80	250 mm	80 mm
6Y112=265X80	265 mm	80 mm
6Y112=280X80	280 mm	80 mm
6Y112=300X80	300 mm	80 mm
6Y112=320X80	320 mm	80 mm
6Y112=340X90	340 mm	90 mm
6Y112=360X90	360 mm	90 mm
6Y112=380X90	380 mm	90 mm
6Y112=400X90	400 mm	90 mm



The ring



The material



The orientation aid