

movido.

For young adventurers.

Thanks to its very low build height and very low weight, **movido** is specially designed for small children. With its hydraulic swing phase control, **movido** supports a more physiological gait pattern right from the start and gives children the freedom of movement they need.

Designed for optimal outcomes

- Low build height
- Low weight
- The 4-axis joint geometry provides stability during the stance phase and enables easy swing phase release
- The hydraulic swing phase control supports different walking speeds
- Comfortable knee extension without a hard (audible) stop
- The high knee flexion angle of up to 154° makes kneeling easier
- The included knee protection cap and dust cover protect the knee joint, clothing and floor when kneeling
- Water resistant



| Technical data | 3R68=1-P movido explore | 3R68=2-P movido play |
|--|------------------------------|-------------------------|
| For | Toddlers and school children | |
| Max. body weight | 55 lbs (25 kg) | 99 lbs (45 kg) |
| Weight | 6.7 oz (190 g) | 8.8 oz (250 g) |
| Distal connection | Tube clamp ø 22 mm | |
| Proximal connection | Pyramid | |
| Build height | 4" (102 mm) | 4 ½" (113 mm) |
| Distal build height to alignment reference point | 3 ½" (88 mm) | 3 ⅞" (99 mm) |
| Proximal build height to alignment reference point | 7/8" (21 mm) | 15/16" (23 mm) |
| Knee flexion angle | 145° | 154° |



Very compact design
Easy kneeling
More physiological gait pattern

3R68 Order form.

Customer

| | |
|------------------|------|
| Customer no. | |
| Company | |
| Street | |
| Postal code | City |
| O&P professional | |

Shipping address (if different from customer address)

| | |
|---------------|------|
| Customer no. | |
| Company | |
| Street | |
| Postal code | City |
| Customer ref. | |

Example for ordering 3R68=1-P

| Article no. | | Size | | |
|-------------|---|------|---|---|
| 3R68 | = | | - | P |

Size

| | |
|---|---------------------------|
| 1 | movido explore (< 55 lbs) |
| 2 | movido play (< 99 lbs) |



The scope of delivery includes the 3R68 knee joint, a knee protection cap (4P10=1 or =2) and a dust cover (4P11=1).

Recommended combinations*



4R110 Lamination anchor with pyramid receiver or 4R60 Socket adapter
4R11 Quickchange
2R41 Tube adapter



F24 Maverick junior

* All components are sold separately and are available Ottobock products that are compatible with the 3R68 **movido**, which helps ensure optimal performance. Practitioners need to select components based upon individual patient criteria.

www.ottobock.com

US: P 800 328 4058 | F 800 962 2549

CA: P 800 665 3327 | F 800 463 3659